---

# A 105/1

<b>de</b>	Gebrauchsanweisung Injektoroberkorb
<b>da</b>	Brugsanvisning Injektoroverkurv
<b>en</b>	Operating instructions Upper injector basket
<b>es</b>	Instrucciones de manejo Cesto superior del inyector
<b>fi</b>	Käyttöohje Suorasuihkuyläkori
<b>fr</b>	Mode d'emploi Panier à injection
<b>hr</b>	Upute za uporabu Injektorska gornja košara
<b>it</b>	Istruzioni d'uso cesto superiore a iniezione
<b>nl</b>	Gebruiksaanwijzing injectorbovenrek
<b>no</b>	Bruksanvisning overkurv/injektor
<b>pt</b>	Instruções de utilização Cesto superior injetor
<b>sl</b>	Navodila za uporabo Injektorska košara
<b>sv</b>	Bruksanvisning Injektoröverkorg

**de** ..... 4  
**da** ..... 19  
**en** ..... 34  
**es** ..... 49  
**fi** ..... 64  
**fr** ..... 79  
**hr** ..... 94  
**it** .....109  
**nl** .....125  
**no** .....140  
**pt** .....155  
**sl** .....170  
**sv** .....185

---

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<b>Notes about these instructions</b> .....	34
Queries and technical problems.....	34
<b>Intended use</b> .....	35
IMPORTANT for the U.S. ....	35
<b>Items supplied</b> .....	36
Scope of delivery A 105/1 .....	36
Disposal of the packing material .....	36
Optional accessories.....	37
<b>Warning and Safety instructions</b> .....	38
<b>Areas of application</b> .....	39
Adjusting the upper basket .....	39
Screwing in nozzles .....	41
Securing injector nozzles .....	41
Test point for measuring water pressure .....	41
Preparing the load .....	42
Before loading the machine and before starting a programme .....	42
Reprocessing hollow instruments .....	42
A 803 holder for transmission instruments.....	43
Installation .....	43
Reprocessing transmission instruments .....	44
AUF 1 holder for transmission instruments .....	44
Optional accessories.....	44
Installation .....	45
Reprocessing transmission instruments .....	45
Replacing filter plates .....	46
A 814 holder for the tips of air scalers and piezo scalers .....	46
Fitting E 499 for suction and rinsing cannulae .....	47
Irrigation sleeves .....	47

### Important warnings

 Information which is important for safety is highlighted in a thick framed box with a warning symbol. This alerts you to the potential danger of injury to people or damage to property. Read these warning notes carefully and observe the procedural instructions and codes of practice they describe.

### Notes

Information of particular importance that must be observed is highlighted in a thick framed box.

### Additional information and comments

Additional information and comments are contained in a simple frame.

### Operating steps

Operating steps are indicated by a black square bullet point.

#### Example:

- Select an option using the arrow buttons and save your choice with *OK*.

### Display

Certain functions are shown in display messages using the same font as that used for the function itself in the display.

#### Example:

Settings  menu.

### Queries and technical problems

In the event of queries or technical problems, please contact Miele. Contact details can be found at the end of the operating instructions for your cleaning machine or at [www.miele.com/professional](http://www.miele.com/professional).

This upper basket can be used to reprocess machine-reprocessable medical devices in special Miele cleaning machines and washer-disinfectors. Follow the operating instructions for your machine as well as the instructions provided by the medical device manufacturer on how to reprocess their items by machine.

The A 105/1 upper injector basket is suitable for reprocessing dental and ENT instruments by machine.

The upper basket has two injector manifolds with a total of 22 holders for nozzles and other irrigation connectors designed for the reprocessing of lumened instruments.

The rest of the upper basket is free and can be equipped with inserts and mesh trays for reprocessing instruments without lumen.

Special reprocessing programmes are required for reprocessing instruments from the above areas of application by machine. These programmes are installed ex-works or can be installed, created or activated by the Miele Customer Service Department on the following machines:

- PG 8581
- PG 8582
- PG 8582 CD
- PG 8591
- PG 8592

Special fittings such as A 803 or AUF 1 are needed for reprocessing transmission instruments used in dental applications.


Disposable items must not be reprocessed.

The washer-disinfector is generally referred to as “the machine” in these operating instructions. Reprocessable medical devices are referred to as “items” if they are not more closely defined.

### **IMPORTANT for the U.S.**

Not all sections in this operating manual reflect the Intended Use for the USA. Please refer to the operating manual of the Miele Dental Washer-Disinfector or contact Miele Service Professional at +1 800-991-9380 with application questions.

### Scope of delivery A 105/1

- Upper injector basket A 105/1
  - Height 202 mm, width 528 mm, depth 527 mm
  - Vertical clearance 215 mm ( $\pm 30$  mm)
  - Injector manifold vertical clearance 185 mm ( $\pm 30$  mm)
  - Spray arm with magnets for spray arm monitoring
- A 838, tool set, consisting of 1 x WAF 9 spanner and 1 x T20 offset screwdriver

### **Disposal of the packing material**

The packaging is designed to protect against transportation damage. The packaging materials used are selected from materials which are environmentally friendly for disposal and should be recycled.

Recycling the packaging reduces the use of raw materials in the manufacturing process and also reduces the amount of waste in landfill sites.

### Optional accessories

- A 800, reusable tubular filter with a filtration grade of  $\geq 0.07$  mm, 405 mm long,  $\varnothing$  22.5 mm
- A 801, threaded adapter
- A 803, holder for dental transmission instruments, height 76.5 mm, width 40 mm, depth 32 mm
- A 813, two-piece adapter with 4 screw threads for various irrigation connectors, height approx. 48 mm,  $\varnothing$  45 mm, 4 nozzles with silicone hoses and 3 blind stoppers
- A 814, holder with screw cap for reprocessing the tips of air scalers and piezo scalers, height approx. 41 mm,  $\varnothing$  16 mm
- A 836, nozzle set for e.g. dental applications, consisting of 22 x A 833, injector nozzle with aperture on the side, length 110 mm,  $\varnothing$  4 mm
- A 837, irrigation sleeve set for e.g. ENT applications, consisting of 22 x A 817, irrigation sleeve for instruments  $\varnothing$  4–8 mm, length 67 mm,  $\varnothing$  11 mm, cap with aperture  $\varnothing$  6 mm
- AUF 1, two-piece adapter for transmission instruments (dental handpieces, for example), height approx. 48 mm,  $\varnothing$  45 mm; silicone adapter (ADS) required for reprocessing
- E 336, irrigation sleeve, length 121 mm,  $\varnothing$  11 mm
- E 442, irrigation sleeves for instruments  $\varnothing$  4 - 8 mm, 121 mm long,  $\varnothing$  11 mm, cap with aperture  $\varnothing$  6 mm
- E 443, irrigation sleeves for instruments  $\varnothing$  8 - 8.5 mm, 121 mm long,  $\varnothing$  11 mm, cap with aperture  $\varnothing$  10 mm
- E 445, 12 caps for irrigation sleeves, aperture  $\varnothing$  6 mm
- E 446, 12 caps for irrigation sleeves, aperture  $\varnothing$  10 mm
- E 447, Luer Lock Adapter ♀ for Luer Lock ♂
- E 499 for reprocessing suction and rinsing cannulae with a variety of connector types as well as Veress cannulae.
- FP 20, 20 filter plates,  $\varnothing$  30 mm

Other optional accessories are available to order from Miele.

## en - Warning and Safety instructions

---

Read the operating instructions carefully before using this upper basket.  
Keep these operating instructions in a safe place.

 Please also read the operating instructions for your washer-disinfector and pay particular attention to the Warning and Safety instructions.

- ▶ The upper basket is only approved for use in the applications stated in these operating instructions. Components such as nozzles may only be replaced with Miele accessories or genuine original spare parts.
- ▶ New mobile units and baskets must be washed in the washer-disinfector without a wash load prior to first use.
- ▶ Inspect all mobile units, baskets, modules and inserts as described in the “Maintenance” section of the operating instructions for your cleaning machine.
- ▶ For your own safety, wear gloves and, if necessary, a surgical mask and protective goggles when sorting the instruments and unloading the machine. To avoid recontamination, wear clean gloves when unloading the machine.
- ▶ Load from back to front and unload from front to back.
- ▶ To clean and disinfect lumened instruments thoroughly from the inside, they must be attached to the relevant holders or nozzles on the injector manifolds. The connectors must not come loose during a programme sequence.  
At the end of every programme sequence, check that the connectors have not come loose.
- ▶ Any hollow items that have become disconnected from their injector nozzles, adapters or holders during processing must be re-processed.
- ▶ The last rinse, as a minimum, must be carried out with DI water.

### Adjusting the upper basket

The A 105/1 upper basket can be adjusted between three positions with approx. 30 mm between each position to accommodate items of different heights.

To adjust the height, the brackets with rollers on the side of the upper basket and the water connector at the back of the basket have to be moved. The roller brackets are each secured to the upper basket by two screws. The water connector consists of the following components:

- 1 plastic connector
- 2 plastic blanking covers
- 3 screws



Adjust upper baskets horizontally only. The baskets are not designed for tilting (one side up, one side down).  
Adjusting the height alters the vertical clearance of the upper and lower baskets.

### To adjust the upper basket

- Remove the upper basket by pulling it out until a resistance is felt and lifting it off the runners.
- Unscrew the roller brackets and the water connector.

### Adjusting the height of the upper basket


#### Upper position

- Move the roller brackets on both sides to the lower position and screw them tight.
- Insert the connector into the opening at the bottom and click it into place by pushing it to the left as far as it will go. Screw the connector tight. Attach the blanking covers to the two openings above, click them into place and screw them tight also.


## en - Areas of application

### Middle position

- Move the roller brackets on both sides to the centre position and screw them tight.
- Insert the connector into the opening in the centre and click it into place by pushing it to the left as far as it will go. Screw the connector tight. Attach the blanking covers to the each of the openings above and below, click them into place and screw them tight also.

### Lower position

- Move the roller brackets on both sides to the upper position and screw them tight.
- Insert the connector into the opening at the top and click it into place by pushing it to the left as far as it will go. Screw the connector tight. Attach the blanking covers to the two openings below, click them into place and screw them tight also.


### Then check

- Replace the upper basket on the rails and push it in carefully to check that the water connector is positioned correctly.


## Screwing in nozzles

### Securing injector nozzles

- Screw the injector nozzles into the required positions.
- Use the 9 mm spanner (WAF 9) to tighten the injector nozzles.

## Test point for measuring water pressure

The connection at the top of the right-hand injector manifold is used as the test point for measuring water pressure. This connection point can be used for testing pressure within the framework of performance and validation testing in accordance with EN ISO 15883.

- To measure the water pressure, replace the existing fixture with a Luer Lock adapter, e.g. E 447.

### Preparing the load

For validation purposes please make sure that you follow the loading instructions given on the template.

⚠ Risk of injury caused by load.

There is a risk of injury when loading and unloading due to possible sharp edges, rims or pointed ends.

To minimise the risk of injury, loading should take place from the back forwards and unloading should take place in the reverse order.

⚠ Risk of infection due to insufficient disinfection.

Load items that become detached from the irrigation connector during reprocessing will not be sufficiently disinfected on the inside.

Any items that become detached from the irrigation connector during reprocessing must be reprocessed again.

### Before loading the machine and before starting a programme

- Check that fittings such as irrigation sleeves and injector nozzles are securely screwed into position.

⚠ Make sure that all screw connectors are fitted with injector nozzles, adapters, irrigation sleeves or blind stoppers to ensure that all fittings in use are supplied with sufficient standardised pressure.

Damaged fittings such as injector nozzles, adapters and irrigation sleeves must not be used.

Fittings not equipped with wash items do not have to be replaced blind stoppers.

- Check that spray arms are rotating freely.
- Check spray arm jets for blockages.
- Check that the basket is correctly connected to the water supply in the washer-disinfector.

### Reprocessing hollow instruments


Before inserting lumened instruments, please make sure that all hollow sections are free of blockages. Otherwise the suds solution will not be able to access all areas properly.

- Attach the lumened instruments to the nozzles.
- At the end of the programme, the insides of the instruments must be dried with sterile compressed air before being cleaned and sterilised as appropriate in accordance with the manufacturer's instructions. Please observe national regulations.

Where necessary, the nozzles supplied can be replaced with different fittings, e.g. holders for dental transmission instruments, irrigation sleeves, etc. Further examples are given on the following pages.

⚠ The wash water used for reprocessing narrow lumen instruments needs to be filtered. The injector manifolds can each be fitted with an A 800 tubular filter for this purpose. The tubular filter is supplied with its own operating instructions.

### Cleaning the A 800 tubular filter

The tubular filter must be cleaned **at least once a week**.

We recommend cleaning the tubular filter after every 10 programme sequences.

- Follow the cleaning instructions in the operating instructions for the tubular filter.

### A 803 holder for transmission instruments

A 803 holders must be retro-fitted for reprocessing dental transmission instruments. The holders are suitable for almost all standard transmission instruments.

## Installation


- Remove a nozzle from the injector manifold and replace it with an A 803 holder.

When the A 803 holders are used, the wash water must be filtered. The injector manifold must be fitted with an A 800 tubular filter. The tubular filter is supplied with its own operating instructions. Follow in particular the instructions on how to perform regular cleaning on the tubular filter.

## en - Areas of application

### Reprocessing transmission instruments

- Check that internal sections of transmission instruments are free of blockages.

- Insert the transmission instrument into the holder. The instrument is then held in position by the spring clamp.
- After reprocessing pull the transmission instrument out of its holder. The spring clamp will automatically spring back to its starting position.
- After reprocessing, the insides of the instruments must be dried with sterile compressed air before being cleaned and sterilised as appropriate in accordance with the manufacturer's instructions. Please observe national regulations.

Before using transmission instruments again following reprocessing, a function check must be carried out, e.g. by spraying into a basin, to ensure they are clear.

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
AUF 1 holders are designed for reprocessing dental transmission instruments.

### Optional accessories

- A 801, threaded adapter
- ADS 1, silicone adapter, white, for holding transmission instruments with a  $\varnothing$  of approx. 20 mm
- ADS 2, silicone adapter, green, for holding transmission instruments with a  $\varnothing$  of approx. 16 mm
- ADS 3, silicone adapter, red, for holding transmission instruments with a  $\varnothing$  of approx. 22 mm
- FP 20, 20 filter plates, porosity 2,  $\varnothing$  30 mm

**Installation**

An A 801 adapter is required.

- Unscrew the upper part ② from the AUF 1 holder.
- Fit a silicone adapter (ADS) ① in its place which has the same diameter as the transmission instrument being reprocessed.

⚠ Fitting a transmission instrument into a silicone adapter that is too big can cause damage to the filter plate.

- Remove an injector nozzle from the injector manifold and replace it with the A 801 adapter ⑤.
- Screw the lower part of the holder ④ into the adapter.
- Insert a filter plate ③ into the lower part.

⚠ Do not insert a filter plate into the holder if the injector manifold has been fitted with an A 800 tubular filter.

- Screw the upper part into the lower part.

**Reprocessing transmission instruments**

Before loading check that the silicone adapters are present and correctly fitted.

- Wipe away any deposits from the area around the silicone adapters that may be present depending on local water conditions using a damp cloth.
- Check that internal sections of transmission instruments are free of blockages.
- Fit the instruments into their holders using the respective silicone adapters.
- After reprocessing, the insides of the instruments must be dried with sterile compressed air before being cleaned and sterilised as appropriate in accordance with the manufacturer's instructions. Please observe national regulations.


Before using transmission instruments again following reprocessing, a function check must be carried out, e.g. by spraying into a basin, to ensure they are clear.

## en - Areas of application

### Replacing filter plates

⚠ Filter plates must be replaced **approx. every 2 weeks** or after **20 reprocessing cycles** have been completed.

- Unscrew the holder completely.


- Dismantle the holder and dispose of the filter plate ③.
- Rinse out the lower part ④ thoroughly.
- Insert a new filter plate.


⚠ Damaged or used filter plates must not be inserted under any circumstances.

- Screw the upper and lower parts of the holder back together and fit the holder back onto the injector manifold.

### A 814 holder for the tips of air scalers and piezo scalers


A 814 holder

A 814 holder with A 813 adapter

The A 814 holder enables the tips of air scalers and piezo scalers to be reprocessed. To increase capacity, the holder can be combined with the A 813 adapter. The holder and the adapter each come with their own operating instructions.

## Fitting E 499 for suction and rinsing cannulae

The E 499 is for reprocessing rinsing and suction cannulae, including suction breakers. It comes with its own instructions.

## Irrigation sleeves

Irrigation sleeves are suitable for reprocessing hollow instruments and instruments with the following external diameters:

- E 336, irrigation sleeve, length 121 mm, Ø 11 mm
- E 442, irrigation sleeves for instruments Ø 4 - 8 mm, 121 mm long, Ø 11 mm, cap with aperture Ø 6 mm
- E 443, irrigation sleeves for instruments Ø 8 - 8.5 mm, 121 mm long, Ø 11 mm, cap with aperture Ø 10 mm

The aperture can be varied using optional replacement caps:

- E 445, 12 caps for irrigation sleeves, aperture Ø 6 mm
- E 446, 12 caps for irrigation sleeves, aperture Ø 10 mm

The diameter of the instruments and the irrigation sleeve apertures or caps must be compatible with each other, e.g. caps for irrigation sleeves with an aperture diameter of Ø 6 mm for dental surgical suction cannulae.

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