EN	DE	FR	ES	PT	IT
EL	NL	DA	FI	NO	SV

VOCO Alfacomp LC

EN Instructions for use MD EU Medical device In accordance with DIN EN ISO 4049

Product description:

Alfacomp LC is a light-curing, highly radiopaque hybrid composite, coloured grey to contrast more strongly with the natural dental hard tissue. Its physical stability and high abrasion resistance were developed to meet the special requirements of restorations bearing masticatory loads. **Alfacomp LC** contains 86 % w/w inorganic fillers (equivalent to 70.1 % by volume) in a methacrylate matrix (BisGMA, UDMA, TEGDMA) and is cured by blue halogen or LED light. It is available in syringes and Caps for direct application.

Shade:

Alfacomp LC is grey coloured.

Indications:

- Class I and II restorative treatment in the posterior region
- Core build-up

Contraindications:

Alfacomp LC contains BisGMA, UDMA, BHT, TEGDMA. **Alfacomp LC** should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group:

Alfacomp LC is suitable for use in all patients without any age or gender restrictions.

Performance features:

The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User:

Alfacomp LC should only be used by a professionally trained dental practitioner.

Use:

Preparation:

Clean the teeth being treated. If applicable, mark occlusal contact points. Bring material to room temperature before use.

Moisture control:

Ensure an adequate moisture control. The use of a rubber dam is recommended.

Cavity preparation:

The cavity should be prepared according to the rules for adhesive restorative treatment; preparation should be minimally invasive to protect healthy dental hard tissue. Then clean and dry the cavity.

Cavity lining:

Suitable protection for the pulp should be applied when working close to it and should be covered with a layer of stable cement if necessary.

Placing the matrix:

For cavities in the proximal region apply a matrix. It helps to use translucent matrices, wedged in the proximal region. Keeping separation to the minimum makes it easier to shape the proximal contacts and attach a matrix.

Bonding material:

Alfacomp LC is applied with a dentine enamel bond, using the adhesive technique. Any light-curing bonding materials can be used. Observe the relevant instructions for use during preparation and processing.

Application:

Apply and cure restorations exceeding 2 mm incrementally in layers.

Light-curing:

Conventional polymerisation devices are suitable for light-curing the material. Curing time is **10 seconds per 2 mm** increment when an LED/halogen lamp with a light output of at least 500 mW/cm² is used.

The light aperture of the light guide should be brought as close as possible to the restoration. Make sure that every area of the restoration is sufficiently exposed to the light. You are advised to expose the accessible proximal boxes to the light orally or vestibularly. If you are using metal matrices, remove them beforehand. Inadequate curing can cause discolouration and sensitivity.

Finishing:

You can finish and polish restorations using cooling, immediately after removing the shaping aids (e. g. fine or extra-fine diamond burs, abrasive discs). Finally, the edge of the restoration, or better still the whole tooth, should be fluoridated.

Warnings, precautionary measures:

– Phenolic substances, especially products containing eugenol and thymol, lead to curing failures in restorative composites. Consequently, the use of zinc oxide/eugenol cements or other materials containing eugenol should be avoided in combination with restorative composites.

– Each **Alfacomp LC** Caps unit should be used for one patient only and should be applied exerting uniform but not excessive pressure.

– Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Constituents (in descending order according to content):

Barium aluminium borosilicate glass, BisGMA, UDMA, silicon dioxide, barium sulfate, TEGDMA, initiators, stabilisers, pigments

Storage instructions and application method: Storage at 4°C-23°C. Seal the syringes immediately after removing the material in order to avoid exposure to light and subsequent polymerisation. Do not use after the expiry date.

Disposal:

Dispose of the product in accordance with local regulations.

Reporting obligation:

Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of **Alfacomp LC** must be reported to VOCO GmbH and the responsible authority.

Note:

The Summary of Safety and Clinical Performance of **Alfacomp LC** can be found in the European database on medical devices (EUDAMED – https://ec.europa.eu/tools/eudamed). Detailed information can also be found at www.voco.dental.

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