Alphalink

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Directions for use

Alphalink is used to adhere metal components directly on the working model, ensuring an absolutely passive fit. The material can be applied on a model at the lab or directly in the patient's mouth.

A high proportion of inorganic fillers guarantees a high final hardness producing the stability required for laboratory preparation.

Indications

- Fixing telescopic or conical crowns to frameworks.
- Fixing pre-fabricated attachments to frameworks.

Constructing the joints

Design retentive and stable bonding surfaces. When placing the secondary component onto the primary component, make sure that there is an even, thin gap between the two components. This gap ensures a passive fit and it is filled with Alphalink. When using prefabricated retaining elements, create undercuts in the joint area. Observe the manufacturer's instructions.

Preparing the joints

After ensuring that both the primary and secondary components fit passively, create undercuts in them and sandblast with aluminium oxide (max. 125 µm).

After steam cleaning, condition the metal surfaces with Sebond Smart metal bonding system.

Dosage

Alphalink is available in the user-friendly dual-compartment cartridge as well as in two individual syringes.

When using the **automix syringe**, please proceed as follows:

- Place a disposable mixing tip onto the dual-compartment cartridge. Discard the first bit of mixed material discharged from the tip.
- Press Alphalink from the automix tip and apply the desired amount directly onto the restoration. Spread evenly.
- As the material will cure inside the used mixing tip, leave the tip on the cartridge as a closure until the next use.

- By using the automix cartridge with the mixing tip, Alphalink is always mixed in the ideal ratio of 1 : 1.

When using **individual syringes**, please proceed as follows:

Two equal lengths of components A and B are dispensed from both syringes by pressing the plunger and are mixed homogeneously.

Working time after starting tomix:

Approx. 3 min.

Setting commences after approx. 4 1/2 min.

Setting completed after approx. 8 min.

These times apply to an ambient temperature of 22 °C (71 °F).

When working inside the patient's mouth:

Working time after starting to mix: 45 sec.

Setting commences after approx. 1 min.

Setting completed after approx. 1.5 min.

Higher temperatures accelerate the setting reaction and lower temperatures retard it.

Application

Alphalink applied to matrix and patrix. The metal components are repositioned on the model and fixed into place.

Ensure that the areas to be bonded are completely filled with Alphalink and that there are not any air bubbles or contamination.

Application directly in the patient's mouth

Apply Alphalink with a spatula to the surfaces to be joined. Set and fix the metal elements to a passive fit.

Important

The sections being joined must not be moved until the Alphalink has set completely.

Please note

- Alphalink is dimensionally stable up to 120 °C (248 °F).
- Alphalink is easy to grind.

Storage

Refrigerate at 3-9 °C (37-48 °F). Do not use after the expiry date.

Side-effects

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us.

Contra-indications/interactions

If a patient has known allergies against or hypersensitivities towards a component of this product, we recommend not to use it or

to do so only under strict medical supervision. In such cases, we will supply the composition of our medical device upon request. The dentist should consider known interactions and crossreactions of the product with other materials already in the patient’s mouth before using the product.

Note

Please supply the dentist with the above information, if this medical device is used to produce a special model. Please also consider the safety data sheets.

Avoid contact of mucous membranes and skin with uncured material.

Troubleshooting

Working time too short	– High ambient temperatures (e. g. summer heat or central heating) decrease the working time.
Working time too long	– If the material is used directly after taking it from the refridgerator, curing will be delayed.
Material does not cure sufficiently	– Substances containing eugenol or oil of wintergreen prevent curing. – Material components are insufficiently mixed; only use original mixing tip.
Material is cured	– After use, the cap was replaced instead of leaving the mixing tip on the cartridge (contamination of component A + B).
Construction does not fit	– Construction was moved before the material had cured. – Material is already too firm: Please join fewer parts at a time.

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