# betasil®

## V A R I O

**müller-omicron**  
D E N T A L

- DE **Gebrauchsanweisung**
- EN **Instructions for use**
- FR **Mode d'emploi**
- ES **Instrucciones de uso**
- IT **Istruzioni per l'uso**
- NL **Gebruiksaanwijzing**
- RO **Instrucțiuni de utilizare**

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### Elastomeric precision impression material on vinyl polysiloxane basis, addition-curing.

Thank you for the confidence you have shown in our company. You have made the right decision with the betasil® VARIO impression system. Please always read these instructions for use through carefully before using the product! betasil® VARIO products are only intended for use in dental applications by correspondingly qualified professionals. The medical device must be employed responsibly and with due caution.

### Product description:

betasil® VARIO is an impression material which satisfies the highest demands. The putty-wash materials offer very high thixotropism, outstanding resilience and hydrophilic properties. All betasil® VARIO materials have a thermo-activating formulation which offers the user a variable processing time and thus the greatest processing comfort. The time in mouth is patient-friendly as it is always equally short. The system is rounded off by the highly contrasting colour system and the pleasant mild mint flavour.

### Indication:

Crowns and bridges, partial crowns, inlays, onlays, implants, transfer impressions, functional impressions.

### Impression technique:

- betasil® VARIO PUTTY / PUTTY SOFT: Preliminary impression material for the putty-wash technique (2-step), situation impression.
- betasil® VARIO PUTTY SOFT / HEAVY: Basic material for the double mixing technique (1-step); situation impression.
- betasil® VARIO MEDIUM / LIGHT / LIGHT fast set; wash material for putty-wash technique and double mixing technique.
- betasil® VARIO IMPLANT: Precision silicone for monophasic impressions (can be used as tray- and / or syringe material).

### Dosage / processing (see also illustration):

- betasil® VARIO putty materials for manual mixing applications  
Remove the required amount of betasil® VARIO putty with the respective measuring spoon in a mass ratio of 1:1 of the catalyst and base. Combine the two components, draw the putty out to form a small strand, then fold it together again and knead it with your fingertips. Repeat this process a few times until you achieve a homogeneous colour, but always complying with the minimum mixing time of 30 seconds. **Close containers tightly immediately after use and do not interchange the lids and measuring spoons.**

### betasil® VARIO materials in a 5:1 full cartridge system

The 5:1 full cartridges are suitable for use in all common mixers for 5:1 cartridge systems. Remove the end cap from the 5:1 full cartridge by gently pushing the tab upwards. When doing so, hold the cartridge with the openings facing upward. Do not drop the cartridge, as it can break and be rendered unusable. Insert the 5:1 cartridge into the device (please follow the manufacturer's instructions on using the device). Then start the device and let material flow out until material is flowing evenly from both openings. Now place the mixing tip (single-use) over the openings and check that it is positioned correctly. Then push the fix cap as far as it will go over the mixing tip and turn it to the right until it locks in place. Start the mixer again and dispense the desired quantity of material into an impression tray or into a separate intraoral impression syringe (the latter is only possible with betasil® VARIO IMPLANT). The used mixing tip remains on the cartridge and serves as a cap. Prior to the next use, loosen and remove the fix cap and remove the used mixing tip, check the opening for a plug that forms in isolated cases, clean if necessary and proceed as usual.

- betasil® VARIO materials in 50 ml safety cartridges  
Insert the cartridge in the 1:1 mixing gun, remove the end cap and dispense material until material is flowing evenly from both openings. Now place the mixing tip (single-use) over the openings (ensure that the colour code of the end cap and mixing tip are identical) and ensure that it is positioned correctly. To close it, turn the tip 90° clockwise. If necessary, use an intraoral tip (single-use). Operate the mixing gun again and dispense material in the desired amount. The used mixing tip remains on the cartridge and serves as a cap. Prior to the next use, remove the used mixing tip, check the opening for a plug that forms in isolated cases, clean if necessary and proceed as usual.
- betasil® VARIO materials in tubes (not available in all countries):  
Place the required amount of betasil® VARIO cat and base from the respective tubes in the same strand length in mass ratio 1:1 on a mixing block or glass plate. Mix both components with a mixing spatula until a homogeneous colour tone is achieved, but always complying with the minimum mixing time of 30 seconds. Depending on the type of material and intended use, apply the mixed material in an impression syringe or tray. **Close the tubes immediately after removing the material; do not interchange the lids.**

### Impression procedure:

When using a suitable tray adhesive, please follow the manufacturer's instructions for use. After inserting the tray with the betasil® VARIO materials into the patient's mouth, press it into place briefly and then hold in place until fully cured. Always check the material intraorally to assess the degree of curing. When using the putty-wash technique (2-step), the first impression must be cleaned and dried in order to ensure cross-linking with the wash material.

### Disinfection:

The impressions can be disinfected with a 2% glutaraldehyde solution. After removing the impression from the mouth, rinse it under running water for 15 seconds. Then disinfect the impression by immersing it in the disinfectant thoroughly. Observing the disinfectant manufacturer's specifications is essential. We recommend the (aldehyde-free) products Dentoprint MD liquid and Dentoprint MD pur from Müller-Omicron GmbH & Co. KG.

Further information can be found at [www.dgzmk.de](http://www.dgzmk.de), and in J. Prosth. Dent. 1999 May, 81(5), 621.

### Electroplating:

The impressions can be electroplated with copper or silver in standard baths.

### Model fabrication:

The impression can be poured after just 60 mins. There are no other time restrictions. All class 3 and 4 gypsums which comply with the requirements in ISO 6873, e.g., Moldano, Kulzer (Type 3) and Fujirock EP, GC (Type 4), can be used. Undamaged impressions can be poured several times.

### Working and safety instructions:

- The products described above must only be used in accordance with their instructions for use. The user has sole liability for any other use which contravenes these instructions for use.
- Close containers and tubes tightly immediately after use and do not interchange the lids and measuring spoons.
- Products marked "single use" are intended for a single use. To avoid cross-contamination, disinfect cartridges and mixing guns before the next use, e.g., with Dentalaprid SD wipes.
- Solutions used prior to impression taking, such as mouth rinses or retraction solutions, may interfere with the setting reaction and must be removed thoroughly.
- Contact between the material and certain types of gloves, e.g., latex gloves, should be avoided as the catalyst can be damaged. We recommend nitril / vinyl or PE-based gloves.
- Do not leave any impression material residues in the oral cavity.

- Avoid contact with clothing, as cross-linked silicones are chemically resistant and leave stains which cannot be removed.
- Do not use betasil® VARIO products with polyether, polysulphides or condensation-curing silicones.
- To prevent possible incompatibilities, do not combine betasil® VARIO products with A-silicones from other manufacturers.
- Storage: Store in a dry place at between 5°C and 27°C (41°F and 80.6°F).
- Shelf life: See date printed on package / label. Do not use after expiry date.
- Processing temperature: Room temperature between 18°C and 25°C (64.4°F to 77°F).

- Allergic reactions:** It is not possible to rule out allergic reactions in hypersensitive people. If acute allergic reactions should occur during the treatment, treatment with the product should be stopped immediately and a doctor should be consulted.

Should a serious incident occur unexpectedly in connection with this medical device, it must be reported to the manufacturer and the competent authority in your Member State.

### Please also observe the warnings in the safety data sheet!

Information at [www.mueller-omicron.com](http://www.mueller-omicron.com)

Errors and omissions excepted.

### Contraindication:

The product must not be used in patients with known allergies to flavourings or other ingredients [see safety data sheet or ask the manufacturer]. This should be clarified before use!

### Warnings:

- Eye contact:** Avoid contact with eyes. In case of contact, rinse eyes immediately with plenty of clean water and consult an ophthalmologist without delay.
- Swallowing:** If swallowed, rinse out the mouth and drink plenty of water. Always contact a doctor.

### Technical data and product features

Technical specifications	betasil® VARIO PUTTY / PUTTY SOFT	betasil® VARIO PUTTY SOFT, 5:1	betasil® VARIO HEAVY	betasil® VARIO IMPLANT	betasil® VARIO MEDIUM	betasil® VARIO LIGHT / LIGHT fast set
Consistency / EN ISO 4823	Kneadable / Type 0	Kneadable / Type 0	Heavy-bodied / Type 1	Medium-bodied / Type 2	Medium-bodied / Type 2	Light-bodied / Type 3
Colour	Light Green / Blue	Blue	Violet-Blue	Blue	Lilac	Pink
Dosage (1:1) / mixing time	30 sec. Automatic	-	Automatic Automatic	Automatic Automatic	Automatic -	Automatic -
Working time incl. mixing time*	Up to 2 min.	Up to 2 min.	Up to 2 min. 15 sec.	Up to 2 min.	Up to 2 min.	Up to 2 min. / up to 1 min. (fast set)
Time in mouth	≥ 2 min.	≥ 2 min.	≥ 2 min. 15 sec.	≥ 2 min.	≥ 2 min.	≥ 2 min. / ≥ 1 min. (fast set)
Total setting time	4 min.	4 min.	4 min. 30 sec.	4 min.	4 min.	4 min. / 2 min. (fast set)
Deformation under pressure (approx.)	2.0%	2.5%	2.0%	2.5%	2.0%	3.0%
Recovery from deformation (approx.)	99.4%	99.5%	99.6%	99.7%	99.6%	99.6%
Linear dimensional change (approx.)	0.20%	0.20%	0.10%	0.20%	0.13%	0.16%
Flavour	Mint	Mint	Neutral	Mint	Mint	Mint

\* At 23°C / 73.4°F. The times specified in the table above are based on a relative humidity of 50 +/- 10%. As a general rule: higher temperatures accelerate setting and lower temperatures delay it.

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