
# BioSonic® UC150

Ultrasonic Cleaner

<b>EN</b>	Owner's Guide	3	<b>NO</b>	Eierveiledning	233
<b>DE</b>	Benutzerhandbuch	27	<b>FI</b>	Käyttöohje	257
<b>FR</b>	Manuel de l'utilisateur	53	<b>EL</b>	Εγχειρίδιο χρήσης	281
<b>IT</b>	Guida per l'utente	79	<b>CS</b>	Příručka vlastníka	309
<b>ES</b>	Guía de uso	105	<b>JA</b>	オーナーガイド	333
<b>PT</b>	Guia do proprietário	131	<b>ZH</b>	用戶指南	359
<b>NL</b>	Gebruiksaanwijzing	157			
<b>DA</b>	Brugsanvisning	183			
<b>SV</b>	Bruksanvisning	207			

1.	Overview	3
2.	Safety Information	4
3.	Explanation of Symbols	5
4.	Unpacking	6
5.	Installation and Setup	6
6.	Filling Instructions	9
7.	Degassing	9
8.	Operation Instructions	10
9.	Cleaning Methods	12
10.	Accessory Use	12
11.	Data Log	12
12.	Maintenance	13
13.	Troubleshooting	15
14.	Technical Details	21
15.	Disposal Considerations	22
16.	Warranty information	22
17.	Order Information Details	25
18.	Figures	380

**Date of issue:**  
09-2018

# BioSonic® UC150

## Ultrasonic Cleaner

### 1. Overview

---

#### **An Essential Part of Your Infection Control System**

Congratulations on your purchase of a BioSonic® Ultrasonic Cleaner - a high quality, extremely efficient and popular ultrasonic cleaning system in dentistry.

Developed specifically for the dental industry, BioSonic units feature a unique product design to optimize time and safety in the office or laboratory.

BioSonic is a complete ultrasonic cleaning system, with a full line of cleaning solutions and unit accessories which provide you with the flexibility to efficiently clean instruments, prostheses, burs and other appliances.

---

#### **What is Ultrasonic Cleaning?**

Ultrasonic Cleaning is a process created by high frequency soundwaves. The soundwaves, enhanced by specially formulated cleaning solutions, create high energy cavitation. During cavitation, millions of tiny bubbles form and then collapse or “implode,” releasing enormous amounts of energy and shockwaves which scour the surface of instruments, appliances and other devices placed in the cleaning solution. This powerful scouring action reaches into minute crevices, which manual brushing cannot reach. The combination of energy and specially formulated solutions make ultrasonic cleaning the most effective method for removing large and microscopic debris.

---

#### **Intended Use:**

The BioSonic® UC150 Ultrasonic Cleaner is intended for cleaning medical instruments by the emission of high frequency soundwaves.

## 2. Safety Information

**NOTE: Failure to comply with these Caution and Warning statements may result in injury to personnel or equipment damage.**


 <p><b>CAUTION!</b></p>	<p>Do not operate the unit without sufficient solution in the tank. Do not use water alone as the carrier bath. Water is not a satisfactory coupling agent unless a wetting agent is used. For best results use BioSonic General Purpose, BioSonic Germicidal Cleaner or BioSonic Enzymatic Ultrasonic Cleaning solutions.</p>
	<p>Do not place objects directly on bottom of the tank because they may block cavitation and cause buildup of heat in the transducer causing premature failure. Use only BioSonic UC150 accessory baskets, cassette racks or beakers to hold items.</p>
	<p>Do not restrict airflow. Provide adequate ventilation by maintaining at least 2 inches clearance at the back.</p>
	<p>Do not subject unit to shock or impact.</p>
	<p>Do not immerse the unit in water.</p>
	<p>Keep the display surface dry for proper functioning of device.</p>
	<p>This product is intended for use only by trained healthcare professionals.</p>
	<p>Use this product in accordance with the instructions in the owner's guide.</p>
	<p>Unspecified or improper use of this product may impair safety protection.</p>
	<p>Before cleaning items of considerable value ultrasonically, consult the item's manufacturer for cleaning precautions.</p>
	<p>The power off selection from the keypad does not turn off main power to the unit. Unplug the device from the power outlet to shut down the device.</p>
	<p>This equipment may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures such as reorienting or relocating the unit or shielding the location.</p>
	<p>Never drain tank by manually discarding solution as this may cause the unit to stop working properly.</p>
	<p>Use only the supplied power cord or equivalent rated for 105°C.</p>
	<p>Do not disassemble the unit. There are no user serviceable parts inside.</p>
	<p>Keep unit covered to reduce airborne contaminants.</p>
<p>Do not place fingers or hands in the ultrasonic tank while in operation.</p>	
<p>Ensure appropriate decontamination in case of spillage of hazardous material on or inside of the equipment through device openings.</p>	

Operation for long periods of time will cause the cleaning solution to get very hot. To avoid scalding, do not place hands in the solution or spill any onto you. Be careful when removing an item from the tank as it may also be hot.

**WARNING!**

### 3. Explanation of Symbols

	Caution: Hot Surface		Settings Menu
	Warning/ Caution: Refer to accompanying documents		Log
	Caution: Risk of Electric shock		Heater Settings
	Temperature Limits		Degas Menu
	Humidity Limits		Date and Time Settings
	Proper Disposal Required		Language and Measures Settings
	Eco Cleaning Mode [15 min.]		'MY' mode settings
	Regular Cleaning Mode [20 min.]		Information
	Strong Cleaning Mode [30 min.]		Screen ON/OFF
	'MY' Cleaning/Custom mode		Manufacturer

## 4. Unpacking

---

**Figure 1 BioSonic UC150 Parts and descriptions**

- Remove your BioSonic unit from the shipping container and inspect the unit for any possible damage which may have occurred during shipping.
- Check the serial number of the unit (found on either the underside or rear of the unit) with the number on the shipping carton. If they are not identical, contact your dealer.
- **Register your product warranty at: <http://warranty.coltene.com>**

## 5. Installation and Setup

---

### General Setup

- Select your desired location for the cleaner. It should be on a solid level surface and in close proximity to a sink or to waste lines to facilitate draining. The hose supplied with the unit will allow you to place the unit on the left side of the sink or connect directly to the waste line.
- Connect one end of the hose to the drain valve (**Figure 4\***) and place the other end in the sink. If desired, have your plumber connect the hose directly to the waste line with the appropriate fittings.

\*see last page for figures

- **CONNECT THE DISPLAY.** Attach the connector on the bottom of the display module to the connector on the top right-hand corner of the unit. Press the display module down carefully until it snaps into the right position.
- **PLUG IN THE UNIT.** Attach the BioSonic UC150 electrical cord at the lower back of the unit and plug into a three-wire grounded AC outlet (in accordance with local regulations) (Figure 9\*). It is recommended that the unit be plugged into a Ground Fault Circuit Interrupter (GFCI) outlet.

---

**Installation instructions for recessing unit:**

- Tools Required:
  - Electric Drill
  - Drill bit to drill a hole large enough to accept your reciprocating saw blade
  - Reciprocating saw and blade

**NOTE:** This Ultrasonic Cleaner must be installed by an experienced equipment manager or licensed contractor. Be sure to follow local building codes.

**NOTE:** Before you begin please read and understand the safety instructions provided by the manufacturer of the tools required for installation.

- Select the desired location for the cleaner. You will need at least 19 cm (7.5 inches) clear vertical space below the top of the counter to accommodate the tank enclosure. The countertop must be flat and level. Make sure that you are not selecting a location which will interfere with drawers or vertical partitions between cabinets.
- Tape the enclosed template down to the counter top in the desired location. Make sure that the "front of the unit" is closest to the user.
- Cut along the "cut line" for the counter top with a jigsaw or other reciprocating saw.
- Clean away dirt or sawdust around the opening and remove scrap wood.
- Attach the rubber gasket on the housing for better seating between bezel and counter top.
- Lower unit into opening. Be sure that the printed BioSonic logo on the housing faces the user. (*Figure 11\**)
- The hose supplied with the BioSonic UC150 will allow you to dispose of the waste in accordance with your local, state and/or federal regulations.

\*see last page for figures

- Opening the drain valve will allow you to empty the solution into a waste storage container for disposal in accordance with your local, state and/or federal regulations.


**WARNING!**

(Failure to comply may result in personnel injury and/or equipment damage.)

Power should not be supplied to the unit until all interconnections between the display and tank enclosure and any accessories are securely connected.


This unit must be properly grounded. Failure to properly ground unit may result in serious injury.

**Power up Settings:**

- **SELECT A LANGUAGE AND TEMPERATURE UNIT:** Select the language for your appropriate region from the drop down list. Press button to accept.

*Figure 2 Setup Screen: Language and Measures*

- **SET THE DATE AND TIME:** Set the date and time with the preferred format by using the drop down list and format buttons. Press to accept. User can choose one of the three date formats:

- DD.MM.YY
- MM/DD/YY
- YY-MM-DD

*Figure 3 Setup Screen: Date and Time*

- **SETUP COMPLETE.** The display will now show the main screen once setup is completed in the language of your choice. Your BioSonic UC150 is now setup and ready to use!

**NOTE:** If modifications need to be made to the setup of the unit, go to settings menu by pressing to change the setup options.

## 6. Filling Instructions

---

- Ensure that the drain handle (**Figure 5\***) is fully rotated counter-clockwise to the closed position. Fill the unit with warm water to the fill line indicator (**Figure 6\***). Ensure that the water level is filled up to the minimum line and does not exceed the maximum line while using accessories like beakers or baskets.
- Pour the required amount of ultrasonic cleaning solution (eg. BioSonic General Purpose UC30, BioSonic Enzymatic Ultrasonic Cleaner UC32) and water into the tank. Note that some solutions require the use of a beaker.
- Change the solution in the tank daily or sooner if it appears soiled or discolored. When beakers are used, change the solution in the beaker after each use.
- Drain unit by fully rotating the drain handle clockwise to the “Open” position and allow the solution to empty into sink or waste system. Rinse tank and fully rotate the drain handle counter-clockwise to the “Closed” position. Never drain tank by manually discarding solution as this may cause the unit to stop working properly.

**NOTE:** Never operate the unit without solution. **Operating a dry tank will void the warranty.** Do not place objects directly on bottom of tank. Always use a cassette basket or beakers.

## 7. Degassing

---

Degassing the solution eliminates large air bubbles that interfere with the cleaning process and is necessary each time new solution is added to the tank.

**NOTE:** Be careful not to introduce unwanted gases into the liquid. Do not pump/pour the solution too quickly into the bath, insert/remove the product being cleaned into the bath too rapidly, or move product while in the bath. Doing so will re-introduce gas into the solution and require a longer degas cycle.

- PRESS THE DEGAS icon on  the main screen.
- Select the type of solution to log solution used or degas only and press  to degas.
- The BioSonic UC150 will begin to cycle on and off during the degassing process. The timer will count down from five minutes.
- DEGAS COMPLETE! Once the timer reaches zero, the words “Degas Complete” will be displayed. Your solution is now ready for use. Press  to return to the main screen.

\*see last page for figures

## 8. Operation Instructions

---

### Mode Selection

- Four modes are available to choose from on the main screen - 3 Preset modes: ECO: for light load time; REGULAR: for medium load time; STRONG: for heavy load time and "MY" mode for customized operation. (Figure 12\*)
- **'MY' mode:** For custom mode settings, select  from settings menu. Set the desired temperature and cleaning time. Ensure that HEATER function is ON if heating is desired. Press  to accept.
- Select mode by pressing appropriate icon for your particular load. Ensure that the lid is closed before and during operation. Press  to start cleaning cycle. Press  to pause. The BioSonic UC150 will continue to clean as the timer counts down the remaining time from the user-selected cleaning time.
- Once the timer reaches zero, "Program Successful" will be displayed on the screen until any button is pressed. Press  to return to the main screen or to abort the process.

**NOTE:** The cleaning cycle can be stopped at any time by pressing start or pause. Press  to return to the main screen.

- The BioSonic Ultrasonic cleaner should always be covered while operating to trap aerosols in the unit. The BioSonic UC150 is equipped with a hinged cover for easy one hand opening. The cover also conveniently lifts off for cleaning in spaces with low cabinets above the unit.
- Do not place additional instruments into a cleaning load mid-cycle. The additional instruments may not be cleaned as effectively as those in the original load. Additional instruments may also impact the cleaning efficacy of the original load by covering soiled instruments and increasing the mass of the load.

---

### Heater Operation

The BioSonic UC150 model is equipped with a heater to accelerate heating of the solution in the tank. Selecting a temperature instructs the heater to activate during the next cleaning cycle and remain active until the temperature is reached.

**NOTE:** The heater feature is not designed to maintain the selected temperature, but rather activate at the start of the next cleaning cycle and deactivate once the set temperature is reached. **Thermostatic control is not a feature of the Biosonic UC150.**

\*see last page for figures

**CAUTION!**

The heating function of the BioSonic UC150 (if enabled) takes place until the set temperature is reached only. If the set temperature is reached before the end of the cleaning cycle, the heater will shut off, but the cycle will continue. Through prolonged use, the natural cavitation process will continue to heat the solution beyond the desired temperature, but at a decreased rate without the heater. When this happens, the unit is still functioning properly.

**Activate the Heat Feature**

In order to obtain the activation code, visit <http://biosonic.coltene.com> or email customer service at [info.us@coltene.com](mailto:info.us@coltene.com) or [info.ch@coltene.com](mailto:info.ch@coltene.com).

**NOTE:** The temperature in the tank is displayed in the top bar of the LCD screen. Pressing the  key to activate this feature will bring up a new screen on the display. The first time this feature is activated, the user must enter the activation code. Press  once the code is entered. The heating function will now be available for use.

**Set the Desired Temperature**

- Press the  button. Set the temperature between 25°C (75°F) and 60°C (140°F).
- Select "HEATER ON" to activate heating function. Press  to accept the temperature.
- A  symbol will be displayed on the top of main screen if heater is selected ON. The heater will now function during the selected cleaning cycle.
- Select "HEATER OFF" if heating is not desired.
- A blinking  symbol on the main screen is displayed while the heating is in progress during a cleaning cycle. Current temperature information will always be displayed at the top of the LCD display.

**NOTE:** The heater will not turn on until a cleaning cycle has been initiated. It is, therefore, recommended that a pre-cycle (with heater activated) be performed when a specific temperature is desired. The cycle can be run without cleaning load and the length of the cycle would depend on the desired temperature.

**NOTE:** The user can change the temperature unit from C to F using the settings menu.

## 9. Cleaning Methods

---

### In Tank

Basket Cleaning: Place parts to be cleaned into the basket or use cassettes. (*Figure 8\**) Then place the basket directly into the main tank containing BioSonic General Purpose, Germicidal, or Enzymatic Cleaner. Follow dosage instructions per cleaning solution label. After cleaning, rinse the basket and instruments thoroughly under running water.

### In Beaker

For specialized item cleaning (e.g. prostheses, burs): Place items into a beaker and add desired cleaning solution (e.g. BioSonic Cement Remover or BioSonic General Purpose) according to the product directions. Suspend the beaker in the tank by using the beaker holder and place the beaker positioning ring around the beaker so that the bottom 2/3rd is suspended in the main bath.

**NOTE:** Beaker holder is to be used in conjunction with the basket (*Figure 7\**). Use caution when removing basket if beakers are left in place.

## 10. Accessory Use

---

- CASSETTE BASKET: The BioSonic UC150 basket cleans large volumes of instruments and increases instrument processing efficiency.

Place items in an instrument cassette or directly at the bottom of the basket. Place basket in tank and clean for desired time. After cleaning, rinse the basket and instruments thoroughly under running water.

- BEAKER HOLDER: The beaker holder holds a 600 ml glass beaker to clean small parts, burs, and prostheses.

Place the Beaker Holder on the basket and insert the basket into the tank. Place positioning rings around the beaker so that the bottom 2/3 of the beaker is suspended in the main bath. Clean for desired time.

## 11. Data Log

---

- VIEW LOG: Select  from main menu (*Figure 12\**) to view cleaner data log.

- DATA DOWNLOAD: Insert a USB stick in the USB port provided on the left of the BioSonic UC150 device display. To download cleaner use data onto the USB device, press  key. Press  key to confirm data download. Monthly data log files will be available in .csv format on the USB file (*Figure 13\**).

\*see last page for figures

## 12. Maintenance

---

### Cleaning

All components should be cleaned each day by wiping with a soft cloth moistened with a commercially available mineral deposit remover followed by cleaning with isopropyl alcohol or a mild disinfectant. It is recommended to drain and dry the tank at the end of each work day.

#### CAUTION!

Do not use petroleum based solvents, iodophors or phenolic based products. Iodophors and phenolics can stain the surface of the unit. Clean up all liquid spills immediately. Rust deposits on tank can be cleaned/ removed with a commercially available stainless steel cleaner.

---

### Stainless Steel Cleaning

#### ✓ Allowable:

- Wash with a clean cloth or rag soaked with warm water and mild soap or liquid dish detergent.
- Always clean in the direction of the original polish lines.
- Always rinse well with clear water (2 or 3 times) after cleaning. Wipe dry completely.
- Specialized household stainless steel cleaner may also be used.

#### × Avoid:

- Do not let dirt accumulate.
- Steel or stainless steel wool or scrapers to remove stubborn debris.
- Harsh or abrasive cleaners.
- Chemicals harmful to stainless steel not limited to but including products that contain chloride, fluoride, iodide or bromide.
- Cleaners that contain bleach.

*In case of questions regarding compatibility of decontamination or cleaning agents please contact the manufacturer or your supplier.*

### Inspecting the Power Cord Periodically

Check the power cord for damage and see that the connectors at both ends of the cord are fully seated.

### Changing the Fuse

- Remove the power cord from the wall outlet first and then from the power cord receptacle on the rear of the unit.

- The fuse compartment is located directly above to the power cord receptacle. To access, take a small, flat-headed screwdriver to the tabs on the fuse holder. Apply slight pressure to release the fuse holder (**Figure 9\***).
- Gently remove the fuses from the fuse holder.
- Examine the fuses (**Figure 10\***). If blown, replace with a fuse that matches the unit's voltage and current. If fuse appears to be in good condition, return the unit to your dealer or Coltène/Whaledent. This is an indication that the unit is defective.
- Plug the power cord into the power cord receptacle on the rear of the unit. Then plug the power cord into the wall receptacle.

**NOTE:** The fuse is a protection device designed to limit the amount of current being drawn or supplied to the unit. If a fuse is blown, this can mean that there was a surge of current to the unit due to an external power surge such as lightning, or more likely due to an internal component failure of the cleaner.

## **Firmware Update Process**

### **Checking the Firmware Version**

- Press the "Setup" icon from the main menu, bottom left icon on the display
- Press the "Information" icon from the setup menu, far right icon on the display
- The unit serial number and firmware version will be displayed
- Press the "Check Mark" icon then the "Home" icon to return to the main menu

### **Firmware Update Process**

- Place the most current firmware update folder "coltene-update" on a blank USB stick
- With the BioSonic UC150 displaying the main menu, place the USB stick in the display's USB port
- Wait a minute and the software update screen will appear on the display noting the current firmware version on the BioSonic UC150 and the current version on the USB stick
- Press the "Update" icon to start the software update process
- Do not touch the unit or the USB stick during the update process until the unit reboots and the main menu is displayed
- The firmware update process is complete and the USB stick can be removed

\*see last page for figures

## 13. Troubleshooting

**NOTE: Consult your dealer or Coltène/Whaledent for service if needed.**

### Error Codes

- The BioSonic UC150 is designed to interrupt a cleaning cycle and display an error code if a problem is detected. The table below lists the BioSonic UC150 error messages and the recommended user actions if encountered.

**NOTE:** The below table is provided as a reference to help users further understanding the BioSonic UC150's error code.

<b>Customer Error Code Table</b>			
<b>Error</b>	<b>Error Name</b>	<b>Description</b>	<b>Recommended User Action</b>
00	No Error	Cycle completed	Displayed in the user log file
01	Aborted by user	Cycle was aborted by the user	This error is displayed in the user log file only to note that the cleaning cycle was not completed. The unit can continue to be used normally.
02	Power interrupted	Power to (or within) the unit was interrupted	If power was interrupted by external means (power outage, unit was unplugged, display was removed) press the "OK" icon on the display screen to go to the main menu and resume the next cleaning cycle.
03	Liquid	Liquid level is too low	Verify that the solution level in the tank is filled between the minimum and maximum fill lined stamped on the inside back of the tank. If solution level is within the fill lines, contact customer service.
04	File missing	Firmware files are missing	Contact customer service for firmware update.
05-09	N/A	N/A	N/A
10	Over temperature Board	Temperature on power PCB assembly above 60° C/140° F	Ensure room temperature and humidity are within the recommended operational environmental conditions and there is at least 2 inches of clearance on all sides of the unit for adequate ventilation. If problem persists, contact customer service.
11	Over temperature Board	Hardware has shutdown due to temperature on power PCB assembly above 65° C/150° F	Ensure room temperature and humidity are within the recommended operational environmental conditions and there is at least 2 inches of clearance on all sides of the unit for adequate ventilation. If problem persists, contact customer service.
12	Over temperature Tank	Tank solution temperature is too high	If the unit has been operated for extended periods of time (e.g. multiple hours without intermission) the bath temperature may rise beyond an acceptable level. Either change the solution or wait for it to cool. If the problem persists, contact customer service.

<b>Error</b>	<b>Error Name</b>	<b>Description</b>	<b>Recommended User Action</b>
13	Over temperature Tank	Tank solution temperature is too high	If the unit has been operated for extended periods of time (e.g. multiple hours without intermission) the bath temperature may rise beyond an acceptable level. Either change the solution or wait for it to cool. If the problem persists, contact customer service.
14	Tank temperature rising	Solution bath temperature is increasing too quickly	Ensure that the tank is filled with solution to within the fill lines. If problem persists, contact customer service.
15	Tank temperature rising	Solution bath temperature is increasing too slowly	Ensure that the tank is filled with solution to within the fill lines. If problem persists, contact customer service.
20	Overcurrent	Abnormal current level	Contact customer service.
21	Overcurrent	Abnormal current level	Contact customer service.
22	Overvoltage	Abnormal voltage level	Contact customer service.
23	Overvoltage	Abnormal voltage level	Contact customer service.
24	Power deviation	Power output is lower than expected	Verify that the input voltage to the unit is within the specifications noted in the instructions for use. Deviations in supply voltage may trigger this error. If the problem persists, contact customer service.
25	Power deviation	Power output is lower than expected	Verify that the input voltage to the unit is within the specifications noted in the instructions for use. Deviations in supply voltage may trigger this error. If the problem persists, contact customer service.
26	Power deviation	Power output is higher than expected	Verify that the input voltage to the unit is within the specifications noted in the instructions for use. Deviations in supply voltage may trigger this error. If the problem persists, contact customer service.
27	Power deviation	Power output is higher than expected	Verify that the input voltage to the unit is within the specifications noted in the instructions for use. Deviations in supply voltage may trigger this error. If the problem persists, contact customer service.
28 - 39	N/A	N/A	N/A
40	Piezo no connection	Unexpected error with the piezo driver	Contact customer service.
41	Only one piezo connected	Unexpected error with the piezo driver	Contact customer service.
42	Piezo short-circuit	Unexpected error with the piezo driver	Contact customer service.
43	Piezo no power	Unexpected error with the piezo driver	Contact customer service.
44	Main switch failure	Main relay switch error	Unplug the unit and contact customer service. Do no attempt to reactivate the unit.
45	Tank Temperature Sensor failure	Temperature sensor not responding	Contact customer service.
46	Voltage Measure failure	Unexpected drive signal measurement	Contact customer service.
47	Current Measure failure	Unexpected drive signal measurement	Contact customer service.
48	Phase Measure failure	Unexpected drive signal measurement	Contact customer service.
49	N/A	N/A	N/A
50-69	Cable connection error	Internal cable failure	Contact customer service.

- **Unit is dead and the LCD screen on the unit doesn't light up**

Check the outlet that the unit power cord is plugged into by plugging a lamp (or similar device) into the outlet to verify that the GFCI outlet or circuit breaker is active and didn't trip. If tripped, verify that the current requirements are not exceeded by any devices plugged into the outlet or total circuit. If the outlet is working, remove the power cord from the IEC connector of the unit to gain access to the unit fuse holder. Follow the instruction for changing the fuse in the maintenance section of this manual.

## **Frequently Asked Questions (FAQ)**

### **1. "How do I test my unit to verify that it's working?"**

In order to verify the performance of your BioSonic UC150, follow the instructions on the BioSonic UC150 Ultrasonic Cleaner Performance Test flyer. A copy of the procedure is included with each device and is also available at [www.coltene.com](http://www.coltene.com).

#### **NOTE:**

- In some areas where the water quality is considered "hard," distilled water may be required.
- To adequately test an ultrasonic cleaner's ability to cavitate, solution must be present in the tank prior to verifying performance.
- Various solutions may inhibit the functioning of a unit (powder or tablet). It may therefore be necessary to dissolve the powder/liquid solutions in warm water prior to being added to the tank or use warm water in the tank and then perform the performance test.

### **2. "How long should I clean my instruments?"**

Listed below are several parameters to consider when determining the duration of any cleaning cycle. Adjustments to one or more of these parameters may be needed to achieve the desired level of cleaning.

**NOTE:** Always consult the manufacturer or suppliers instructions for use before cleaning a device.

The below guidance cannot be used to replace the cleaning instructions provided by the manufacturer of the instruments; they are provided to help users get the most out of each cleaning cycle. If soil is still visible after the completion of a cleaning cycle, modification to one or more parameters may be required, and then additional cleaning cycles may be required until all soil is, at a minimum, visibly removed.

- Size of the load being cleaned

The more instruments placed into an ultrasonic bath, the less effective the cleaning cycle will be. It is not recommended to place instruments in the basket so that they overlap, as open surface areas will be easier to clean. Filling a basket with clear space between instruments should lead to a more effective cleaning cycle.

- Amount and type of soil to be removed

Some contaminants have a stronger bond to the items being cleaned and may need a longer cleaning cycle. Also, the amount of soil on items being cleaned has an effect on the time required to complete a cleaning cycle.

- Type of cleaning solution used

Not all kinds of cleaning solutions are the same in regards to cleaning effectiveness. Some cleaning solutions work well on some contaminants, but not as well on other contaminants.

- Type of water

Not all water is the same and may vary in level of contaminants and hardness which can affect the level of cavitation in the solution. This can increase or decrease the time of a cleaning cycle.

- Solution temperature

The temperature of the solution can have an impact on the effectiveness of the cleaning cycle. The formation of cavitation events increases with temperature and the corresponding reduction in the fluid surface tension. Also, some solutions are more effective at a higher temperature, while others are more effective at a lower temperature. Consult the instructions of the solution for any requirements or recommendations from the manufacturer.

- Amount of soil in the solution

Repeated use of the same solution bath will lead to decreased performance as the bath becomes saturated with soil. Continue reading below for guidance on frequency of solution change.

- How long instruments have been soiled prior to cleaning

Soil should not be allowed to dry on instruments prior to cleaning; they should be processed immediately following a procedure. If they cannot be cleaned immediately following a procedure, it is recommended they be placed in a soaking solution so the soil does not harden on the instruments, which would increase the time of the cleaning cycle. The BioSonic ST16 Dental Soaking Tray is an example of a device designed for this purpose.

### **3. “How often should I change the solution bath?”**

How often the solution bath is changed depends on the solution used, frequency of use, and size of cleaning loads. It is recommended to change your solution every 8 hours at a minimum, or whenever the solution becomes visibly soiled or dirty.

### **4. “I’ve activated the heater option, but the temperature doesn’t remain stable. Why doesn’t the temperature remain at a set point?”**

Please refer to the section of the manual that describes the heater operation in detail. The BioSonic UC150 heater option only provides additional energy to the solution to increase the temperature more quickly than with cavitation alone. It is not designed to maintain a set temperature or cool the solution to a lower temperature.

### **5. “I’ve activated the heater option, but the solution isn’t heating.”**

Please refer to the section of the manual that describes the heater operation in detail. The BioSonic UC150 heater is only activated when the heater is turned on through the temperature menu and a cleaning cycle has been activated. The heater will not turn on during a degas cycle or when the unit is idle.

### **6. “Why is there a high pitch sound coming from the Ultrasonic Cleaner?”**

When a unit is initially filled up with a fresh solution mix there are lots of bubbles in the solution. If the unit is turned on right after fresh solution is added to the tank, it may make a screeching sound. Additionally, when placing a basket with instruments into the solutions, air is introduced and dissolves in the solution, which may create a screeching sound. The noise is present in all BioSonic cleaners, but is a normal occurrence and is not an indication of any issues with your device, simply that degas may be required.

To degas, please refer to the section of this manual that describes degassing in detail.

### **7. “What is degassing and why do I need to do this?”**

Before cavitation can become effective in an ultrasonic cleaner the gas trapped in the liquid must be removed. If it is not removed, the cavitation bubbles will fill with this gas, cushioning the implosion force. To degas, refer to the section of this manual that describes degassing in detail.

### **8. “What solution should I be using?”**

The purpose of the solution is to break the bonds between the instruments and the soil. Water alone has no cleaning properties and has high surface tension, which inhibits the cavitation bubbles from developing. Adding cleaning solution reduces fluid surface tension and contains ingredients designed to optimize the ultrasonic cleaning process.

Liquid solution allows for instantaneous solution preparation where tablets and powder solutions must be fully dissolved prior to beginning the cleaning process. It is recommended that

tablets and dry solutions be dissolved in hot water prior to placing them in the tank to prevent them from settling on the bottom of the tank until the ultrasonic action dissolves them.

The BioSonic system offers a full line of cleaning solutions which offer the flexibility to efficiently clean and disinfect instruments, prostheses, and other items. It is recommended that BioSonic solutions be used with your BioSonic cleaner for consistent and effective performance.

Visit [www.coltene.com](http://www.coltene.com) or the order information details section of this manual for more information about BioSonic solutions.

### **9. “What can be cleaned using my Ultrasonic Cleaner?”**

Ultrasonic cleaning is effective for most instruments such as:

- Dental burs which are reprocessable
- NiTi endodontic files
- Jointed instruments, e.g. scissors
- Stainless steel syringes
- Serrated instruments
- Plastics, cemented glass syringes, mirrors etc. may not be successfully processed using ultrasonics.

**NOTE:** Always consult the manufacturer’s instructions regarding the cleaning of instruments/ devices to ensure suitability for ultrasonic cleaning!

### **Service**

When calling customer service please have the following information available:

- Office name and phone number
- Contact person in case the service department needs to talk to the end user of the equipment
- Contact email
- Doctor’s name
- Model number of the unit
- Serial number of the unit
- Dealer’s name
- Detail description of the problem including any error codes displayed
- Photos and short videos of the unit problem are very helpful
- Copy of your purchase invoice for warranty validation
- Solution type (liquid, dry, or tablet) and the brand name
- Any cassettes that are being used, and if so, the cleaning cycle time

Customer service may require the end user to troubleshoot the unit during the call to determine if the issue requires a service technician. Please be prepared and ready to access the unit.

If service is required, a Coltene Customer Service Representative will issue the customer a complaint number and return shipping instructions.

For service, contact your dealer or Coltene customer service.

## 14. Technical Details

### Environmental Conditions:

- Intended for indoor use.
- Maximum altitude 2000 m.
- Recommended Storage and transport temperature: 0°C (32°F)  45°C (113°F)
- Temperature during operation: 5°C (41°F)  40°C (104°F)
- Relative humidity: 50%@40°C (104°F)  80%@31°C (88°F)
- Main supply voltage fluctuations not to exceed  $\pm 10\%$  of the nominal voltage.
- IEC-60664-1 Pollution degree 2.

### Electrical Information:

- Power Supply Input: 100 V – 240 V AC / 50 ~ 60 Hz
- Transducers: 2
- Operating Frequency: 48 kHz - 50.5 kHz
- Fuse Specifications: 2 Fuses FST, 5 × 20 mm, Slow Blow, 250 VAC, 5A, IEC 60127-2 Sheet 3

Max Power @ Input Voltage			
Input Voltage	Max Power without Heater*	Max Power with Heater*	Max Current
100 V	160 W	425 W	4.3 A
115 V	160 W	510 W	4.4 A
230 V	160 W	510 W	2.2 A
240 V	160 W	540 W	2.3 A
*average over a sweep period			

- EMC Class A: Suitable for use in all establishments other than domestic and those directly connected to a low-voltage power supply network which supplies buildings used for domestic purposes.
- Class I electrical equipment per IEC 61140.

- Transient over voltages category II per IEC 61010.

**Mechanical information:**

<b>Tank Capacity</b>	6.6 L (1.75 US gal)
<b>Use Capacity</b>	5.7 L (1.5 US gal)
<b>Overall Dimensions</b>	L × B × H = 34.5 × 27.6 × 35.5 cm (13.6" × 10.8" × 14") (w/ display)
<b>Tank Dimensions</b>	28.7 × 16.0 × 17.8 cm ( 11.3" × 6.3" × 7")
<b>Above countertop surface</b>	34.5 × 27.6 × 17.2 cm (13.6" × 10.8" × 6.8")
<b>Below countertop surface</b>	32.9 × 26 × 18.3 cm (13" × 10.2" × 7.2")

## 15. Disposal Considerations

The BioSonic UC150 is an electric appliance requiring proper disposal according to local laws and regulations.

*This appliance is labeled in accordance with European directive 2012/19/EU concerning used electrical and electronic appliances (waste electrical and electronic equipment - WEEE). This guideline determines the framework for the return and recycling of used appliances as applicable throughout the EU. The symbol on the product or on the documents accompanying the product indicates that this appliance may not be treated as household waste. Instead it shall be handed over to the applicable collection point for the recycling of electrical and electronic equipment.*

## 16. Warranty information

The BioSonic UC150 has been developed specifically for use in dentistry and intended to be operated only by dental professionals in accordance with the instructions contained in this guide. However, notwithstanding anything contained herein, the user shall at all times be solely responsible for determining the suitability of the product for the intended purpose and the method of its use. Any guidance on application technology offered by or on behalf of the manufacturer or distributor, whether written, verbal or by demonstration, shall not relieve the dental professional of his/ her obligation to control the product and to make all professional judgments regarding its use.

This product has been carefully manufactured and has met stringent quality assurance requirements. Our products are manufactured from new parts or new and serviceable used parts. Regardless, our warranty terms apply. If, in normal use, it is found defective in material or workmanship within a period of two years from date of purchase, it will be repaired at Coltène/

Whaledent's option, or replaced. Tampering with any of the components, misuse, negligence, alteration, water damage, accident or lack of reasonable or proper maintenance and care will void this warranty. Failure to use and/or maintain the product in accordance with the instructions contained in this manual will void this warranty.

**THE FOREGOING IS IN LIEU OF ANY EXPRESSED OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION, ANY WARRANTIES OR MERCHANTABILITY OR FITNESS FOR USE.**

**COLTÈNE/WHALEDENT AG ASSUMES NO RESPONSIBILITY FOR ANY INCONVENIENCE, LOSS, INJURY OR DIRECT, INDIRECT OR CONSEQUENTIAL DAMAGE ARISING FROM THE POSSESSION OR USE OF THIS PRODUCT.**

Any claim for damage or breakage in transit should be made at once against the carrier. If factory service is required, be certain to properly pack your equipment and return prepaid and insured to the factory.

#### **CERTIFICATE OF LIMITED WARRANTY**

Coltène/Whaledent (C/W) warrants this product to be free from defects in workmanship and materials for a period of twenty-four (24) months from the date of original purchase. During such warranty period, C/W undertakes to repair or replace this product, at its option, if such product does not function properly under normal use and service, which malfunction is due solely to a defect in workmanship or material.

This warranty does not cover any of the following:

- Damage caused by accident, abuse, misuse, neglect, alteration, transit, failure to follow manufacturer's instructions, lack of reasonable and proper maintenance and care, or Acts of God.
- Normal wear and tear (i.e.: bulbs, fuses, etc.).
- Minor imperfections that do not materially affect the product's performance.
- Use of the product by anyone other than a licensed dental practitioner or qualified dental hygienist.

This warranty shall become null and void if any repair or servicing of the product is made or attempted by any person not authorized by C/W, or if any part not authorized by C/W is used in any repair or servicing.

**This warranty is in lieu of any and all other warranties relating to the product, expressed or implied, including, without limitation, implied warranties of merchantability and fitness for a particular purpose. In no event shall C/W be liable for direct, indirect, tort, exemplary, punitive, special, incidental, or consequential damages arising from the possession or use of the product, including but not limited to loss of use of the product, lost revenues or**

**profits, or the cost of substitute equipment. The limitations provided for herein shall inure to the benefit of C/W's suppliers and subcontractors.**

This Certificate of Limited Warranty is valid only if it has been completed with date of purchase and dealer's name.

If this product becomes defective during the warranty period, it must be delivered together with this Certificate of Limited Warranty and a written warranty claim to the C/W dealer from which it was originally purchased. All warranty claims must include a detailed written description of the claimed defect. The purchaser will be responsible for all costs and risks of transit, and we therefore suggest that the product be carefully packed and that it be shipped to the dealer insured for its full value.

Register your product warranty by either option 1 or 2:

**#1. Please return the following information to**

Coltène/Whaledent AG, Feldwiesenstrasse 20, 9450 Altstätten / Switzerland

Name and address of purchaser:

Distributor's name and address:

Phone:

Email:

Model Number:

Date of purchase:

Serial number:

**#2. Register your product warranty at: <http://warranty.coltene.com>**

## 17. Order Information Details

	REF No*	Description	Quantity
UC150 Units	60021987	UC150 Cleaner Unit USA Plug Type A	1
	60021988	UC150 Cleaner Unit Euro Plug Type F	1
	60021989	UC150 Cleaner Unit UK Plug Type G	1
	60021990	UC150 Cleaner Unit AUS Plug Type I	1
Accessory kit	60022042	Basket, Beakers (2), Covers(2), Positioning rings(2), Beaker holders(2)	1
Accessories	60022040	Cassette Basket (26.8 × 15.8 × 11.5 cm 10.4 × 6.2 × 4.4")	1
	60022041	Beaker Holder	1
	UC53	600 ml Beaker with Cover and Positioning Ring	1
	60022043	UC150 Replacement Drain Filters	3
	60022044	UC150 Replacement Drain Hose	1
	60022045	UC150 Replacement Lid	1
	60022046	UC150 Gasket for sealing	1
Ultrasonic Solutions	UC30	General Purpose Cleaning Solution Concentrate in MeterDose® Dispensing Bottle	473 ml/16oz
	UC31	BioSonic General Purpose Cleaner plus Super Rust Inhibitor	15 ml × 20
	UC32	Enzymatic Ultrasonic Cleaner Concentrate in MeterDose Dispensing Bottle	236 ml/8oz
	UC34	Plaster & Stone Remover Powder Concentrate in Jar	840g/1.85 lbs
	UC38	Germicidal Ultrasonic Cleaning Concentrate in MeterDose Dispensing Bottle	236 ml/8oz
	UC39	Cement Remover Concentrate in MeterDose Dispensing Bottle	473 ml/16oz
	UC42	Germicidal Ultrasonic Cleaning Concentrate in MeterDose Dispensing Bottle	473 ml/16oz
	60021960	UC40 Universal Ultrasonic Cleaner Concentrate, 1L dosed bottle*	1L
	60021961	UC40 Universal Ultrasonic Cleaner Concentrate, 5L Refill can*	5L

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#### Coltène/Whaledent Inc.