## Mode d’emploi Bio-Ink® Flow

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Color Indicator Paste for checking sore spots on dentures and partial dentures

## Instruction for Use Bio-Ink® Flow

### General description and purpose of the medical device:

Color Indicator Paste for checking sore spots on dentures and partial dentures

#### Medical Product - For professional use only

#### Indication

Bausch Intra Oral Ink = Bio-Ink® - Flow green is indicated for marking sore spots on the gum. Bio-Ink®-Flow is a highly viscous paste, containing green food dye.

It can be applied directly to the oral mucosa with the single use needles or the attached Magic-Brush™ applicator. The sore spot on the gum is gently marked. After replacing the denture, a green contact spot can be detected on the lower side of the denture showing the area to be adjusted.

Sore Spots also known as decubitus can occur frequently with new dentures or partial dentures. Also a degeneration of the alveolar ridge can cause sore spots as well as unbalanced occlusal forces on the denture. The adaption of the oral mucosa to the new denture often causes inflammatory lesions.

Because of its thixotropic properties the green paste can easily be applied with Magic-Brush™ applicator in spots, areas and lines. Also the Palatal Vibrating Line and the mucobuccal fold can be marked. If necessary the denture can be adjusted. Before application the areas to be marked on oral mucosa have to be dried.

#### Contraindication

A direct application onto highly inflamed areas or purulent lesions of the oral mucosa areas should be avoided.

#### Side Effects

Possible side effects could be allergic reactions.

#### Notes

- The color paste can be applied directly with the disposable cannulas (single use applicator needles) and spread with the disposable applicator (e.g. Magic-Brush™, ) to avoid cross-contamination.
- Alternatively, the paste can be applied to a disposable applicator (e.g. Magic-Brush™, Easy Brush, ...) and then applied to the surface to be tested.
- The mucous membranes, palate and mandible must be dried with clean pulp or cotton roll previously.

#### Application

Carefully dry the gum, the oral mucosa and the alveolar ridge with a clean cotton roll to avoid smearing of the color indicator paste on the gingiva.

There should be no saliva on the surface to be tested. The denture should be cleaned and dried as well. The color paste could also be applied with single use applicator needles onto the gingiva.

#### Effect / Diagnosis

The color marking of the sore spot will be transferred to the denture when the denture is seated into position. When the denture is removed the spot is visible on the surface of the denture. If necessary the denture can be adjusted by the dentist.

#### Cleaning

After diagnosis, the color can simply be removed with water or any detergents.

#### Reprocessing of multiple use components

- In principle, the following should be observed: „Requirements for hygiene in the processing of medical devices Recommendation of the Commission for hospital hygiene and infection prevention (KRINKO) at the Robert Koch Institute (RKI) and the Federal Institute for Drugs and Medical Devices (BfArM)“, Bundesgesundheitsblatt 2012 -55: 1244-1310.
- The instrument must be reprocessed immediately after each treatment (within a maximum of 2 h)
- Thoroughly clean the brush holder by hand with a small brush with firm bristles.
- Note: The cleansing should be done in a water bowl without further additives below the water surface in order to achieve adequate cleansing of the instrument parts whilst avoiding protein fixation as well as to protect the environment from contamination due to splash water.

- Intermediately rinse the brush holder with water (minimum drinking water quality).
- Place the parts in a standard cleaning and disinfecting bath. Examples: Becht Bechtol Futura, Dürr Dental ID 213 instruments Disinfection, Pluradent instrument bath Pluline, Schülke & Mayr Lysetol V
- Refer to „List of disinfecting agents and methods tested and approved by the Robert Koch Institute“ or disinfectant list VAH.
- Note: Formaldehyde containing cleaning and disinfecting agents may only be used after sufficient cleaning in order to avoid protein fixation.
- Note: The instructions for use of the manufacturer of the cleaning and disinfecting agent must be strictly observed. In particular, the concentrations to be used and the reaction times must be observed!
- Note: For machine cleaning, the manufacturer’s instructions for use must be followed carefully!
- Final rinse of the instrument parts with water (minimum drinking water quality, recommended: demineralized water with microbiological quality according to drinking water).
- Dry.
- Visual inspection for corrosion, damaged surfaces, chipping, mold damage and impurities. Damaged instruments must be discarded (For a limited restriction of reuse, see the chapter „Reusability“). If residual contamination still exists, the entire cleaning process must be repeated with all steps (Cleaning, intermediate rinsing, disinfection, final rinsing and drying).
- The instrument must be free of any residue and dry before further processing.
- Maintenance of the instruments is not required.
- The instrument is packaged and sealed in a sufficiently large single sterilization bag according to EN 868-2 ff ISO 11607 (suitable for steam sterilization). Observe the instructions of the manufacturers of the sterilization bags and the sealing machines and the current normative requirements.
- Sterilization must be carried out in a validated process with moist heat in an autoclave according to DIN EN 13060 type B or DIN EN 285 resp. ANSI AAMI ST79. Observe the manufacturers operating instructions of the autoclave.
- Sterilize the instrument with humid heat (saturated water vapour) using a pre-vacuum procedure for 5 minutes at 134°C.
- After sterilization, the instrument must be stored dry and dust-free in the closed sterilization package.
- The recommended storage life for sterile medical devices is described in standard DIN 58953-8 and depends on external influences and effects during storage, transport and handling.
- **Reusability**
- Frequent reprocessing restricts the life of the plastic instrument. The instrument must be disposed of in case of cracks in the plastic.
- It is not permitted to reuse single use items (Magic Brush™ and applicator needles). Safe use cannot be guaranteed if these products are used again, since there is an infection risk and / or the safety of the products is not liable.
- The use of damaged or soiled instruments is the responsibility of the user.
- In case of non-compliance, all liability is excluded.

#### Disposal instructions

Bio-Ink® - Flow consists of a green food color and natural waxes, which are approved for food, cosmetics and pharmaceutical applications. All components can be disposed of with the usual practice waste.

**Note: Serious incidents occurring in connection with the product must be reported to the manufacturer and the competent authority of the Member State.**