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**Manufacturer Information on Reprocessing of Medical Devices according to
DIN EN ISO 17664:2004**

Manufacturer: HAGER & WERKEN GmbH & Co. KG
P.O.B. 10 06 54, 47006 Duisburg, Germany
T +49 (203) 99269-0, F +49 (203) 299283
www.hagerwerken.de

Product: Cotton-Clip II
REF 605 351
valid as of 02/2010 (Rev: 1)

Warning notices:	no special notices
-------------------------	--------------------

Limitation of reprocessing: Frequent reprocessing has only got a small impact on this instrument. The instrument shelf-life is mainly determined by wear and damage by use.

Place of Handling: Clean surface with a disposable cloth or paper towel.

Storage and Transport: Storage and Transport.

Preparations for Cleaning: n/a

Cleaning: Only use suitable cleansers according to the manufacturer instructions. After cleaning, remove cleanser thoroughly under running water. Do not exceed 93 °C for mechanical cleaning.

Disinfection: Use suitable disinfecting solution according to the manufacturer instructions. Remove disinfecting solution thoroughly under running water.

Maintenance: No special requirements.

Packaging: Standardized packing material for sterilisation can be used (recommendation).

Sterilization: Vapor sterilization at 134 °C for 5 minutes or vapor sterilization at 121 °C for 15 minutes.

Control/Functional Check: Sight check on damages, wear, deformation.

Storage: No special requirements.

The a.m. instructions have been validated as SUITABLE for preparation of a medical device and its re-use by the medical device manufacturer. It is the reprocessor's responsibility that the actual reprocessing achieves the required results with the equipment and materials applied and personnel involved in the reprocessing site. Normally, validation and routine control are necessary. Furthermore, the reprocessor should evaluate any deviation from the provided instructions regarding efficiency and possible adverse consequences.

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Fabricante: HAGER & WERKEN GmbH & Co. KG