 **dentona**[®]
NEW WAY DENTAL TECHNOLOGY

new
way
dental
technology

CE

Gebrauchsinformationen Working Instructions

dentona[®]-tray LC

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Lagerung

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dentona®-tray LC | Light curing custom tray

dentona®-tray LC is an easy moulding and strong light cure individual tray material suitable for:

- individual trays (for jaws with full and partial dentition)
- function trays (for toothless jaws)
- base plate for wax pattern
- individual and functional bite plates and occlusal registrations

The material is ready for use and easy to handle. It can be used within sufficient time to mould, under day and room light. dentona®-tray LC can be cured in all typical polymerization devices: halogen light (about 420 - 480 nm) or UV light (wavelength of 350 - 400 nm).

The working characteristics during shaping and after curing are excellent. dentona®-tray LC shows a smooth and non sticky surface.

Application

Before the working process starts, store the trays under room temperature. Mark the line of the tray with a pencil. If necessary, put space-stoppers onto the tray. In case no space-stoppers are needed, isolate the plaster model. Tear off the fleece completely, place the tray onto the model and bring the tray in shape. Cut off the edges, as drawn before. The grip can be formed from surplus material.

Curing

Polymerization takes place in a suitable light oven with halogen light with a wavelength of 420 - 480 nm or UV-light with a wavelength of 350 - 400 nm. Curing time depends on the type of the used light oven. Curing takes place on the plaster model inside the light curing-device, then take the tray off the model and cure the reverse side of the model:

Lumos from dentona®:

polymerization time 5 - 6 minutes each side

Device with minimum 4 UVA-tubes:

polymerization time 10 minutes each side

Longer curing times improve the physical properties of the material (higher strength). The cured tray is now ready for use and can be manipulated. The inhibition layer, which occurs by the oxygen inhibition, can easily be wiped off with an alcoholic-moisturized napkin.

Storage

Store at room temperature. Protect it from light.