Detaseal® function

DE

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Detaseal[®] function

GB
Instruction for use

Precision impression material, silicone based, addition curing, plastic phase, heavy bodied consistency

Fig. 1


Fig. 2

Fig. 3

1. Preparation of the tray

Depending on use, select a suitable impression tray (custom or prefabricated tray of metal or hard resin). **trayloc[®] A** is recommended for secure adhesion of the material to the impression tray. Coat the tray with a thin layer of **trayloc[®] A** and allow to dry for 3 min. (Fig. 1).

2. Mixing and Dosing

Extrude **Detaseal[®] function** from the proportioning syringes. Maintain a 1:1 proportion by using strands of base and catalyst of the same length. Knead between the fingers until a homogeneous mixture is obtained (Fig. 2). The tray can now be filled. Application for shaping the individual margin is facilitated by using the enclosed application syringe (Fig. 3). The tip of the application syringe can be shortened to increase the thickness of the strand.

3. Disinfection

After removing from the mouth, rinse the impression under running lukewarm water. If desired, impression can then be disinfected by immersing into 2% glutaraldehyde or other usual solutions for 15 minutes.

4. Model casting

The impression should not be cast for 30 min. after removing from the mouth. The dimensional stability is guaranteed for 14 days. Recommended model materials include class III and IV dental gypsums.

5. Cleaning of trays


Remove set material mechanically with a blunt instrument. **trayloc[®] A** is dissolved by placing in acetone or isopropanol. Use only in well ventilated rooms. Then clean and disinfect tray as usual.

Important working hints

- Do not process with condensation curing silicones.
- Cured impression materials are chemically resistant – avoid getting them on clothing.
- Do not leave any impression material residues in the mouth.
- Latex gloves and latex-contaminated surfaces can influence the setting of **Detaseal[®] function**.
- For a high resolution correction, an addition curing correction material, e.g. **Detaseal[®] hydroflow lite/Xlite** is recommended with the static impression procedure.
- Do not interchange the caps of base and catalyst.
- Please follow the instructions of the safety data sheet!

Further information:

Silicone based materials are proven a million times. On condition of a proper application, undesired effects are not to be expected. However, reactions of the immune system like allergies, irritations, cannot be absolutely excluded. In case of doubt, we recommend to make an allergy test before the application of the material.

Indications for use:

- Special impression material for functional impression taking with a prolonged setting phase
- Functional border moulding

Technical Data:

DIN EN ISO 4823 - Type 1

- **Mixed volume:** 2 x 80 ml (dispenser)
- **Mixing ratio:** 1:1
- **Colour code:** Base: red Catalyst: yellow
- **Mixing time:** approx. 30 sec.
- **Working time:** approx. 2 min. 30 sec.*
- **Setting time:** approx. 6 min.*
- **Time in the mouth:** approx. 3 min. 30 sec., incl. 2 min. 30 sec. plastic phase
- **Strain in compression:** approx. 9.5 %
- **Recovery from deformation:** approx. 99.5 %
- **Linear dimensional change:** approx. 0.2 %
- **Application:** at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity
- **Storage:**

* from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity. Increased temperature accelerate, decreased temperatures retard a. m. times.

Ordering information:

Detaseal[®] function
Standard packing **02521**
2 x 80 ml dispenser
base + catalyst
1 moulding syringe

trayloc[®] A
bottle with brush **03098**
17 ml
Vaporizer of 15 ml **03392**

application syringes **02802**
10 pcs.

DETAX

GmbH & Co. KG
Carl-Zeiss-Str. 4 · 76275 Ettlingen/Germany
Telefon 07243/510-0 · Fax 07243/510-100
www.detax.de · post@detax.de

Caution: Federal U.S. law restricts this device to sale by or on the order of a dentist (or trained specialist personell).

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