Detaseal® hydroflow mono 5:1

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Detaseal® hydroflow mono 5:1

GB
Instructions for use

VPS Precision silicone impression material, medium consistency – medium bodied – normal setting

Fig. 1

1. Preparation of the impression tray

Select a tray type suited for desired application (standard-size or custom, in metal or hard plastics). To avoid teeth contact with the tray, place stops from wax or light-curing resin into the tray at non-relevant areas of the impression. For optimal adhesion of impression material, we recommend applying a thin layer of **trayloc® A** to tray (Fig. 1). Let dry for approx. 3 minutes.


Fig. 2

2. Mixing and dispensing

Disengage and remove the tamper-proof seal on the cartridge by raising the locking lever (Fig. 2). Insert the cartridge into the dispenser. Start the device and dispense material until material flows uniformly. Discard the first inhomogeneous material, then attach the dynamic mixer (yellow). Fix the mixer by the applying and turning the yellow rotating ring clockwise (Fig. 3).


Fig. 3

Then fill the tray with the desired quantity of **Detaseal® hydroflow mono**, holding the tray at an angle and pressing against it slightly. While dispensing, leave the mixer in the material (Fig. 4). Leave the filled mixer on the cartridge holder as a closure. Before the next application remove the used mixer by undoing the rotating ring and check the cartridge openings for blockages. Attach a new mixer and proceed in the normal way. When starting up, discard the initial material emerging.


Fig. 4

3. Procedure

At the beginning introduce a sufficient quantity of **Detaseal® hydroflow mono** to the impression tray. Then syringe around the relevant areas with **Detaseal® hydroflow lite** or **Detaseal® hydroflow Xlite** intraorally. For this purpose intraoral tips or application syringes (Figs. 5 and 6) can be used. Position the loaded tray in the patient's mouth, applying slight pressure. Fix the impression without pressure until the material has set completely (Fig. 7).


Fig. 5

4. Disinfecting

After removing tray from the mouth, rinse impression under running, lukewarm water. If desired, impression can then be disinfected by immersing into 2% glutaraldehyde or other usual solutions for 15 minutes.


Fig. 6

5. Pouring the model

We recommend waiting at least 30 minutes before pouring the model. The impression will retain its dimensional stability for up to 14 days. All class III and IV dental plasters and standard modelling resins can be used.

Fig. 7

6. Cleaning the tray

After setting, material can be removed with a blunt instrument. Remaining film of **trayloc® A** is dissolved by immersion in acetone or isopropanol (be sure to use only in well-ventilated rooms). Tray can then be cleaned and disinfected in the usual manner.


Important working hints

- Substances such as retraction liquids can adversely affect the setting reaction and should therefore be thoroughly rinsed off before taking the impression.
- Do not combine with condensation curing silicones.
- Contact with latex examination gloves or other latex material may impair the setting reaction of **Detaseal® hydroflow mono**.
- Be sure to remove all traces of impression material from patient's mouth.
- Cured impression materials are chemically inert – spots on clothing should be avoided.
- For the one-step putty wash technique we recommend **Detaseal® hydroflow lite** or **-Xlite** as precision wash material for single-stage impression-taking.
- To ensure trouble-free use we only recommend using dynamic mixers from DETAX following the mixer manufacturer's instructions
- In impression-taking situations with deep undercuts and wide open interdental spaces the usual blocking-out procedures must be performed.
- Please follow the instructions of the safety data sheet!

Further information:

Silicone based materials are proven a million times. On condition of a proper application, undesired effects are not to be expected. However, reactions of the immune system like allergies, irritations, cannot be absolutely excluded. In case of doubt, we recommend to make an allergy test before the application of the material.

Only for dental use by qualified personnel.

Indications for use:

- Single-mix impression technique
- Double-mix impression technique
- Situation impressions
- Over-ring and fixing impressions
- Impressions for indirect denture relining
- Transfer impressions
- Tray material for the sandwich and double mix impression technique

Technical data:

- DIN EN ISO 4823 - Type 2
- Mixed volume:** maxi cartridge: 380 ml, base + catalyst
 - Mixing ratio:** 5:1
 - Product colour:** base: yellow catalyst: white
 - Working time:** approx. 1 min. 30 sec.*
 - Setting time:** approx. 3 min. 30 sec.*
 - Time in the mouth:** approx. 2 min.
 - Strain in compression:** < 3,0 %
 - Recovery from deformation:** > 99,8 %
 - Linear dimensional change:** < 0,2 %
 - Application:** At 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity.
 - Storage:** at ambient temperature

* from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity. Increased temperatures accelerate, decreased temperatures retard these times.

Ordering information:

- Detaseal® hydroflow mono**
Jumbo packing 03373
maxi cartridge of 380 ml, base + catalyst
10 dynamic mixers, yellow
1 rotating ring, yellow
- Standard packing 02754
2 cartridges of 50 ml
6 mixing cannulas, green
- Multi packing 02757
4 cartridges of 50 ml
12 mixing cannulas, green
- Detaseal® hydroflow putty**
Jumbo packing 03283
maxi cartridge of 380 ml, base + catalyst
10 dynamic mixers, yellow
1 rotating ring, yellow
- Standard pack 02727
2 jars of 250 ml
2 measuring scoops
- Detaseal® hydroflow heavy**
Jumbo packing 02729
maxi cartridge of 380 ml, base + catalyst
10 dynamic mixers, yellow
1 rotating ring, yellow
- Detaseal® hydroflow lite**
Standard packing
2 cartridges of 50 ml
12 mixing cannulas, yellow
regular set 02730
fast set 02733
- Detaseal® hydroflow Xlite**
Standard packing
2 cartridges of 50 ml
12 mixing cannulas, yellow
regular set 02741
fast set 02744
- trayloc® A**
17 ml bottle 03098
with brush 03392
15 ml spray bottle

DETAX

GmbH & Co. KG
Carl-Zeiss-Str. 4 · 76275 Ettlingen/Germany
Telefon: 0 72 43/510-0 · Fax: 0 72 43/510-100
www.detax.de · post@detax.de

Caution: Federal U.S. law restricts this device to sale by or on the order of a dentist (or trained specialist personell).

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