
Rinn® Eezee-Grip®

Universal Sensor Holder / Porte-captteur
universel / Porta-sensore universale / Universal-
Sensorhalter / Soporte de sensor universal /
Universele sensorhouder / Uniwersalny uchwył
czujnika / Suporte de sensor universal / Βάση
αισθητήρα γενικής χρήσης / Anturin yleispidike /
Universal sensorholder / Universal sensorholder /
Universalsensorhållare

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
ENGLISH

Indications for use: Eezee-Grip Universal Sensor Holder is indicated for holding digital sensors for intraoral periapical dental radiographs.

Contraindications: None known.

Warnings: Sterilize device before use. DO NOT use phenol-based glutaraldehyde or ultrasonic cleaners. DO NOT chemiclave, dry-heat sterilize or cold sterilize. Only sterilize with steam sterilization according to the provided sterilization instructions. DO NOT exceed 273°F (134°C). It is the user's responsibility to qualify alternative methods of disinfection and sterilization. Failure to follow properly validated sterilization processes and approved aseptic techniques may result in cross contamination. Any method of sterilization will shorten the life of plastic parts.

Precautions: Follow sterilization instructions carefully to prevent damage to plastic during sterilization. Always use a sensor cover to protect sensor.

Adverse Reactions: None known.

Residual Risks: Failure to follow the validated reprocessing instructions provided may result in exposing the patient to cytotoxic residuals and/or cross-contamination.

Step-by-Step Instructions:

1. Insert digital sensor according to photographs.
2. Position holder in patient's mouth and align x-ray tube.
3. Take exposure.
4. Adjust as needed for the next desired position.

Reprocessing: Disassemble instrument(s) and separate all parts prior to reprocessing.

Cleaning: Manual / Mechanical (Ultrasonic)

1. Rinse the disassembled device under potable running water while brushing with a soft-bristled brush to remove gross soil.
2. Prepare an enzymatic instrument cleaning solution and immerse the device per the manufacturer's instructions.
3. Rinse the device under potable running water for at least 30 seconds and dry using a disposable lint-free cloth.
4. Visually inspect the device for visible soil. If soil is seen, repeat the cleaning procedure.

Cleaning and Disinfection: Automated – Instrument Washer-Disinfector
Use only properly maintained, calibrated, and approved washer-disinfectors according to ISO 15883.

1. Rinse the disassembled device under potable running water while brushing with a soft-bristled brush to remove gross soil.
2. Follow manufacturer instructions for instrument washer solution, observing concentration rates and contact times as well as the instrument washer's recommended cleaning and/or thermal disinfection cycle (1 minute at 90°C).
3. Visually inspect the devices for visible soil and repeat the cleaning and disinfection procedure if necessary.

Sterilization: Ensure that the device is thoroughly dry prior to sterilization. Enclose the device into an appropriately sized sterilization pouch intended for steam sterilization. Place the pouched device into the sterilizer per the manufacturer's instructions. Sterilize the device according to one of the recommended methods listed below.

Sterilize metal and plastic components in separate pouches to avoid damage to the components. DO NOT exceed 273°F (134°C).

Steam Sterilization Cycle	Minimum Temperature	Minimum Full Cycle Exposure Time	Minimum Dry Time
Pre-vacuum	270°F (132°C)	3 min	16 min
Steam-Flush Pressure-Pulse	270°F (132°C)	3 min	16 min
Gravity Steam Sterilization	270°F (132°C)	10 min	15 min

Alternate Method: Place the non-pouched device into the sterilizer and sterilize according to one of methods listed above. Non-pouched device must be used immediately following sterilization.

Storage: To maintain sterility, the device should remain in the sterilization pouch until ready for use.

Any serious incident that has occurred in relation to this device should be reported to the manufacturer and the Competent Authority of the Member State.

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©2019 Dentsply Sirona
MDSS GmbH