Etching Gel

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English

Product description

DMG Etching Gel is a universal etching gel (37% phosphoric acid gel) with optimal flow and thixotropy characteristics. Due to its blue color DMG Etching Gel is easily distinguishable from the tooth and provides reliable conditioning results.

Indications

Etching of enamel and dentine

Basic safety instructions

- For dental use only!
- Keep out of reach of children!
- In case of contact with the soft tissue of the mouth, eyes or skin rinse immediately with plenty of water, and, if necessary, consult a physician.

Instructions for use

- Make sure that the etched surface is not touched or contaminated with saliva until the next treatment step (e. g. with the application of a bonding agent). If contamination occurs after drying re-etch for approx. 10 seconds.

Recommended use

1. Clean the tooth surface prior to preparation. Remove any residue with water spray. Prepare the cavity according to the general rules for adhesive dentistry. Bevel enamel margins. We recommend the use of a rubber dam. Protect pulp in deep cavities with a suitable underfilling, such as a small amount of calcium-hydroxide and a thin layer of glass ionomer cement.
2. Apply DMG Etching Gel on the cavity's enamel margins (the use of the Mikrotip is explained below), and let set for 20 – 60 seconds. When using the Total Etch method start the application at the enamel and move to dentine so that the enamel is etched for 20 to 60 seconds but the application time on the dentine does not exceed 15 seconds. The etched enamel should have a chalky white appearance. If this is not the case, the etching process must be repeated. Please follow the manufacturer's instructions for the materials used in combination with the DMG Etching Gel.
3. Rinse with water for at least 15 seconds. Then dry with oil-free and water-free air.

Storage and shelf life

- ▶ Do not store above 25 °C/77 °F.
- ▶ Do not use after expiration date.

Composition

o-phosphoric acid in water (37%), pyrogenic silicic acid

Available product sizes

2 tips @ 2 ml etching gel, 6 Luer-Lock-Tips	REF 212873
---	------------

Accessories

50 Luer-Lock-Tips, White	REF 109942
--------------------------	------------

Federal law (USA) restricts this device to sale by, or on the order of a dentist, or other practitioner licensed by law of the state in which he or she practices to use or order the use of this device.

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