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FANTESTIC® POST

CE 0482

Hersteller

R-dental Dentalerzeugnisse GmbH - Winterhuder Weg 88 - D-22085 Hamburg
Telefon +49 (0)40 - 22 75 76 17 - Fax +49 (0)40 - 22 75 76 18 - www.r-dental.com

FANTESTIC POST

Bestell-Nr.: FPS2010 (zylindrisch, 1,0 mm, inkl. Bohrer FPB2010)
Bestell-Nr.: FPS2012 (zylindrisch, 1,2 mm, inkl. Bohrer FPB2012)
Bestell-Nr.: FPS2014 (zylindrisch, 1,4 mm, inkl. Bohrer FPB2014)
Bestell-Nr.: FPS2017 (zylindrisch, 1,7 mm, inkl. Bohrer FPB2017)

Bestell-Nr.: FPS2510 (konisch, 0,80-1,25 mm, inkl. Bohrer FPB2510)
Bestell-Nr.: FPS2512 (konisch, 1,00-1,45 mm, inkl. Bohrer FPB2512)
Bestell-Nr.: FPS2514 (konisch, 1,20-1,65 mm, inkl. Bohrer FPB2514)
Bestell-Nr.: FPS2517 (konisch, 1,40-1,85 mm, inkl. Bohrer FPB2517)

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Description and indication

FANTESTIC® POST is a root canal post based on glass fibre in a resin matrix with optimal biomechanical and biocompatible characteristics. It is slightly radioopaque (comparable to a titanium post). FANTESTIC® POST has on a biomechanical point of view a dentinlike elasticity modulus (E-modulus). Therefore occlusal loads of the post build-up will be transferred non-traumatically within the tooth. FANTESTIC® POST is conically tapered and cylindrically shaped.

Characteristics and composition

FANTESTIC® POST consists of silanized glass fibres in a resin matrix, which guarantees an optimal internal adhesion and chemical adhesion to the adhesive composite. The non-metallic material guarantees a good biocompatibility.

Recommended procedures**1. Isolation**

A dry area and use of rubber dam (Order-No.: KDL31XX, XX = depends on strength) or latex-free dam (Order-No.: KDL3500) is recommended. Prepare the root canal lege artis. Do not use eugenol containing products!*

2. Adaptation of the post

Take the FANTESTIC® POST which corresponds to the geometry of the prepared root canal with a forcep and prepare its length for an adaptation within the endometrially measured root canal (localization of the apex). The post should penetrate about 6-7 mm within the root canal. It can be shortened with a diamond instrument.

Caution: The post is silanized and should not be touched with fingers. If needed clean the post with ethanol. Before adapting the post a silanization e. g. with FANTESTIC® CerBond (Order-No.: FCF4100) is recommended.

Further procedure according to**3. / 3 a. Total-Etch-Technique (Etch & Rinse) or 4. Self-Etch-Technique****3. Total-Etch-Technique (Etch & Rinse)**

Apply FANTESTIC® ETCH (Order-No.: EGS3000) into the root canal on dentin (etch dentin for about 10 s). Then rinse about 20 s with water and dry with oilfree air. The dentin should stay wet!

3 a. Application of a dual-cure adhesive like FANTESTIC® FLOWSIVE respectively FLOWSIVE 2 each with FANTESTIC® FLOWSIVE ACTIVATOR

Equivalent amounts of FANTESTIC® FLOWSIVE (Order-No.: FFF3100) respectively FANTESTIC® FLOWSIVE 2 (Order-No.: FFF3105) will be mixed for dual-cure each with FANTESTIC® FLOWSIVE ACTIVATOR (Order-No.: FAF3150) light protected for about 15 s. Apply the mixture by using MICRO BRUSHES (Order-No.: MPB3981) or a paper stick into the root canal and agitate intensively for about 30 s and aerate for about 15 s in a weak air stream. It will auto-cure within about 3 min (37 °C, 99 °F). Continue the work immediately.

4. Application of the adhesive system FANTESTIC® FLOWSIVE SE (Self-Etch-Technique)

FANTESTIC® FLOWSIVE SE (Order-No.: FFF5155) Liquid A and Liquid B will be mixed 1:1 within the enclosed mixing mould light protected for about 10 s. The homogenous mixture at this point is self-etching (SE). Then the mixture will be applied with a brush or a paper stick into the root canal and agitate intensively for about 30 s and aerate for about 15 s in a weak air stream. It will auto-cure within about 5 min (37 °C, 99 °F).

5. Application of the luting material

For luting the silanized post the automatically mixed FANTESTIC® FLOWFILLCEM DC (Order-No.: FFS326X, X depends on shade) or FANTESTIC® CORE DC (Order-No.: FCX3X6X, X depends on system and shade) composite or the dual-cure luting cement FANTESTIC® DirectCem DC (Order-No.: FDS426X, X depends on system and shade) will be applied within the working time of about 1 min 30 s directly into the root canal. Then the post will be set. The self-cure system will auto-cure within about 3 or 4 min (depends on luting material) after the working time. On the core side the luting material can be light-cured optionally for about 40 s. In the case the cavity is circular open, the core can be built-up directly with FANTESTIC® CORE DC.

If using cold composites (however do *not* store the luting cement DirectCem DC within a refrigerator) the working time (attention for prolonged hardening time!) can be prolonged.

Precautions to application, product life and/or storage

Avoid contact of the etching gel with soft tissues, skin or eyes. If there is an accidental contact, flush immediately with numerous amounts of water. Avoid contact of the FANTESTIC®-luting materials with soft tissues, skin or eyes. Do not use the FANTESTIC®-luting materials in case of allergic reactions.

Pay attention to the material safety data sheets (MSDS >> www.r-dental.com).

*Do not store FANTESTIC®-luting materials in proximity of eugenol containing materials nor let it come in contact with eugenol containing materials. Eugenol can impair hardening.

Attention: For a self-curing application with luting composites always a self(dual)-curing adhesive must be applied! This does not apply to the self-adhesive luting cement FANTESTIC® DirectCem DC.

Shelf life see label. Do not use after expiry date. Keep out of the reach of children!

Information for reprocessing**FANTESTIC® POST drill**

This information applies to the drill of each system kit.

Drills delivered in non-sterile condition have to be processed that means prepared prior to first use.

Recommended procedures of reprocessing the drill

After the surgical-invasive treatment the drill should be placed immediately into a drill disinfection bath.

Manual cleaning procedure

1. Remove mechanically dried residues immediately after treatment from the drill (especially the hollow spaces) using a nylon brush.
2. Then rinse the drill with running water.
3. Place it into an ultrasonic device filled with a detergent/disinfectant.
4. Take the drill off, rinse it with running water and dry it with compressed air.
5. Examine the drill for being clean. If necessary repeat the procedure.

Sterilisation procedure

1. Place the pre-cleaned suitably packed drill onto a sterilisation tray.
2. The drill must be placed mechanically protected on the tray.
3. Take a steam sterilisation device using a fractionated vacuum process that complies with the provisions of DIN EN 13060 with validated processes.
4. Choose a fractionated pre-vacuum (type B).
5. Make sure to use a sterilisation process temperature of 134 °C (273 °F).
6. The process time should be at least 5 min, the drying time at least 10 min.
7. Follow completely the instructions for use of the sterilisation device manufacturer.
8. After the process protect the drill from dust, moisture and recontamination during transport and storage.

FANTESTIC is a registered trademark of R-dental Dentalerzeugnisse GmbH

Caution: Pay attention to the legal provisions of your country for sterilisation.

FANTESTIC, Cal de Luxe, ZM-Bond and KwikkMix are registered trademarks of R-dental Dentalerzeugnisse GmbH.

Caution: For dental use only.

The described product must be used according to the instructions for use and is determined for dental use only. There is no liability for deflections resulting from inappropriate use or not intended purposes. The user is responsible for testing the suitability of the product for any purposes not explicitly stated in the instructions for use.

If there are damages to this product within the warranty period, R-dental is obliged to replace or to repair this product. There is no further claim. Outside of the liability by law there is no liability of R-dental for any loss or damages arising from this product, independently whether these damages are direct, indirect, special, incidental or consequential.

© 2013 R-dental Dentalerzeugnisse GmbH. All rights reserved.