# DE

# EN

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## Specific advices:

- only for dental use
- FAVORIT powder and liquid should be used only in combination or with components formerly known as Kulzer Carboxylate powder and liquid
- the best adhesion occurs when applied to dry tooth structure
- protect the powder from moisture
- protect the liquid from sun light
- store between 2–28 °C
- do not use after expiration date
- additional information and material safety data sheet (MSDS) on [www.favodent.de](http://www.favodent.de)

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**Composition:**

Powder: zinc oxide, magnesium oxide, feldspar

Liquid: polymeric acrylic acid, water

**Indications:**

- permanent fixations of crowns, bridges, inlays
- cavity liner for protection against thermal and chemical irritants

**Handling and dosage:**

In order to mix the powder and liquid correctly, a glass plate with a temperature of 18–20°C (64–68°F) and a spatula (classically made of tantalum or agate) should be used.

Ensure that spatula and glass plate are clean.

**1. Dosage:**

- For luting: 1 g liquid + 2,9 g powder
- cavity linings: 1 g liquid + 4,7 g powder

**1. Mixing:**

Quickly add the portions of powder to the correct amount of liquid.

Once mixed, carboxylate cement has a thicker consistency than phosphate cement, this does not affect its ease of application and, when used as a luting cement, the flow properties remain excellent.

The ratios quoted for use as a cavity liner result in a condensable consistency.

Mixing time: 30 seconds

**2. Working time:**

- For luting: 2 minutes
- For cavity linings: 1 minute

**Setting-time:**

- For luting: 3–4 minutes
- For cavity linings: 2–2.5 minute

The cement can be trimmed once the setting time has elapsed.