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Hazard statements

H317 May cause an allergic skin reaction. H413 May cause long lasting harmful effects to aquatic life.

Precautionary statements

P261 Avoid breathing dust/fume/gas/mist/vapours/spray. P273 Avoid release to the environment. P280 Wear protective gloves. P302+P352 IF ON SKIN: Wash with plenty of water. P333+P313 If skin irritation or rash occurs: Get medical advice/attention.

inserted in situ within the above mentioned processing time, and is to be hold in position until reaching the flexible phase (approx. 2 min) under constant pressure. Afterwards the temporary can be removed along with the impression, but it is recommended to let it remain inside the impression for about 5 min (end of polymerisation).

b) Indirect manufacturing

For manufacturing the temporary in a dental laboratory a precise situation model of the state before the preparation is required, as well as a prep-model. A thermoforming blank (e. g. Drufofen H) is placed onto the model in the area of the teeth that need to be prepared; subsequently a matrix is cut out and filled with Fixtemp C&B directly from the double cartridge with the mixing tips. 30 seconds after beginning of mixing, the Fixtemp C&B-filled thermoforming blank is set on the insulated preparation model (e. g. with Isolat Film) and fixed safely at the remaining teeth. After 2 minutes, the flexible phase of Fixtemp C&B is reached, and the thermoforming blank can be removed carefully along with the temporary. Until the material has ultimately hardened after 5 minutes, the temporary is recommended to remain inside the thermoforming blank.

Final Treatment

- The inhibition layer on the surface, which has been created by the aerial oxygen, can be removed from the not-abraded areas with the help of organic solvent (e. g. isopropyl alcohol).
- The temporary is to be treated with fine carbide-trephines after it has hardened completely (at least 5 minutes after begin of mixing). Do not breathe in abrading dust! Use mask or suction!
- Apply light-curing high-gloss lacquer (e. g. Plaquit) or polish manually with rubber-polisher respectively caprine-hair brush.

Cementation

Cements containing Eugenol can interfere the curing composites when fixing. Insert temporary with usual temporary fixing cements (e. g. Fixtemp cement).

Disposal

P501 Dispose of contents / container in accordance with local / regional / national / international regulations.

Important advice

- Due to technical filling an unbalanced mixing of components can occur when using new cartridges. For this reason it is recommended to dispose the initially mixed pea-sized amount.
- In case of reusing a started double cartridge the outlets have to be cleaned of possible blockages before attaching the mixing tips.
- Disinfection with usual material is possible without any problems.
- This product is only to be used in the above mentioned application area. Only professional and introduced persons are allowed to use this product.

Keep away from children!

Only for dental use!

Contraindications

The material should not be used for allergies to any of the ingredients and / or contact allergies.

Side effects

So far no systematic side effects are known.

Disinfection of the Injector 4:1 / 10:1

The disinfection can be carried out with a common immersion disinfection and is also autoclavable.

Storage conditions

Store at room temperature that means not over 25 °C. Do not store in refrigerator!

Durability:

see date of expiry
Fixtemp C&B is not to be used when expiry date has passed.

Caution

Fixtemp C&B contains methacryl acid ester. In case that allergic reactions according to this component are known, do not use!

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