# FLEXISTONE®

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**FLEXISTONE® cat** **03165**  
Endurecedor, botella de 30 ml

# FLEXISTONE®

Insulating silicone material for use in dental technology,  
condensation curing, low viscosity, high final hardness

Fig. 1

Fig. 2

Fig. 3

## 1. Mixing ratio and mixing

Fill mixing bowl with **FLEXISTONE®** to the second ring marking (approx. 12 g), add 35 – 40 drops of catalyst (fig. 1) and mix homogeneously by means of a spatula (fig. 2).

## 2. Application

Apply **FLEXISTONE®** in an even layer of approx. 2 mm thickness to the acrylic or porcelain teeth, beginning with the interdental spaces. Before the material starts to set, provide undercuts with an instrument (fig. 3). Thus a firm adhesion to the plaster is guaranteed. No retention crystals are necessary.


## 3. Deflasking

After polymerisation deflask as usual and remove **FLEXISTONE®**.

## Important working hints

- Avoid thin-tapered borders – cutting of a small shoulder is recommended to prevent separation in the flask-counter.
- Avoid contact of catalyst with eyes, may cause irritations. In case of accidental contact immediately flush under running water for at least 15 minutes and seek immediately medical care.
- Please follow the instructions of the safety data sheet!

Only for dental use by qualified personnel.

**GB**  
Instructions for use

## Indications for use:

- Insulating material, C-silicone based, for the application with flask pressing technique. For the insulation and protection of acrylic or porcelain teeth during investing. Can be easily applied in one single procedure onto acrylic teeth and wax, forming a thin, very hard insulating layer for fixing of teeth during pressing.

## Technical data:

- **Mixing ratio:** approx. 12 g / 35 – 40 drops
- **Colour code:** **FLEXISTONE® paste:** green  
**FLEXISTONE® cat:** colourless
- **Mixing time:** approx. 45 sec.
- **Working time:** approx. 3 min.\*
- **Setting time:** approx. 10 min.\*
- **Final hardness:** approx. 75 Shore A
- **Application:** at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity
- **Storage:**

\* from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5% rel. humidity. Increased temperatures accelerate, decreased temperatures retard these times.

## Ordering information:

**FLEXISTONE®**  
Standard packing **03009**  
Paste, bottle of 190 ml (350 g)  
Catalyst, bottle of 30 ml  
1 Mixing bowl  
1 Stand

Large packing **03010**  
Paste, bottles of 5 x 190 ml (5 x 350 g)  
Catalyst, bottles of 3 x 30 ml  
2 Mixing bowls  
1 Stand

**FLEXISTONE® cat** **03165**  
Catalyst, bottle of 30 ml

# DETAX

GmbH & Co. KG  
Carl-Zeiss-Str. 4 · 76275 Ettlingen/Germany  
Telefon: 07243/510-0 · Fax: 07243/510-100  
www.detax.de · post@detax.de

Made in  
Germany

08/2017

# FLEXISTONE®

recommended by IVOCLAR  
for the SR-Ivocap-System

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