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Instructions for Use

Y-Holder for Articulating Silk and Occlusion test films BK 140

Bausch
WE MAKE OCCLUSION VISIBLE**Manufacturer**

Dr. Jean Bausch GmbH & Co. KG
Oskar-Schindler-Str. 4
50769 Cologne - Germany
Tel: +49 221 70936-0
Fax: +49 221 70936-66
E-mail: info@bauschdental.de
Internet: www.bauschdental.de

Date of issue: 2019-12-06

1 General Description and purpose of the product

The Y-holder is a laboratory instrument, made of stainless steel, for easier handling of wide articulating silks or occlusion test foils in the articulator.

The Y-holder is not intended for use on humans and is not a medical device!

The Y-holder is sterilizable.

Indication

Laboratory-Instrument.

Contraindication

Not intended for use on humans!

Side effects

N/A.

2 Notes

- Before use, the Y-holder (instrument) has to be cleaned, disinfected and sterilized, if necessary.
- The instrument can be sterilized with steam (moist heat, 134°C).
- The instrument can be used in combination with articulating silks or occlusion test foils, produced by Dr. Jean Bausch GmbH & Co. KG.
- Disposal: The instrument must be disposed of with the usual laboratory waste.

3 Package content

- 1 Y-Holder (BK 140)
- Instructions for Use

4 Preparation

The instrument should be cleaned before use. This applies in particular to initial use after delivery, as the instrument is delivered non-sterile.

5 Application

- Insert the occlusion test material (e.g. silk or foil) into the holder so that it protrudes over the clamping jaws on all sides.

- Close the clamping jaws by pushing the locking slider on the handle forward.
- Perform the occlusion test on the model in the articulator.
- After the occlusion test has been completed, open the Y-holder by pushing the locking slide back and dispose of the used occlusion test material with the usual laboratory waste.
- If necessary, prepare the Y-holder for reprocessing.

6 Disposal Instructions

The instrument has to be disposed of with the normal laboratory waste.

7 Symbols

- Manufacturer
- Production Date
- „Follow manufacturers instructions“
- Reference number
- Unique Device Identifier
- Conformity with the relevant EU regulations

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