Dreve Design Rev.070319-3990

- Remove all bubbles with a cleaned object.
- Processing temperature:

RapidShape DII printer	23 ± 3 °C
Asiga MAX UV	35 ± 3 °C
- Select the product specific parameter settings for FotoDent® IBT for your process.
- When the building process is finished a direct post treatment is recommended.
- After the platform is taken up a dripping off time of approx. 10 minutes is recommended.
- Resin coated parts get clean with isopropanol (97 %) within a time of about 2 x 6 minutes. For a particular effective cleaning with isopropanol the use of explosion-protected ultrasonic units is recommended.
- For post-curing of the parts, please select the parameter settings for your process. You can find more qualified printers on our website www.dreve.dentamid.de.
- Post-curing:

PCU LED N ₂	RapidShape DII
	12 minutes, 40% light output using nitrogen
	Asiga MAX UV
	10 minutes, 80% light output using nitrogen
- FotoDent® flash

RapidShape DII / Asiga MAX UV
12 minutes + 1 minute gas pre-flow
- Otoflash G171

RapidShape DII / Asiga MAX UV
4000 flashes with program II using nitrogen and UVB-blocking tray
- Impurity due to the construction or a break of the material and impurity due to operation mistakes cannot be excluded. Thanks to the low viscosity however it is possible to filtrate FotoDent® IBT. It is recommended to regularly take out the container of the production unit, to homogenize and filtrate the content. Upcoming entrained bubbles get eliminated through a holding time of approx. 1–2 hours with FotoDent® IBT.

Possible dangers

Please note the safety instructions for processing FotoDent® IBT.

Hazard statements: **H315** Causes skin irritation. **H319** Causes serious eye irritation. **H317** May cause an allergic skin reaction. **H411** Toxic to aquatic life with long lasting effects.

Precautionary statements: **P261** Avoid breathing dust/fume/gas/mist/vapours/spray. **P280** Wear protective gloves/eye protection/face protection. **P273** Avoid release to the environment. **P305+P351+P338** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. **P333+P313** If skin irritation or rash occurs: Get medical advice/attention. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Contents: Aliphatic difunctional Urethane Acrylate, (5-ethyl-1,3-dioxan-5-yl) methyl acrylate, 2-phenoxyethyl methacrylate, Diphenyl(2,4,6-trimethylbenzoyl) phosphinoxide

Contra indication

Intended use of FotoDent® IBT is for manufacturing dental bonding trays only; an allergic reaction might occur in case of any allergy to (meth)acrylates.

Lot number / Expiry date

The lot number and the expiry date are indicated on each FotoDent® IBT packaging. In case of claims please always indicate the lot number of the product. Do not use the product after the expiry date.

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Intended use

FotoDent® IBT is a material based on (meth-)acrylate resin for DLP systems with 385 nm LED for the manufacturing of dental indirect bonding trays. RapidShape DII series and Asiga MAX UV are appropriate printer. For the post curing of FotoDent® IBT built parts, the light curing units PCU LED N₂, FotoDent® flash and Otofash G171 are recommended, which are optimal balanced to the production process.

Indications for use

For manufacturing dental indirect bonding trays to place brackets.

Please note

- Due to the flexibility of the FotoDent® IBT material it is recommended to built the object directly on the platform without support structures in an horizontal orientation. Therefore, the occlusal area of the indirect bonding tray has to be designed as a flat surface without chewing surface structure. Details can be found in the instruction for use of your design program.
- Due to different bracket shapes, please be aware of sufficient fitting during designing IBT. If not the brackets may fall out of lower jaw IBT during placing IBT.
- To avoid damage of material quality do not expose the liquid material to light under any circumstances. Variations from the described manufacturing process may impair biocompatibility, lead to modified mechanical characteristics and/or color variations of the FotoDent® IBT material.

Processing

- During the application of FotoDent® IBT we recommend wearing personal protective equipment e. g. suitable gloves, safety goggles etc.
- Please see further steps within printer manual.
- Carefully pour FotoDent® IBT into the pre-determined container of the production unit.

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