Dreve Design Rev.200820-3939

• For post-curing of the parts, please select the parameter settings for your process. You can find more qualified printers on our website www.dreve.dentamid.de.

• Post-curing **385 nm**:
PCU LED N₂
Rapid Shape DII / Asiga MAX UV 5 minutes, 33 %
light output (incl. nitrogen)

FotoDent® flash
Rapid Shape DII / Asiga MAX UV 6 minutes + 1 minute
gas pre-flow

RS Cure (Gas)
According to respective parameter set for FotoDent® guide
385 nm.

Post-curing **405 nm**:
FotoDent® flash
Rapid Shape D 10 minutes + 1 minute
gas pre-flow

• Impurity due to the construction or a break of the material and impurity due to operation mistakes cannot be excluded. Thanks to the low viscosity however it is possible to filtrate FotoDent® guide. It is recommended to regularly take out the container of the production unit, to homogenize and filtrate the content. Upcoming entrained bubbles get eliminated through a holding time of approx. 1–2 hours with FotoDent® guide.

To avoid damage of material quality do not expose the liquid material to light under any circumstances. Variations from the described manufacturing process may impair biocompatibility, lead to modified mechanical characteristics and/or color variations of the FotoDent® guide material.

Possible dangers

Please note the safety instructions for processing FotoDent® guide.

Hazard statements: H319 Causes serious eye irritation. **H317** May cause an allergic skin reaction. **H412** Harmful to aquatic life with long lasting effects. (REF D35600 only) **H413** May cause long lasting harmful effects to aquatic life. (REF D35650 only)


Precautionary statements: P261 Avoid breathing dust/fume/gas/mist/vapours/spray. **P273** Avoid release to the environment. **P280** Wear protective gloves/eye protection/face protection. **P305+P351+P338** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. **P333+P313** If skin irritation or rash occurs: Get medical advice/attention. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Contains: 2-hydroxyethyl methacrylate, hydroxypropylmethacrylate, monoester with propane-1,2-diol, Diphenyl(2,4,6-trimethylbenzoyl) phosphin oxide

Lot number / Expiry date

The lot number and the expiry date are indicated on each FotoDent® guide packaging. In case of claims please always indicate the lot number of the product. Do not use the product after the expiry date.


Glossary

	Manufacturer		Protect from sunlight
	Batch Code		Temperature limitation
	Item number		Medical device
	Expiry date		Use by only professionals
	CE marking of conformity according to the requirements for medical devices		Please note instructions for use

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Intended use

FotoDent® guide 385 nm or 405 nm are materials based on (meth-)acrylate resins for DLP systems with 385 nm / 405 nm LED for the manufacturing of dental drill guides for tooth-supported and /or teethridge-supported application. For the post curing of FotoDent® guide built parts, the light curing units PCU LED N₂ and FotoDent® flash are recommended, which are optimal balanced to the production process.

Indications for use

Manufacturing dental drill guides for placing dental implants.

Processing

- Alignment of drill guide: vertical, supports, connection of supports to a baseplate.
- During the application of FotoDent® guide we recommend wearing personal protective equipment e. g. suitable gloves, safety goggles etc.
- Please see further steps within printer manual.
- Carefully pour FotoDent® guide into the pre-determined container of the production unit.
- Remove all bubbles with a cleaned object.
- Processing temperature:
RapidShape D and DII series $23 \pm 3^\circ\text{C}$
Asiga MAX UV $35 \pm 3^\circ\text{C}$
- Select the product specific parameter settings for FotoDent® guide for your process.
- When the building process is finished a direct post treatment is recommended.
- After the platform is taken up a dripping off time of approx. 10 minutes is recommended.
- Cleaning and post curing incl. supports and baseplate to avoid warpage.
- Resin coated parts get clean with isopropanol (97 %) within a time of about 2x6 minutes. For a particular effective cleaning with isopropanol the use of explosion-protected ultrasonic units is recommended.