# Granisit® XF Speed

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DIN EN ISO 15912, Typ 2, Klasse 2

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DIN EN ISO 15912, type 2, class 2

## Description:

Granisit® XF Speed is a phosphate bonded, graphite free and very fine grain precision investment for partial denture frameworks. It is usable for the conventional and speed casting technique and is suitable for casting all types of dental alloys.

## Casting models:

Smear Vaseline on the inside of the overlapping section of the duplicating sleeve to prevent cracking. Measure the required amount of investment and the correct concentration of expansion liquid according to the mixing ratio for casting models. Mix manually with a spatula before vacuum mixing for 60 sec. Use small vibration for pouring the contours. An instrument is useful for this. Switch off the vibrator and pour the rest of the model, let it harden shock-free in the form for approx. 30 minutes.

## Technical data:

Mixing liquid:	SILADENT type 100 expansion liquid
Mixing ratio, Powder : liquid:	100 g : 20 ml
Mixing under vacuum:	60 sec.
Working time:	4 - 5 min.
Working temperature:	approx. 22 °C (Room temperature)
Initial setting time:	7 - 10 min.
Setting expansion:	approx. 1.00% (75% diluted type 100 liquid)
Preheating temperature range for placing in the furnace:	850 - 900° C
Place in the hot furnace:	25 min. after investing

## Concentration of the expansion liquid:

20 ml liquid / 100 g powder

## Frameworks with clasps

Liquid ratio Typ 100 / distilled water: 75 % / 25 %

100 g	180 g	400 g
15 ml / 5 ml	27 ml / 9 ml	60 ml / 20 ml

## Fixed / removable work

Liquid ratio Typ 100 / distilled water: 80 % / 20 %

100 g	180 g	400 g
16 ml / 4 ml	29 ml / 7 ml	64 ml / 16 ml

### Investing:

Measure and mix Granisit® XF Speed as described above. For the investing use the same liquid concentration as for the model. A small reduction can optimize the adhesion between model and investing. Pour with minimum vibration until the pattern is completely covered. Switch off the vibrator, fill the ring and do not vibrate again. Put the muffle immediately after investing in a pressure chamber at approx. **2 bar for 25 minutes** and let it harden vibration free.

### Preheating with the speed casting technique:

25 minutes after investing, place the mould with the sprue hole pointing downwards in a furnace preheated to the correct temperature of 850° C – 900° C). SILADENT recommend a maximum preheating temperature of 900° C. This temperature is perfectly adequate for both casting the mould and final thermal expansion. Moulds are ready to cast approx. 60 minutes after the final temperature is reached. During the first 15 minutes after inserting the moulds, the furnace fan and extractor should remain switched off and the furnace door should not be opened because of the risk of instant combustion.

### Preheating with the conventional casting technique:

Granisit® XF Speed can be used in the conventional casting with following heating rates:

	Temperature	Heating rate	Holding times
1. Heating rate	300° C	5-8° C / min.	30-60 min.
2. Heating rate	600° C	5-8° C / min.	30-60 min.
1. Final temperature	850 - 900° C	7-9° C / min.	30-60 min.

### Please note:

Technical recommendations are based on tests and findings from work in our development laboratory and can only be regarded as guidelines. SILADENT products are subjected to strict quality controls. We reserve the right to make technical changes.

### Warning:

**Investments contain quartz and cristobalite. Do not breathe in dust. Danger of lung diseases (silicosis or cancer). Use a suitable dust mask!**

**Further questions:** Contact SILADENT technical services (Tel.: +49 (0) 53 21-37 79 25/26) or our sales representatives.

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