greenbite apple

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greenbite apple

GB
Instructions for use

Bite registration material, A-silicone based, addition curing fast setting – extreme high final hardness


Fig. 1

1. Mixing and Dosing

Extrude material from the cartridge with Automix2-system gun (Fig. 1). Place cartridge into mixing gun and remove closure cap by rotation. Eject a small amount of material (bleeding) until equal flow is being extruded from both orifices (Fig. 2). Mount the mixing cannula by using the guides on the cannula and the cartridge and lock in place by rotating in the opposite direction (Fig. 3). Feed material applying even pressure. Prior to use, dispense a small amount of material and visually check that the base and catalyst paste are mixed omogeneously (Fig. 4). Only then dosing can be performed individually. Leave mixing cannula on cartridge after use, thus sealing the cartridge.


Fig. 2

2. Instruction for use

Extrude material from cartridge with mixing cannula or with attached contouring tip onto teeth (Fig. 5). At terminal occlusion application from vestibular is possible. Let registration harden in terminal occlusion and remove from the mouth after 1 min.


Fig. 3

3. Disinfection

After removal from the mouth rinse registrate under lukewarm water. Disinfection for 15 min. in 2 % glutaraldehyde or other usual solutions.


Fig. 4

Important working hints

- Remove any retraction astringents or mouth wash solutions prior to the registration – could affect setting reaction.
- Do not intermix with condensation curing silicones.
- Avoid staining of clothing. Hardened impression materials cannot be removed.
- Remove all residues of registration material from the patient's mouth.
- For reducing the pressure needed during injection, we recommend using the green mixing cannula.
- Please follow the instructions of the safety data sheet!


Further information:

Silicone based materials are proven a million times. On condition of a proper application, undesired effects are not to be expected. However, reactions of the immune system like allergies, irritations, cannot be absolutely excluded. In case of doubt, we recommend to make an allergy test before the application of the material.

Only for dental use by qualified personnel.

Fig. 5

Indications for use:

- Bite registration with crown and bridge impressions and for orthodontic indications, key for intraoral registration

Technical Data:

- DIN 13903
- **Mixed volume:** 50 ml (cartridge)
- **Mixing ratio:** 1:1
- **Colour code:** base: green catalyst: yellow
- **Mixing time:** omitted (Automix2-system)
- **Working time:** approx. 30 sec.*
- **Setting time:** approx. 1 min. 30 sec.*
- **Time in mouth:** approx. 1 min.
- **Final hardness:** approx. 45 Shore D
- **Linear dimensional change:** approx. 0.2 %
- **Application:** At 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5% rel. humidity
- **Storage:**

* from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity. Increased temperatures accelerate, decreased temperatures retard a. m. times.

Ordering information:

greenbite apple
Standard packing **02641**
2 x 50 ml cartridges
8 mixing cannulas, green
8 contouring tips
for mixing cannulas green

Ecopacking 02647
8 x 50 ml cartridges
32 mixing cannulas, green
32 contouring tips
for mixing cannulas green

Mixing cannulas 02705
green, 48 pcs.

Contouring tips 02648
for mixing cannulas green
12 pcs.

milkbite
Standard packing **03380**
2 x 50 ml cartridges

greenbite colour
Standard packing **02207**
2 x 50 ml cartridges

glassbite
Standard packing **02571**
1 x 50 ml cartridge

DETAX

GmbH & Co. KG
Carl-Zeiss-Str. 4 · 76275 Ettlingen/Germany
Telefon: 072 43/510-0 · Fax: 072 43/510-100
www.detax.de · post@detax.de

Caution: Federal U.S. law restricts this device to sale by or on the order of a dentist (or trained specialist personell).

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