ENGLISH

Harmonize™: NanoHybrid Universal Composite

DESCRIPTION

Harmonize is a light-cured, esthetic, resin-based nanohybrid dental restorative designed for the placement of direct restorations. Harmonize is designed to be used after placement of a dental adhesive, such as a product from within the OptiBond™ family. Harmonize is available in Unidose composite compules and syringe delivery systems.

INDICATIONS FOR USE

Harmonize, a nanohybrid dental restorative, is indicated for direct placement in all cavity classes in anterior and posterior teeth. Additional indications include: repair of enamel defects, repair of provisional, repair of porcelain restorations, minor occlusal build-ups, core build-ups, and incisal abrasions.

CONTRAINDICATIONS

None known.

WARNINGS

- Do not use this product on patients with known methacrylate allergies. Uncured methacrylate resin may cause contact dermatitis and damage to the pulp.
- Avoid contact with skin, eyes, and soft tissue. Wash thoroughly with water and soap after contact.
- Get medical attention for eyes.
- Do not take internally.

PRECAUTIONS

Kerr recommends the following precautions when using Harmonize or any other composite material:

- Provide the patient with protective eyewear
- Use a rubber dam to isolate the oral cavity (Kerr recommends Kerr OptiDam™)
- Avoid fluid contamination during placement
- Use of protective gloves and a no-touch clinical technique is strongly recommended. If the composite is touched with commonly used gloves, acrylics may penetrate the glove. Dispose of glove immediately and wash hands.
- Protect the material from ambient light that may initiate polymerization

The Unidose tips are designed for SINGLE PATIENT USE ONLY to prevent cross-contamination between patients.

ADVERSE REACTIONS

None known.

INSTRUCTIONS FOR USE

Preparation of the Tooth

- Clean the tooth with pumice and water to remove any debris or stain that may interfere with the shade selection process.
- Select the desired shade. A composite mock-up or VITA-matched shade guide may be used to check the color of the tooth to be restored and any adjacent teeth. In many cases, a single shade of Harmonize in enamel opacity will be sufficient to restore the tooth. The shade of the restorative material is affected by the thickness of the composite. Shade selection and mock-ups should match the thickness of the tooth to be restored.
- After shade selection, isolate the tooth with your preferred isolation method.
- Prepare the tooth using your preferred rotary instruments and using contemporary cavity designs. Harmonize offers excellent blending (chameleon effect) which may be further enhanced by the use of slight bevels on preparation margins. If replacing a restoration, be sure to remove all existing restorative material from the operative site.
- If a pulp exposure has occurred, place a small amount of calcium hydroxide (Kerr recommends the use of Life) followed by the placement of a conventional glass ionomer base.
- Complete the bonding procedure. Harmonize is compatible with all methacrylate dental adhesive systems. OptiBond products have been shown to produce the strongest and most durable restorative interfaces.
- If desired, place a liner with a low-viscosity flowable composite, such as Herculite™ XRV Ultra Flow, and light cure according to manufacturer's instructions.
- Harmonize features optimal, non-sticky handling. It permits the use of hand instruments to define most of the restoration anatomy, which will help minimize the time and effort associated with rotary instrumentation. For **Unidose delivery**: Inject the restorative material from the Unidose applicator into the cavity, starting with the deepest portion. For **syringe delivery**: After extruding material onto pad, back off syringe at least two full turns to relieve pressure on the material.
- Apply composite layers up to 2mm in height for enamel/translucent (incisal) shades or 1.5mm in height for dentin shades, taking care to properly adapt the material to the cavity surface.
- Cure each layer. Refer to the table at the end of these instructions for recommended curing times of the various shades at various light intensities. Avoid exposure of Harmonize to intense optical light.
- After curing, complete the finishing and polishing procedure. Adjust occlusion of the restoration. Kerr Trimming and Finishing Carbides, such as the TDF™ (Trim Define Finish) and OA™ (Occlusal Anatomy) series burr, for finishing and adjustments and Kerr Opti1Step Polisher one-step, multi-use polisher for polishing. Kerr Opti1Step Polisher has been found to produce equivalent polish to other multi-step polishers tested in half the time.

***Recommended Curing Times**
Adequate curing is essential for restorative function. Uncured or undercured composite material is far more likely to experience some form of clinical failure. Kerr recommends regular inspections of your curing lights to ensure they are functioning and providing the expected irradiance. Harmonize is intended to be cured by a halogen or LED light with a minimum light intensity of 600 mW/cm2 and a light wavelength output within the 400-520 nm range. Please refer to the table below for curing recommendations. All times refer to a single occlusal surface cure. In all other posterior restoration, Kerr recommends additional 10 second cures on both the buccal and lingual surfaces of the tooth.

	Dentin Shades	Enamel/Translucent (Incisal) Shades
Incremental Depth	1,5 mm	2,0 mm
OptiLux™ (or light with output 600 – 1000 mW/cm²)		20 Seconds
Demi™ Ultra/ Demi™ Plus (or light with output >1000 mW/cm²)		10 Seconds

Disinfection and Storage

Syringes may be safely disinfected. Use an intermediate level disinfectant, such as Cavicide and Caviwipes, or as recommended by the Centers for Disease Control. Refer to manufacturer's instructions for recommended time, dilution, and solution temperature. Harmonize should be stored at ambient temperature. Avoid excess heat. If product is stored in a refrigerator allow the product to reach room temperature prior to use. Every effort should be made to protect the composite from visible light prior to use.

Shelf Life and Disposal

The expiration date is marked on packaging. Do not use restorative after expiration date. For proper disposal always follow Federal or State laws.

LIMITED WARRANTY - LIMITATION OF KERR'S LIABILITY
Kerr's technical advice, whether verbal or in writing, is designed to assist dentists in using Kerr's product. The dentist assumes all risk and liability for damages arising out of the improper use of Kerr's product. In the event of a defect in material or workmanship, Kerr's liability is limited, at Kerr's option, to replacement of the defective product or part thereof, or reimbursement of the actual cost of the defective product. In order to take advantage of this limited warranty, the defective product must be returned to Kerr. In no event shall Kerr be liable for any indirect, incidental or consequential damage. EXCEPT AS EXPRESSLY PROVIDED ABOVE, THERE ARE NO WARRANTIES, BY KERR, EXPRESS OR IMPLIED, INCLUDING WARRANTIES WITH RESPECT TO DESCRIPTION, QUALITY, OR FITNESS FOR A PARTICULAR PURPOSE.

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