Hera

Heraenium® Sun

DE Gebrauchsanweisung	2	HU Használati utasítás	95
GB Instructions for use	9	LV Lietošanas instrukcija	102
FR Mode d'emploi	16	LT Naudojimo instrukcija	109
ES Instrucciones de uso	24	PL Instrukcja obsługi	116
IT Istruzioni per l'uso	31	RU Инструкция по применению	123
PT Instruções de uso	38	UA Інструкція по застосуванню	131
NL Gebruiksaanwijzing	45		
SE Bruksanvisning	52		
DK Brugervejledning	59		
NO Bruksanvisning	66		
FI Käyttöohjeet	73		
GR Οδηγίες Χρήσης	80		
CZ Návod k použití	88		

KULZER
MITSUBI CHEMICALS GROUP

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Heraenium® Sun Instructions for use

Type 4 alloy (as per EN ISO 22674 and EN ISO 9693)

Heraenium Sun is a cobalt chromium basic alloy for fabricating crowns and bridges for veneering with high expansion veneer-ceramics or veneer-composites. Heraenium Sun bonds well to veneer-ceramics and exhibits optimum physical properties as well as a wide range of indications. Heraenium Sun is easy and straightforward to process.

1. Waxing up: Coping patterns must be at least 0.4 mm thick – this ensures that the coping remains at least 0.3 mm after trimming. To ensure that the ceramic forms to a uniformly thick layer, copings and pontics must be waxed up to approximate the anatomy of natural teeth. Avoid waxing up sharp line angles and undercuts. The minimum wall thickness must also be 0.3 mm for veneering with composites. Use retention beads (\varnothing 0.4 mm).

2. Attaching sprues and investing: We recommend using a reservoir bar sprue system for bridgework. 5 mm long and 3.5 mm diameter feeder sprues are attached to the palatal or lingual surfaces of the wax pattern at an angle of 45°. Each bridge unit requires its own feeder sprue. Large molar copings or solid pontics require two. The feeder sprues must always be attached to the thickest section of the pattern. They are connected by a 5 mm diameter horizontal sprue (bar). 4 mm diameter sprues are used for single connecting sprues. The higher the volume of the casting, the thicker the sprues should be. Fine grained,

phosphate bonded, graphite free and plaster free investment materials for crowns and bridgework such as Heraest Onyx or Moldavest exact are optimised for producing precision fit castings and perfect cast surfaces.

3. Casting: The correct amount of Heraest Sun required for casting is calculated by multiplying the weight of the wax pattern by the density of the alloy divided by the density of the wax (approx. 0.93 g/cm³).

Density of Heraest Sun: 8.2 g/cm³

Formula: Weight of pattern x 8.2 g/cm³ ÷ 0.93 g/cm³.

Casting with induction heated vacuum /pressure casting machines:

Induction heated vacuum/pressure casting machines by Kulzer are the most suitable for melting and casting. During pre-melting, the ingots are melted until the last ingot has sunk into the molten alloy and the edges or dull film of the ingot are no longer visible. Melting is then interrupted, and the casting chamber is allowed to resume normal pressure and opened. The casting ring (mould) is then inserted and the main melting process commenced. The casting process is started **5** seconds after the dull film has disappeared from the surface.

Casting with a motor driven centrifugal casting machine and melting with a casting torch:

Ignite the oxyacetylene gas. Open both valves fully and adjust the flame by reducing the acetylene flow until approximately 3 mm long blue cones appear on the orifices of the casting torch. Pressure settings for the flame: **Acetylene approx. 0.7**

bar, Oxygen approx. 1.9 bar. These settings will produce the neutral flame required for melting Heraenium Sun. Place the Heraenium Sun ingots in the preheated crucible. Hold the torch approximately 40 mm above the top edge of the crucible. Move the torch around the ingots in circles to heat them uniformly until they melt and fuse. Once the ingots have fused, start the centrifugal casting machine immediately. After casting, bench cool the casting ring. Devest the castings very carefully. Remove the investment material with care. Sandblast with aluminium oxide abrasive (110–125µm).

4. Trimming and preparing for ceramic facings: After casting, the copings are trimmed with cross cut tungsten carbide cutters in one direction only, but without exerting high pressure. Do not use diamond burs. Do not use rubber polishers on the surfaces to receive ceramic. Avoid sharp line angles and undercuts. Sandblast and evaporate the surfaces with aluminium oxide abrasive (110–125 µm, 4 bar).

Metal dust is hazardous to health.

When trimming and sandblasting, use a suitable extraction system and/or facemask (type FFP3-EN 149-2001).

5. Oxide firing: 950°C (1742°F) for 10 minutes. The oxide layer must be stained uniformly. Carefully remove the oxide layer on the facing surfaces with aluminium oxide abrasive (110–125 µm, 4 bar) using a non-recycling sandblaster, and then clean (e.g. evaporation). When using either HeraCeramSun ,HeraCeram NP Primer or HeraCeramSun PreOpaque no oxid firing is necessary.

6. Veneering with ceramics

- a) Veneering with HeraCeramSun: The HeraCeramSun processing instructions from Kulzer apply.
- b) Veneering with other ceramics: When veneering with other ceramics, the other ceramic manufacturers' instruction for use must be observed.
- c) Stress-relief cooling: stress-relief cooling is recommended for large-span and more massive frameworks: setting a cooling period with opening of the ceramic furnace (2–4 minutes) or leaving the object on the firing tray in the radiated heat from the open firing chamber.

7. Acrylic facing: Process the facing material according to the manufacturer's instructions.

8. Soldering: Special "Heraenium Sun solder 1" is recommended for soldering Heraenium Sun frameworks. Solder according to the current solder instructions for use.

9. Laser welding: Special CoCr laser welding wire is ideal for welding Heraenium Sun frameworks. \varnothing 0.5 mm, 150 cm. Weld according to the current instructions for use of the laser welding wire.

10. Recasting: The particular physical and chemical qualities of the alloy can be guaranteed only if new material is used. We recommend using only new material for ceramic facings.

11. Applying Blendgold Neu: Blendgold Neu enhances the shades and aesthetics of ceramic. Condition the surfaces as described in sec-

tion 4 and 5 before applying and firing a very thin coat of ceramic layer (i.e.: HeraCeram NP-Primer, HeraCeramSun Pre-Opaque, Opaker etc.). Extrude Blendgold Neu paste out of the syringe directly onto a clean brush and apply a thin coat of Blendgold Neu onto the the first fired ceramic basic layer. Heat the ceramic furnace to 400°C (752°F), and dry the Blendgold Neu in the heat radiating from it. The temperature of the object should reach approximately 180°C (356°F). Raise the temperature to 820°C (1508°F) at a rate of 55°C (131°F) / min. and fire without a vacuum and holding time. Remove the object on the firing tray from the furnace, and allow it to bench cool. Continue processing the ceramic as usual. Simply applying a coat of Blendgold in the cervical region may be adequate. Blendgold Neu cannot be directly applied to Heraenium Sun. If a warm, yellowish gold colour is needed, we recommend using HeraCeramSun intensive opaquer gold.

12. Chemical composition in mass percent

Co: 43,0; **Fe:** 27,0; **Cr:** 23,45; **W:** 2,5; **Mo:** 2,0; **Si:** 1,0; **Mn:** 0,8; **N:** 0,15; **C:** 0,1

13. Technical data

Type 4 alloy following casting and ceramic firing according to EN ISO 22674 and EN ISO 9693.

Density: 8.2 g/cm³

Investment material: Phosphate-bonded e.g. Heravest Onyx, Moldavest exact

Preheating of casting ring: 950°C (1742°F)

Crucible: Ceramic NPM

Melting range: 1290°C–1380°C (2354°F–2516°F)

Casting temperature: 1500°C (2732°F), 5 seconds after dull film disappears

Oxide firing: 950°C (1742°F) / 10 min (oxide firing is not required for facings with HeraCeramSun, using HeraCeram NP-Primer or HeraCeram Pre-Opaque)

CET 25°C–500°C (77°F–932°F): 16.2 $\mu\text{m}/\text{m}^*\text{K}$

Hardness after casting: 265 HV 10

Hardness after firing: 280 HV 10

Yield point 0.2%: after casting 430 MPa
after firing 490MPa

Tensile strength after firing: 800 MPa

Elongation at rupture after firing: 18%

Elasticity modulus: 250 GPa

Contraindications: In cases of hypersensitivity (allergy) to the constituents of the alloy, discontinue its use.

Side effects: In individual cases, hypersensitivity reactions (allergies) and electrochemically induced local dysaesthesia have been reported, such as changes in taste and irritation of the oral mucosa. Interactions with other dental alloys: Galvanic effects can occur under proximal or antagonistic contact with dentures of different alloys. If lasting, electrochemically induced, local dysaesthesia arises from contact with other alloys, the restorations must be replaced with other materials.

14. Notes

CoCr base alloys can contain less than 0.1 mass percent trace elements, especially nickel.

Please consult the material safety data sheets or national regulations for disposal.

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Heraenium® = registered trademark of Kulzer GmbH
Kulzer GmbH
Kulzer, LLC