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hydroC Activator
Presentación normal **02245**

Tubo, 60 ml

Precision impression material, silicone based, condensation curing

Fig. 1

1. Tray preparation

Select a tray type suited for desired application. For optimal adhesion of impression material, we recommend use of **sili Adhesive liquid** or **spray**. Apply a thin layer of adhesive film onto the tray and let dry for 60 seconds (Fig. 1). Take preliminary impression with **hydroC putty**.

Fig. 2

2. Preliminary impression

Fill measuring scoop evenly to top with **hydroC putty** (Fig. 2). Smooth **hydroC putty** onto mixing pad with a spatula and make a series of crosswise grooves (Fig. 3). Add recommended amount of **hydroC Activator** (Fig. 4). Mix thoroughly, first using spatula and then by kneading between fingers until a uniform colouring is achieved (approx. 30 sec.). To facilitate the flowing out of excessive light bodied material from the impression, cutting of escape vents and removal of interdental septi in the impression with the use of the special instrument **Deta-Cut** is recommended. Clean the surface of the preliminary impression prior to the corrective impression.


Fig. 3

3. Corrective impression

Extrude **hydroC wash** and **hydroC Activator** in equal lengths and thickness onto mixing pad (Fig. 5). Mix both components for approx. 30 sec. with spatula until a homogeneous uniform colouring is achieved (Fig. 6). Load preliminary impression. A normal impression syringe or disposable syringe can be used for bubble free application of **hydroC wash** (Fig. 7).

Fig. 4

4. Disinfecting

After removing tray from mouth, rinse impression under lukewarm water. Then impression can be disinfected by immersing into 2% glutaraldehyde for 15 min.


Fig. 5

5. Model casting

Casting between 30 minutes and 72 hours after removal of the impression. All class III and IV dental plasters and standard model-casting acrylics can be used.


Fig. 6

6. Electro plating

The impression can be electro plated using all commercially available copper or silver plating baths.

Fig. 7

7. Cleaning the tray

After setting, material can be dislodged with a blunt instrument. Remaining adhesive film can easily be removed with **sili Solvent** spray (FCHC-free; use only in well ventilated rooms). Tray can then be cleaned and disinfected in the usual manner.

Important working hints


- We recommend using **hydroC wash** for corrective impression after making a preliminary impression with **hydroC putty**.
- Overdosing of Activator and increased temperatures accelerate, underdosing and decreased temperatures retard the setting time.
- Activator tubes should be closed tightly immediately after use.
- Cured impression materials are chemically inert – spots on clothing should be avoided.
- Avoid contact of Activator with eyes, may cause irritation. If necessary, flush eyes with plenty of water and seek medical advice immediately.
- Don't leave any residual material in the patient's mouth.
- Please follow the instructions of the safety data sheet!

Avoid contact with eyes and skin!
Wear suitable eye protection and gloves!


Other information:
Silicone impression materials are proven a million times over, unwanted effects can not be expected in proper use. However, immune reactions, e.g. allergies or irritations can not be excluded basically. In the case of doubt, we recommend performing an allergy test before application.

For use by trained specialists.

hydroC putty / hydroC Activator

hydroC wash / hydroC Activator

Indications for use:

- **hydroC putty**
Preliminary impression material for two-stage impressions, putty for double mix technique, situation impressions, check bite
- **hydroC wash**
Corrective material for double impression and double mix technique, relin impressions

Technical data:

hydroC putty
DIN EN ISO 4823 – Type 0

hydroC wash
DIN EN ISO 4823 – Type 3

■ **Dosage:**
hydroC putty / hydroC Activator

1 scoop = 13.5 g /

8 cm = 0.33 g

hydroC wash / hydroC Activator

10 cm = 4.5 g /

10 cm = 0.42 g

■ **Colour code:**
hydroC putty:

pastel green

hydroC wash:

rosa-opaque

hydroC Activator:

blue

■ **Mixing time:**

approx. 30 sec., each

■ **Working time:**

approx. 1 min. 15 sec.*

each

■ **Setting time:**

hydroC putty:

approx. 4 min. 45 sec.*

hydroC wash:

approx. 4 min.*

■ **Time in the mouth:**

hydroC putty:

approx. 3 min. 30 sec.

hydroC wash:

approx. 2 min. 45 sec.

■ **Strain in compression:**

hydroC putty: 3.9 %

hydroC wash: 9.3 %

■ **Recovery from deformation:**

hydroC putty: 98 %

hydroC wash: 98.7 %

■ **Linear dimensional change:**

hydroC putty: 0.44 %

hydroC wash: 0.95 %

■ **Application:**

At 23 °C ± 2 °C /

73 °F ± 4 °F, 50 ± 5 %

relative humidity

■ **Storage:**

15 °C / 59 °F

25 °C / 77 °F

* from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity. Increased temperatures accelerate, decreased temperatures retard a.m. times.

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