Icon vestibular

Inhalt | Contents

DE – Gebrauchsinformation	4
EN – Instructions for use	9
FR – Mode d'emploi	14
IT – Istruzioni d'uso	20
ES – Instrucción de uso	25
PT – Instruções de uso	31
NL – Gebruiksaanwijzing	36
DA – Brugsanvisning	42
SV – Användarinstruktioner	47
PL – Instrukcja użycia	52
RU – Информация по применению	58
TR – Kullanım kılavuzu	64
RO – Instrucțiuni de utilizare	69
CS – Návod k použití	74
EL – Οδηγίες χρήσης	79
ZH – 使用说明	85
JP – 取扱説明書	89
KO – 사용 설명서	94

Produktbeschreibung

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English

Product description

Icon smooth surface is an innovative product for the micro-invasive treatment of lesions. It was especially developed for smooth surface infiltration and is particularly suitable for the treatment and masking of white spots.

Intended purpose

- Micro-invasive treatment of non-cavitated enamel lesions.
- Aesthetic treatment of mild to moderate fluorosis with an intact surface.
- Aesthetic treatment of idiopathic or trauma induced spots.

Limitations of use

Sufficient drying and protection of gingiva is not possible

Indications

- Non-cavitated cariogenic lesions
- Fluorosis
- Idiopathic or trauma induced spots

Contraindications

Do not use the material where there are known allergies to any of the components or contact allergies.

Patient target group

People treated in the course of a dental procedure.

Intended users

Dentist

Notes for use

- Light-curing units should have an output of 450 nm and should be checked regularly.
- The light intensity should be a minimum of 800 mW/cm²
- Place the light unit as close to the material as possible.
- The etched surface must not be touched or contaminated with saliva or moisture until the treatment resumes. If contamination occurs after drying re-etch for approx. 10 seconds and repeat drying with Icon-Dry.
- For best treatment results a sufficiently dry working area is critical. Therefore, appropriate measures for

complete isolation must be taken. A liquid rubber dam can also be used for this purpose. Do not use rubber dams made from thermoplastic elastomers such as Flexi-Dam (Coltène Whaledent / Hygenic).

- Discard syringes after use.
- The syringes included in the treatment unit contain sufficient material for treating two to three smooth surface lesions.
- Multiple smooth surface lesions can be treated at the same time.
- For topical white spots an area of 2 mm beyond the lesion site is etched. At the discretion of the dentist the entire smooth surface area ought to be etched and infiltrated in cases of large area white spots which occur after bracket removal.
- The overall aesthetic result may be improved by means of a whitening treatment before infiltration, especially in the case of fluorosis treatments. The manufacturer's instructions must be taken into account.

Recommended use

1. Preparation

- Before the start of treatment, clean the affected tooth and adjacent teeth. Remove all cleaning residue with water spray.


- Ensure that the working area is sufficiently dry. Therefore, appropriate measures such as a rubber dam or liquid dam must be taken. Make sure that the lesion to be treated remains fully accessible.

2. Uncovering the lesion body


- Pretreatment measures: To ensure successful infiltration, the complete hypermineralized surface layer must be removed prior to applying Icon-Infiltrant in order to gain access. Similarly, access to deeper white spots must also be gained first.
- Optional use of suitable auxiliaries (e.g., medium grit abrasive disk, microabrasion or sandblasting) prior to etching.

Application of Icon-Etch:

- Screw the Smooth Surface-Tip onto the Icon-Etch syringe.

- Apply an ample amount of Icon-Etch onto the lesion site by turning the syringe shaft carefully, let set for 2 minutes and activate it by moving. Remove excess material with cotton roll.
- Suction off Icon-Etch and rinse with water for at least 30 seconds. Then dry using oil-free and water-free air.

3. Result preview and drying

For best treatment results it is necessary to dry the lesion to its full depth.

At this stage a preview of the final result is shown.

Application of Icon-Dry:

- Screw the application tip onto the Icon-Dry syringe, apply an ample amount of material onto the lesion and allow to set for 30 seconds. Dry carefully with oil-free and water-free air.


- When wetted with Icon-Dry, the whitish-opaque coloration on the etched enamel should diminish. If this is not the case, repeat step 2.
- After the visual check dry the lesion with oil-free and water-free air.

4. Infiltration

Do not apply Icon-Infiltrant under direct operating light as this may cause the material to set prematurely.

Application of Icon-Infiltrant

- Screw a new Smooth Surface-Tip onto the Icon-Infiltrant syringe.
- Apply an ample amount of Icon-Infiltrant by turning the shaft.

- Allow Icon-Infiltrant to infiltrate for 3 minutes, occasionally activate by moving and top up if necessary.

In case of deeper and larger defects the aesthetic result can be improved by extending the exposure time.

- Remove excess material with a cotton roll and, if necessary, dental floss.
- Light-cure Icon-Infiltrant for a minimum of 40 seconds from all sides.
- Repeat the previously described treatment and allow penetration for 1 minute
- Remove the rubber dam.
- Use polishing cups, or similar, for the surface finish.

Interactions

Do not use together with substances which inhibit polymerization. This is especially valid for materials which contain eugenol, such as temporary cements.

Residual risks/side effects

- In some cases, contact allergies with similar composite products have been reported.
- Contact of Icon-Etch with the oral mucosa causes a white surface coloration. This will subside after a few days.

Warnings / precautions

- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact immediately wash the affected area thoroughly with soap and water.
- Avoid eye contact! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

Caution! Icon-Etch can cause chemical burns

- Avoid contact with soft oral tissue, eyes, and skin. In case of inadvertent contact rinse with copious amounts of water, and, if necessary, consult a physician.
- Limit the application to the surface to be treated. Avoid contact of the etching gel with gingiva, exposed dentin, root element and adjacent teeth.
- Use only the enclosed syringe and application tip to apply the etching gel in the patient's mouth.

- Discard the applicator tip after use. Replace syringe cap.

Performance characteristics of the device

Penetration coefficient: approx. 147 [cm/s]

Storage / disposal

- Store in a dry place at 2 to 25°C/36 to 77°F!
- Protect against exposure to direct sunlight!
- Do not use after the expiration date!

Composition

- Icon-Etch: Hydrochloric acid (approx. 47 %), water (approx. 33 %), glycerol (approx. 15 %), pyrogenic silicic acid, surface-active substances
- Icon-Dry: Ethanol (approx. 99%)
- Icon-Infiltrant: Methacrylate-based resin matrix (approx. 99%), initiators, additives

Packaging

Starter Pack	
2 treatment units, each including: 1 Syringe @ 0.45 ml Icon-Etch 1 Syringe @ 0.45 ml Icon-Dry 1 Syringe @ 0.45 ml Icon-Infiltrant 6 Vestibular-Tips 1 Luer-Lock-Tip	REF 220343

Economy pack	
7 treatment units, each including: 1 Syringe @ 0.45 ml Icon-Etch 1 Syringe @ 0.45 ml Icon-Dry 1 Syringe @ 0.45 ml Icon-Infiltrant 6 Vestibular-Tips 1 Luer-Lock-Tip	REF 220238

Icon-Etch Refill	
3 Syringes @ 0.45 ml Icon-Etch 15 Vestibular-Tips	REF 220384

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DMG Chemisch-Pharmazeutische Fabrik GmbH