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CE 0297

Prothesenkunststoff Methylmethacrylat -cadmiumfrei-
Typ II Klasse 1 gem. ISO 1567,
gem. Medizinproduktegesetz Kl. 2a CE

Indikation:

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Working Instructions:

Indication:

For dentures with the pressure and injection technique.

Material:

Hot polymerising material based on methylmeth-
acrylate. Powder and liquid available.
For short-time and normal polymerisation.

Characteristics:

High mechanical properties, excellent stability of colour,
very good flowing ability, natural shading.
Important

Preparation:

Invest the wax modulation in half-hard plaster.
Thoroughly scald flask to make it clean and brush
alginat isolating medium Hedent Inkolant onto the
plaster surfaces.

Mixing ratio:

10 g powder: 4 ml liquid.
Free dosage also possible.

Mixing:

Place required quantity of liquid into mixing bowl and stir
in powder, until no more powder can be absorbed. Then
mix thoroughly with a wide spatula. After 10 minutes the
dough is workable for up to 45 minutes. for a longer
working time the mixing ratio is 3 vol. parts powder : 1
Vol. liquid. Here with the dough is - at normal room
temperature - workable for 2-3 hours. The high
mechanical properties are not changed by this.

Pressing:

Place acrylic dough with excess in the flask which is
cooled down to approx. 23 °C. Slowly close flask. k keep
pressure after closing for 10 minutes. A trial pressing is
not necessary due to the good flowing capability of
Inkodon.

Polymerization:

Inkodon was developed especially for the short time
polymerisation. Place flask into boiling water (per flask
approx. 3l) and Polymerise for 20 min. Put thick pieces
into boiling water turn off heat source and wait for
additional 15 minutes. Thereafter get the water to boil
and polymerise for 10 minutes. Of course all other
polymerising techniques are possible with Inkodon, too.
Cooling down as usual.

Cooling down:

Cool down the flask in room temperature.

**Use only original Hedent Inkodon powder and
liquid. Proceed only with that instruction. In case of
different handling the polymerisation cannot be
sufficient and the pateint could notice irritation of
skin and taste. Also the mechanical values can be
lost. To trace back, always the LOT number has to
be noted for each work. (European law of medicine
products)**

Caution Perform liquid monomer contains:

Methylmethacrylate. Highly flammable. Harmful by
inhalation and in contact with skin and if swallowed.
Keep bottle always closed. Keep bottle in a well
ventilated place. Keep away from children and also
keep away from sources of ignition. No smoking. Bottle
has to be closed always. After contact with skin, wash
immediately with plenty of water and soap. Do not apply
into drains. Do not throw to refuse if bottles are not
empty. Pay attention for the rules of toxic waste
Storage temp max 25 °C 77 °F Flashpoint 10 °C, 40 °F

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Hedent GmbH,