# Inkotherm-press 0297

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### Inkotherm press:

GB

### Cold polymerizing material for model casting prosthesis and relinings. Stability of colour

### Working Instruction:

Inkotherm - press  
Universal acrylic resin with normal working time for the model casting technique.

### Indication:

For completion of model casting prosthesis for relinings borders and repairs of dentures.

### Material:

Cold polymerizing material based on methylmethacrylate. Available powder and liquid. Colour stability due to a catalyst system without tertiary amines.

### Characteristics:

Inkotherm-press is colour stable and simple in application. The high mechanical properties correspond to those of a good hot curing material. Natural shading, safe and quick curing.

### Important

Use only original Hedent Inkotherm powder and liquid. Proceed only with that instruction. In case of different handling the polymerisation cannot be sufficient and the patient could notice irritation of skin and taste. Also the mechanical values can be lost. To trace back, always the LOT number has to be noted for each work. (European law of medicine products).

### Preparation:

For Model casting dentures fix the teeth which are standing in wax, by a plaster or silicone precast of Inkoform. The precasts, the model casting frame and the teeth are to carefully cleaned by scalding off the wax remainders. The plaster model must be well watered, so that during pressure polymerizing in the Hedent Inkovac PPA, no air is being pressed from the plaster into the resin which would form bubbles.

### Isolating:

Apply a thin layer of Inkolant-Alginate-Insulation to the model and plaster precasts.

### Mixing ratio:

Mix 10 g Powder: 7ml liquid. Free dosage is possible.

### Working time:

Setting time: 1/2 minute  
Flowing: 2-3 minutes, there after rigid.  
Plastic modelling phase: 3 Minutes.  
Time indicated for a room temperature of 23° C. Higher temperatures shorten, lower temperatures extend the time indicated

### Mixing:

Place liquid into mixing bowl and quickly stir in the powder until a creamy consistency, capable of dripping, is reached.

### Processing:

The mixed resin is to be filled into the precasts during the flowing phase. At the beginning of the plastic phase the resin is rigid, does not longer flow out of the precast and can be modelled.

### Polymerization:

Polymerization takes place after max. 8 minutes in a pressure polymerizing unit - Hedent inkovac PPA - for 6-8 minutes at a pressure of 2.5 – 6 bar and a water temperature of 45-55 °C.

### Important:

The alginate insulation must be well dried in order that no white discolouring occurs at the basis areas. The application of Hedent activator guarantees that no white discolourings occur.

### Packages:

500 g Powder -	250 ml Liquid
1000 g Powder -	500 ml Liquid
5000 g Powder -	1000 ml Liquid

**Use only original Hedent Inkotherm press powder and liquid. Proceed only with that instruction. In case of different handling the polymerisation cannot be sufficient and the patient could notice irritation of skin and taste. Also the mechanical values can be lost. To trace back, always the LOT number has to be noted for each work. (European law of medicine products)**

### Caution Perform liquid monomer contains:

Methylmethacrylate. Highly flammable. Harmful by inhalation and in contact with skin and if swallowed. Keep bottle always closed. Keep bottle in a well ventilated place. Keep away from children and also keep away from sources of ignition. No smoking. Bottle has to be closed always. After contact with skin, wash immediately with plenty of water and soap. Do not apply into drains. Do not throw to refuse if bottles are not empty. Pay attention for the rules of toxic waste Storage temp max 25°C 77°F Flashpoint 10 °C, 40 °F

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**Hedent GmbH - Germany**