Intensiv UniglossCellbrush

DE

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GB

Clinically tested and scientifically approved by the University of Zurich, Switzerland

Product description

- Unique cellulose brush hyper filled with ultra-fine diamond particles for high gloss polishing of aesthetic restorations.
- The cellulose filaments change from hard to ultra soft with water absorption.
- Unique cylindrical shape, with more than 600 extra-long filaments.
- Through water absorption, the softened filaments release the ultra-fine diamond particles.
- The free ultra-fine diamond particles, adapt to the surface structure and significantly increase polishing performance compared to currently available results with other standard brushes.
- The brush has a low-vibration and is atraumatic for the soft tissues.
- Stable polymeric mandrel for conrtriangle pieces.
- Autoclavable instrument before first use.
- Single use recommended (due to water and debris absorption, possible deformation).

Indications

- Final high-gloss polishing of all aesthetic restorations.
- Maintenance of existing restorations.
- For concave and convex morphologies.

Practical Tips

- Autoclave before first use.
- Use Intensiv UniglossCellbrush after pre-polishing with rotating diamond instruments (15µm yellow).
- Use UniglossCellbrush on the dental surface starting with hard filaments for 15 seconds under light pressure (2 N) and low speed (5000 RPM).
- Subsequently add water spray for few seconds.
- The filaments stiffness gradually changes to ultra soft.
- Adjust spray if there is too much water.
- Continue polishing for about 15 seconds more without pressure only with wiping movements to achieve final high gloss.
- Apply Intensiv UniglossCellbrush without paste.
- Single use application recommended (due to water and debris absorption, deformation).
- Intensiv UniglossCellbrush can be safely disposed, as it is made entirely of recyclable materials.

Kavo: Speed calculation example

RPM	"Green ring" reduction 5,4:1 KAVO INTRACompact 0767 LHC	"Blue ring" 1:1 KAVO INTRACompact 2068 LHC
Max.	100% (7'400 RPM)	100% (40'000 RPM)
5'000	ca. 67%	ca. 12%

Sterilization

- Sterilize before the first use.
- Sterilization must be done according to a validated method. If possible, use pre- or post-vacuum autoclave Class B with vacuum drying.
- For sterilisation processes refer to the ISO 17664 standard; the following procedures are recommended:
 - Cycles of 134°C:
 - Tmin. = 134°C Tmax. = 138°C
 - Pressure = 3.15 bar (abs)
 - Time = 4 minutes
 - Cycles of 121°C:
 - Tmin. = 121°C Tmax. = 125°C
 - Pressure = 2.10 bar (abs)
 - Time = 16 minutes

Risk Warnings

- Single patient use is recommended due to water and debris absorption and deformation.
- Additional sterilization cycles may also change the polymeric structure and increase the risk of fracture of the shaft, or increase the risk of detachment of filaments.
- Temperatures above 140°C should be avoided because otherwise the plastic shaft tends to fracture during subsequent use.
- Avoid contact pressure of over 2N, because of fracture risk of the plastic shaft.
- To ensure traceability of the instruments during their entire application, it is recommended to keep the packaging.

Ref. 2200/6 (box of 6 pcs.)
 Ref. 2200/12 (box of 12 pcs.)
 Ref. 2200/18 (box of 18 pcs.)

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