Luxatemp

Automix Plus

Inhalt | Contents

Gebrauchsinformation - DEUTSCH	4
Instructions for use - ENGLISH	11
Mode d'emploi - FRANÇAIS	17
Istruzioni d'uso - ITALIANO	24
Instrucciones de empleo - ESPAÑOL	31
Instruções de uso - PORTUGUÊS	38
Gebruiksaanformatie - NEDERLANDS	45
Brugsanvisning - DANSK	52
Användarinstruktioner - SVENSKA	59
Instrukcja użycia - POLSKI	66
Информация по применению - РУССКИЙ	73
使用说明书 - 中文	81
사용 정보 - 독일어	87
使用説明書 - 日本語	94

Produktbeschreibung

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Product description

Luxatemp -Automix Plus is a self-curing composite for the fabrication of temporary crowns and bridges, inlays, onlays and veneers. The 2-component material is based on multifunctional methacrylates and is automatically mixed.

Indications

- Fabrication of temporary crowns, bridges, inlays, onlays, partial crowns and veneers.
- Fabrication of long-term temporaries.

Contraindications

Do not use the material in cases of allergies to any of the ingredients or, as can occur in rare cases, in the event of contact allergies.

Basic safety instructions

- Only for professional dental use!
- Keep away from children!

- Avoid contact with skin and eyes. In the event of inadvertent contact, rinse thoroughly and immediately with plenty of water and consult a physician if necessary.
- Wear a face mask and goggles when finishing the temporary, otherwise sanding dust may get into the eyes and respiratory tract.

Application system

- Automix cartridge: see “Using the Automix cartridge” on page 14

Note: When using a cartridge for the first time, express a pea-sized amount of material and discard it.

Timing

0:00–0:45 min	Insertion in the mouth
2:00 - 3:00 min	Removal from the mouth
6:00 - 7:00 min	End of setting Finishing of the temporary

Note: The times provided apply to a room temperature of 23 °C and a normal relative air humidity of 50%. Higher temperatures shorten, lower temperatures increase these times.

Recommended use

1. Before preparation of the crown or bridge, or a planned tooth extraction, take an impression using alginate, silicone or a thermoplastic impression material.

Alternatively, the temporary can be made using a laboratory-prepared thermoforming tray or, in the case of single crowns, with a strip crown. Strictly follow the corresponding manufacturer's instructions.

Note: In a silicone impression, undercuts should be blocked out and, if necessary, pressure relief channels cut in. To create a connection with the abutment teeth in the case of a gap in the molar region, cut a channel into the impression.

2. Dry the prepared teeth and block out any undercuts in the preparation. Then lightly

lubricate the prepared teeth and surrounding tissue plus any composite core build-ups using Vaseline or a similar separating medium.

3. Select the desired Luxatemp-Automix Plus colour and place the corresponding cartridge in the Automix-Dispenser (see "Using the Automix cartridge" on page 14).

Note: To avoid bubbles, it is important that the end of the mixing tip always remains immersed in the material during application.

4. Fill the mould or thermoforming tray with Luxatemp-Automix Plus. First, apply the material to the occlusal surfaces of the mould or thermoforming tray and fill up to the gingival areas, only slightly overfilling in the process.
5. Insert the impression or thermoforming matrix, filled with Luxatemp-Automix Plus, over the prepared teeth within 45 seconds after start of mixing, using moderate pressure and fix in place.

Note: The setting reaction has to be monitored intra-orally (e.g. using a scaler), because the cur-

ing is significantly affected by the temperature in the mouth and the temporary can only be removed during its elastic phase.

6. The temporary must be removed from the mouth during the elastic phase, i.e. in the period from 2:00 to 3:00 minutes after start of mixing.

Note: To further improve the precision of large temporaries, they can be put back into the mould after removal from the mouth until they reach final hardness.

7. Before finishing, remove the oxygen inhibition layer from the surface of the temporary (e.g. with alcohol).

Caution! When finishing the temporary, sanding dust may get into the eyes and respiratory tract.

- ▶ Wear a face mask and goggles.
- ▶ Do not inhale the sanded dust and use an evacuator.

8. After about 6:00 to 7:00 minutes from the start of mixing, the temporary can be finished using cross-toothed grinders, flexible disks, etc.
9. Check occlusion and grind if necessary using suitable instruments.
10. If necessary, the surface of the temporary may be glazed with a light-curing high-gloss varnish (e.g. Luxatemp-Glaze & Bond) or polished manually with a rubber polisher or goat's hair wheel.
11. Clean the prepared tooth, dry it and thoroughly remove the separating medium.
12. Cement the temporary.

Note: Any temporary cement can be used for cementing the temporary. DMG recommends using automatically mixing luting cements, e.g. TempoCemNE (non-eugenol zinc oxide cement).

Cements containing eugenol may inhibit the polymerization of resin-based luting cements for the permanent cementation.

13. Check occlusion and adjust accordingly.

Customizing and repair of temporaries

Note: When customizing and repairing a freshly-made temporary, it does not need to be roughened and a bonding agent is not required (proceed with step 2).

Caution! When working on the temporary, sanding dust may get into the eyes and respiratory tract.


- ▶ Wear a face mask and goggles.
- ▶ Do not inhale the sanded dust and use an evacuator.

1. Repair of an older temporary: Roughen the temporary slightly using a grinder or a sandblaster at the appropriate points. Then apply a bonding agent (e.g. Luxatemp-Glaze & Bond) according to manufacturer's instructions.
2. Apply Luxatemp-Automix Plus to all surfaces.
3. Press the parts together for 3 minutes.
4. After 6:00 to 7:00 minutes, the appropriate points can be finished as usual.

Note: To speed up the polymerization, the finished temporary can be placed in warm water at 50 °C (120 °F) for a few minutes.

Alternatively, the temporary can also be customized or repaired using the light-curing, flowable composite LuxaFlow Star. Follow the accompanying instructions for use!

Using the Automix cartridge

Insert cartridge

1. Push up the lever [A] at the back of the Automix-Dispenser and pull back the slide [B] completely.
2. Lift up the plastic lever [C], insert the cartridge and lock it by pushing down the plastic lever [C].

Note: Make sure that the notches on the cartridge line up with those on the Automix-Dispenser.

Attach mixing tip

1. Remove the cartridge cap [D] or used mixing tip by turning it 90° counterclockwise, and discard it.
2. Attach a new mixing tip [E].

Note: Make sure that the notches on the mixing tip line up with those on the cartridge.

3. Turn the mixing tip clockwise 90° to lock it into place.

Apply material

- ▶ When using a cartridge for the first time, express a pea-sized amount of material and discard it.
- ▶ Mix the material in the mixing tip by pulling the trigger [F] on the Automix-Dispenser and then apply directly.

Note: Leave the used mixing tip on the cartridge after use. It serves as a cap.

Remove cartridge

1. Push up the lever [A] at the back of the Automix-Dispenser and pull back the slide [B] completely.
2. Lift up the plastic lever [C] and remove the cartridge.

Disinfection of the Automix-Dispenser

The Automix-Dispenser is reusable after disinfection, although it should be replaced if visible signs of wear appear.

Disinfection can be carried out using a commercially available dipping disinfectant. DMG

recommends using only disinfecting solutions listed by the RKI (Robert Koch Institute). The Automix-Dispenser can also be autoclaved.

Side effects

There are no known systemic side effects to date.

Composition

Glass filler in a matrix of multifunctional methacrylates; catalysts, stabilizers and additives. Free of methyl methacrylate and peroxides. Total filler volume: 44 w. % = 24 vol. % (0.02 to 2.5 μm)

Storage and shelf life

- ▶ Store in a dry place at room temperature (15 to 25 °C/59 to 77 °F)!
- ▶ Do not store in the refrigerator!
- ▶ Protect against exposure to direct sunlight!
- ▶ Leave the used mixing tip on the Automix cartridge after use. It serves as a cap.
- ▶ Do not use after the expiration date!

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DMG Chemisch-Pharmazeutische Fabrik GmbH