# Micro

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## DIN EN ISO 15912, Typ 2, Klasse 1

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## DIN EN ISO 15912, type 2, class 1

### Description:

Micro is a phosphate-bonded, graphite-free precision investment for use in CrCo and crown and bridge work. It is recommended for non-precious alloys.

### Casting models:

Smear Vaseline on the inside of the overlapping section of the duplicating sleeve (REF 101708) to prevent cracking. Measure the required amount of investment and the correct concentration of expansion liquid according to the mixing ratio for Casting models. Mix manually with a spatula before vacuum mixing for 30 - 40 seconds. Use increased vibration for pouring the contours. An instrument is useful for this. Switch off the vibrator and pour the rest of the model. When duplicating models with silicone, it is unnecessary to harden and soak the models. When duplicating models with gel, reduce the powder : liquid ratio to 100 g : 13 - 14 ml. The liquid concentration for the model and muffle should be reduced to 50 %.

### Technical Data:

Mixing liquid:	SILADENT type 100 expansion liquid
Mixing ratio, Powder : liquid:	100 g : 16 ml
Mixing under vacuum:	30 - 40 sec.
Working time:	3 - 5 min.
Working temperature:	approx. 22 °C (Room temperature)
Initial setting time:	6 - 8 min.
Setting expansion:	approx. 0.75% (Type 100 / 75 %)

### Use of pressure:

Investing under pressure is not necessary due to Premium's extremely homogenous grain distribution. If investing under pressure is preferred, the pressure should be released after 6 minutes so that it does not impede the setting expansion.

### Investing:

Place the crepe sleeve (REF 102601) around the investment model. Measure and mix Micro according to the above instructions. For investing use the same liquid concentration as for the casting model. A small reduction of the liquid concentration may improve the connecting of the model and investing. Pour using minimum vibration until the pattern is completely covered. Switch the vibrator off, fill the ring and do not vibrate again. A setting under pressure is not necessary. The preheating can start after 60 minutes.

### Crown and bridge work:

Micro is suitable for investing without a casting ring. Use SILADENT silicone sleeves. Mix the investment according to the instructions and fill the ring quickly using minimum vibration. Micro is also suitable as a die material for firing porcelain using the direct method.

### Recommended preheating curve:

Temperature	Rate of increase	Holding time
Up to 270°C	2°C per minute	40 minutes
270 - 580°C	5°C per minute	30 minutes
580°C - casting temp.*	9°C per minute	30 minutes

\*Recommendation: CrCo or fixed/removable appliances 920 - 950°C

### Concentration of the expansion liquid:

16 ml Liquid / 100 g powder

#### Framework with clasps:

Liquid ratio type 100 / distilled water: 75 % / 25 %

100 g	180 g	400 g
12 ml / 4 ml	22 ml / 7 ml	48 ml / 16 ml

#### Fixed / removable work:

Liquid ratio type 100 / distilled water: 85 % / 15 %

100 g	180 g	400 g
14 ml / 2 ml	25 ml / 4 ml	54 ml / 10 ml

#### Crown & bridge work, precious alloy:

Liquid ratio type 100 / distilled water: 60 % / 40 %

100 g	180 g	400 g
10 ml / 6 ml	17 ml / 12 ml	38 ml / 26 ml

#### Crown & bridge work, semi precious alloy:

Liquid ratio type 100 / distilled water: 75 % / 25 %

100 g	180 g	400 g
12 ml / 4 ml	22 ml / 7 ml	48 ml / 16 ml

#### Crown & bridge work, non precious alloy CrCo:

Liquid ratio type 100 / distilled water: 100 % / - %

100 g	180 g	400 g
16 ml / - ml	29 ml / - ml	64 ml / - ml

\*\* SILADENT type 140 expansion liquid can be used for higher expansion.

### Firing porcelain on Micro:

Duplicate the die with Adisil® and mix at least 50 g of Micro with a 75% to 90% concentration of expansion liquid (must be tested with porcelain). Heat in a porcelain furnace under vacuum to 650°C and hold for 4 minutes before opening the furnace. The dies do not have to be cooled slowly. Glazing the die gives a smoother surface. Fire according to the porcelain manufacturer's instructions.

### Please note:

Technical recommendations are based on tests and findings from work in our development laboratory and can only be regarded as guidelines. SILADENT products are subjected to strict quality controls. We reserve the right to make technical changes.

### Warning:

**Investments contain quartz and cristobalite. Do not breathe in dust. Danger of lung diseases (silicosis or cancer). Use a suitable dust mask!**

**Further questions:** Contact SILADENT technical services (Tel.: +49 (0) 53 21-37 79 25/26) or our sales representatives.

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