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Fig. 7

## 1. Preparation of the denture

- **Removal of existing relinings** (fig. 1): The old, existing relining has to be removed **completely**. Grind off the denture base approx. 1–2 mm, depending on the situation. Incorporate in the areas of transition a shoulder of approx. 2–3 mm in height with a depth of minimum 1 mm. Round off the borders and clean the denture thoroughly.
- **Degreasing** (fig. 2): Degrease by means of pure medical alcohol and let dry.
- **Applying the adhesive for adhesion** (fig. 3): Apply the adhesive **uniformly and completely** onto the grinded denture areas to which the material has to adhere. Avoid contact with eyes! Wear suitable eye protection! **Let the adhesive dry for approx. 1 minute**. A continuous film of adhesive is essential for satisfactory reliner adhesion. After use close the bottle immediately. Saliva, humidity or monomer liquids prevent the adhesion of **mollosil®** to the denture base.

## 2. Direct chairside relining method

- **Coating the denture:** Prepare denture as described under point 1. Mix homogeneously equal lengths of **mollosil®** base material and catalyst for 30 seconds. Apply bubble-free material with spatula onto the vestibular, lingual and dorsal borders. Then coat **completely** the remaining denture with a layer of minimum 2 mm. **All denture areas in contact with the gingiva have to be coated.**
- **Functional movements** (fig. 4): Seat the coated denture into the patient's mouth and let him close in normal occlusion. Setting is finished after approx. 5 minutes. Then remove the relined denture and trim and polish as described under point 4.

## 3. Indirect method in the laboratory

Invest denture to be relined into a flask with counter. For easier removal of the denture, prior to pouring the counter coat teeth with silicone e.g. **Flexistone® Plus**. Remove denture from the flask and prepare as described under point 1. Coat the model with an usual alginate based insulation. Mix homogeneously equal lengths of **mollosil®** base material and catalyst for 30 seconds. Apply bubble-free material with spatula uniformly onto the denture area to be relined (the layer of **mollosil®** should be approx. **1 mm more** than the base material has been grinded off). Close flask with counter. Press with light pressure for approx. 30 minutes or place flask in a pressure pot at 40–45 °C/104–113 °F for approx. 10 minutes in the water bath. Then deflask the relined denture and trim and polish as described under point 4.

## 4. Final trimming and varnishing

- **Smoothing of uneven spots:** Remove excess material with a sharp scalpel or fine scissors (fig. 5). Then smooth any rough areas or ridges with suitable grinding instruments (e.g., **Molloplast® cutters, pre-polishers** and **grinding sleeves** with approx. 15000 r.p.m. under low pressure (fig. 6).
- **Polish** (fig. 7): **Lustrol** gloss varnish is applied only for smoothing of the trimmed functional borders. Avoid contact with eyes and skin! Wear suitable eye protection and gloves! Open the aluminium cap with a spatula, squeeze both liquids onto a mixing pad in a mixing ratio of 1:1, spatulate to a homogeneous mixture and brush on the mixed liquid. Let dry at room temperature for approx. 5 minutes.

## 5. Recommendations for daily care of relined denture

Patients should be advised to keep dentures clean by brushing several times daily with a soft-bristled toothbrush and neutral soap and rinsing thoroughly. Dentures should not be soaked in cleaning baths for more than 10 minutes.


## Important working hints

- Do not interchange tube closures!
- Absolutely firm bonding is obtained only by applying adhesive to roughened surfaces!
- Adhesive must cover entire area of contact completely!
- For proper denture adhesion, all surfaces of the denture in contact with gingiva must be coated evenly with **mollosil®**!
- Setting times must be strictly observed!
- Please follow the instructions of the safety data sheets!

Further information:

Silicone based materials are proven a million times. On condition of a proper application, **undesired effects** are not to be expected. However, reactions of the immune system like allergies, irritations, cannot be absolutely excluded. In case of doubt, we recommend to get tested the allergenic potential **before the application** of the material.

Only for dental use by qualified personnel.

## Indications

### for use:

- Suited for all polymethyl-methacrylate based synthetics
- Fast and easy chairside preparation of permanently soft relinings; chairside or indirect use in the dental laboratory
- Correction of pressure areas in the upper or lower jaw
- Post-damming of the A-line in case of adhesion problems or with flabby ridges
- To steady the position of upper and lower dentures by improving adhesion
- Elastic covering of sharp-ridged alveolar processes
- Cushioning of implants and preparation of over-dentures for wearing during healing

## Technical data:

DIN EN ISO 10139 – 2 Type A

- **Mixed volume:**  
2 x 50 ml (tubes)  
2 x 30 ml (tubes)
- **Mixing ratio:** 1:1
- **Colour code:**  
base: dark red  
catalyst: clear blue
- **Mixing time:**  
manually approx. 30 sec.
- **Working time:**  
approx. 1 min. 30 sec.\*
- **Setting time:**  
approx. 7 min.\*
- **Time in the mouth:**  
approx. 5 min.
- **Recovery from deformation:** approx. 99.5%
- **Linear dimensional change:** approx. 0.2%
- **Application:**  
at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity
- **Storage:**

\* from initiation of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity. Increased temperatures accelerate, decreased temperatures retard these times.

## Ordering information:

- mollosil®**  
Standard packing **03005**  
base, tube of 50 ml  
catalyst, tube of 50 ml  
adhesive, bottle of 5 ml  
Lustrol gloss varnish,  
bottles of 2 x 6 ml  
2 pipettes
- Starter kit** **02274**  
base, tube of 30 ml  
catalyst, tube of 30 ml  
adhesive, bottle of 5 ml  
Lustrol gloss varnish,  
bottles of 2 x 5 ml  
2 pipettes
- mollosil® adhesive** **03007**  
bottle of 5 ml
- Lustrol gloss varnish** **03008**  
bottles of 2 x 6 ml

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