## Multilayer PMMA Disc

Denture Teeth | Temporary Crown and Bridge

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#### Mode d'emploi

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#### Indications for Use

Dentsply Sirona MultiLayer PMMA Discs are used for the manufacturing of temporary crowns and bridges, as well as full arch and tooth segments for dentures using a CAD/CAM system.

#### Contraindications

- May not be used for patients with hypersensitivity to methacrylate or any of the other ingredients.

#### Warnings

- Possible cross-reactions or interactions of this material with other materials already present in the oral environment must be taken into consideration by the dentist.

#### Precautions

- Keep product dust away from eyes.
- After use, wash your hands.
- Do not swallow the product.
- Do not inhale dust particles during grinding.
- Use local vacuum suction and suitable mouth/face protection during manual machining at the workplace.

#### Adverse Reactions

- If properly processed and used adverse effects, are highly unlikely. However, reactions of the immune system (such as allergies) to the substances contained in the material or localized paraesthesia (such as taste disturbances or irritation of the oral mucosa) cannot be completely ruled out.
- Residual dust during milling operation may cause irritation.

#### Symbols on product labels

	Manufacturer
	Batch Code
	Reorder Number
	Manufacture Date
	Expiration Date
	Consult Instructions for Use
	Caution
	Keep Away from Sunlight
	Medical Device
	U.S. Prescription Only
	CE European Conformity Marking

Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.

#### Technical Specifications

Meets applicable requirements of ISO 22112: 2017[E] Dentistry - Artificial Teeth for Dental Prostheses and ISO 10477: 2018[E] Dentistry - Polymer-based Crown and Veneering Materials.

#### Requirements for Milling Denture Teeth

	Recommended Option	Permissible Options
Machine	Dentsply Sirona inLab MC X5	Any 5-axis milling machine with PMMA milling capability (e.g. imes-icore® CORITEC 350i® or DGSHAPE DWX-52D®)
CAM Software	Dentsply Sirona inLab CAM Software	Any CAM software with a milling strategy for PMMA arches or bridges
Material Selection	For Dentsply Sirona inLab CAM machine: Composite	For other equivalent machines: PMMA
Files for Import	.stl(s) of tooth arch(es) or tooth segments as designed by Dentsply Sirona inEos CAD software	.stl(s) of tooth arch(es) or tooth segments as designed by equivalent CAD software (e.g. 3shape® or exocad®)
Import Item or Appliance Type	For Dentsply Sirona inLab CAM: Tooth Arch	For other equivalent CAM software: Full-Contour Bridge, Tooth Arch, Bridge, or a similar selection
Milling Strategy	Premade strategy based upon material and item type selection	Not Applicable
Disc Size	Thickness based upon height of teeth, usually 10mm-20mm	Not Applicable
Milling Tools	Selected as required by milling strategy	Not Applicable

\*imes-icore, 3shape, exocad, and DGSHAPE are not registered trademarks of Dentsply Sirona.

#### Step-by-Step Instructions for Milling Denture Teeth

**Note:** The following instructions support the Dentsply Sirona inLab MC X5 5-axis milling machine. Only 5-axis milling machines with PMMA milling capability should be used.

##### A. Procedure for Fabrication: Design

- Create Tooth Arch or Tooth Segment design(s) in Dentsply Sirona inEos CAD software or equivalent.

##### B. Procedure for Fabrication: Mill

- Create a new job, selecting MC X5 as the machine.
- Select New Item and import the tooth arch or tooth segment .stl or .stls as item type tooth arch and Material Misc Composite.

**NOTE:** In equivalent software for other 5-axis milling machines, some of these steps may be in a different order or omitted. Refer to the machine manufacturer's instructions or consult the reseller for additional information.

- Select imported items and move forward to Select Blank.

- Select Add Blank, selecting Misc and Composite for Manufacturer and Material respectively. Ensure the height matches the disc used for milling. Select this blank and move forward.

- Use Positioning Tools to place the tooth arch(es) or tooth segments in the disc as desired. Occlusal Surfaces should be directed similarly. Place in vertical space according to desired shading from shaded disc. Move forward.

- Support Sprues will autogenerate. Modify as desired, then move forward.

- Under Analyzing Tools, ensure that Production Simulation is selected. The tooth arch(es) or tooth segments should appear mostly green. The recommended Detail Level is Very High. The recommended Machining Mode is Normal.

- Select Check Instruments then ensure that the proper toolset is listed in the software and is placed properly in the machine.

- Select Start Production. Place the MultiLayer PMMA Disc in the machine with the lighter side facing up and secure. Press the Play button on the machine, or select Start Production on the computer.

- After milling is complete, remove the disc from the machine and clean (vacuum or compressed air is sufficient), remove the tooth arch(es) or tooth segments from the disc, and move to bonding.

**NOTE:** In equivalent software for other 5-axis milling machines, a pre-installed milling strategy will be used based upon the material and item type selections. Tools should be used according to the pre-installed milling strategies.

**NOTE:** In equivalent software for other 5-axis milling machines, some of these steps may be in a different order or be omitted. Refer to the machine manufacturer's instructions or consult reseller for additional information.

##### C. Procedure for Fabrication: Bond (Printed & Milled)

#### Printed Denture Base Procedure - Lucitone Digital Print

**NOTE:** Use of personal protective equipment is recommended.

- Plug in the electric hot plate and set the temperature to 90° C.

**Note:** If the electric hot plate is not digital, set to 100°C. Allow the electric hot plate to preheat for at least ten (10) minutes.
- Dry-fit the tooth arch(es) or tooth segments. Care should be taken to avoid chipping the gingival margins.
- Remove the foam insert from the metal container.
- Use tweezers to place the tooth arch(es) or tooth segments into the metal container so that the teeth are positioned with the ridge lap down. Shake the bottle of Lucitone Digital Fuse™ Step 1 - 3D Tooth Conditioning Agent to ensure contents are thoroughly mixed. Add enough liquid to ensure full submersion of the ridge lap of the tooth arch(es) or tooth segments. Immediately place the lid on the metal container.

**CAUTION:** Using an inadequate amount of Fuse Step 1 will compromise tooth bond strength.
**CAUTION:** Keep the metal container free of any solid buildup to avoid compromising tooth bond strength.
- Place the tooth-filled metal container inside the metal ring of the electric hot plate until the 40° C rectangle on the temperature strip turns green (normally in four (4) to six (6) minutes).

**CAUTION:** Tooth bond strength will be compromised if 40° C is not reached or reached too fast.

- Carefully remove the metal container from the electric hot plate. Remove the lid from the metal container and use tweezers to retrieve the tooth arch(es) or tooth segments. Allow tooth arch(es) or tooth segments to air dry for minimum of two (2) minutes. Do not wipe tooth arch(es) or tooth segments.

- Place the lid back on the metal container immediately. Do not return the metal container to the electric hot plate until ready to treat the next set of tooth arch(es) or tooth segments.

- Tooth Arch Technique:** Apply a small amount of Lucitone Digital Fuse™ Step 2 – 3D Denture Bonding Resin shade Original into all tooth pockets. Insert the tooth arch into pockets and apply pressure to the occlusal surface. Hold in place while removing excess resin material from the buccal side using 4x4 gauze. Proceed to step 9.

**Tooth Segments Technique:** Apply a small amount of Lucitone Digital Fuse Step 2 – 3D Denture Bonding Resin shade Original into appropriate tooth pockets for the tooth segment. Be careful not to place resin in adjacent pockets. Insert the teeth into the tooth pockets and apply pressure to the occlusal surface. Hold the tooth segment in place while removing excess resin material from the buccal side using 4x4 gauze. Proceed to step 9. Note: Steps 8, 9 and 10 should be repeated for each tooth segment.

**Note:** Lucitone Digital Fuse Step 2 shade Original may be used to Fuse teeth to all Lucitone Digital Print base shades.

- While continuing to hold the tooth arch or tooth segment in place, use the UV light to tack-cure the buccal side of two (2) teeth for ten (10) seconds. Repeat this step to tack-cure all remaining teeth.
- Remove the excess resin from the lingual side using 4x4 gauze, repeating this around the arch on the buccal of the remaining teeth. While continuing to hold teeth in place, tack-cure for 10 seconds. Repeat this step to tack-cure all remaining teeth.
- Use a Flocked Applicator Tip to apply two (2) thin layers of the Lucitone Digital Fuse™ Step 3 - 3D Denture Sealer around the margins and the interproximal areas of the teeth where Lucitone Digital Fuse Step 2 – 3D Denture Bonding Resin was applied, repeating this around the arch on the lingual of the remaining teeth.

The denture is now ready for post-cure. Refer to the Dentsply Sirona Lucitone Digital Print Denture Instructions for Use section D for step-by-step cure instructions.

#### Milled Denture Base Procedure

**NOTE:** Use of personal protective equipment is recommended.

- Dry-fit the tooth arch or tooth segments. Care should be taken to avoid chipping the gingival margins.
- Using an inverted cone or round bur (no larger than ¼" diameter), adjust the pocket where undercut exist.

**Tip:** The color map from the CAM software and the digital preview of the design will guide the technician to the adjustment area. Reference the digital preview of the design for proper tooth placement.
- Prepare bonding of the tooth arch or tooth segments using Lucitone® HIPA. Measure one small scoop (11 grams) of Lucitone HIPA acrylic powder and (8ml) of Lucitone HIPA Liquid; add powder to liquid and mix for fifteen (15) seconds. Refer to Lucitone HIPA Instructions for Use for additional information.
- Fill all pockets with Lucitone HIPA acrylic, for tooth arch or tooth segments, to the top of the gingival margins.
- Place the tooth arch or tooth segments in order of the dry-fit. Tooth arch or tooth segments will displace excess bonding acrylic from the pockets. Make sure to press the tooth arch or tooth segments into the pockets before the material swells to a point where there is not a distinct "pop" when seating the tooth arch or tooth segments. Note: Lucitone HIPA has a working time of three (3) minutes.
- Clean as much excess acrylic as possible from the tooth arch or tooth segments and base surfaces using a cold cure monomer (Lucitone HIPA monomer should not be used for this step).
- Submerge the denture in a pressure pot to cure the acrylic per the Lucitone HIPA Instructions for Use: 30 mins / 113 F (45° C) / 20 psi.
- Minor contouring may be necessary. Finish and polish the denture using conventional techniques.


Requirements for Milling Temporary Crowns & Bridges
> Minimum thickness: 1mm occlusion and 0.6mm cervical
> Minimum connection diameters: Anterior: 3.6mm; Posterior: 4.4mm
> The provisional can be fixed with all cements for temporary fixation.

#### Step-by-Step Instructions for Milling Temporary Crowns & Bridges

##### A. Procedure for Fabrication

- Place the disc on the machine following the equipment manufacturer's instructions.

**Note:** Shade variation depends on crown/bridge placement within the disc - see figures 1, 2 and 3.

- After milling, remove the disc from the equipment and remove the sprues.
- Finish and polish using conventional techniques.

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