### Mentions de danger

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#### Intended use

Dental sealing lacquer for coating resin temporaries.

#### Indications for use

Varnish for final coating for resin temporaries to produce a smooth and shiny surface.

#### Processing

The area to be lacquered is to be readily treated, fat-free and dry. Due to the thin layer NanoVarnish lacquer is not suitable as polish replacement. The areas to be coated are to be smoothed in advance. For a homogeneous and even coating 10 MicroBrushes accompany the packaging. These ensure an optimum distribution onto the resin surface, even in case of small quantities of lacquer. **Attention:** The sensitive photo-initiator of NanoVarnish lacquer can even respond to the table lamp. The lacquer is to be applied rapidly and in one step and polymerised in-between, if necessary. Only take out the requested quantity and close the bottle immediately.

#### Polymerisation

The lacquer can be cured in all current light-curing machines (spectral area 300–400 nm). The curing time is approx. 2–5 min. (e. g. EyeVolution®), however, it can vary according to the performance of the different machines. All coated areas have to be polymerised completely. Leaching monomers may provoke reactions on the mucous membrane.

#### Important advice

- Do not inhale and do not get in contact with skin or eyes.
- Due to the legal product liability we cannot quote binding polymerisation times for units of other manufacturers.
- Sensitive persons can show allergic reactions to some product components.
- The product is exclusively designed for the a.m. field of application and is only to be used by qualified and instructed persons.
- The general information of this product is based on detailed application practice and intensive research work. Due to further development we reserve the right of modification.

#### Keep away from children!

**Only for use in the dental technique!**

#### Storage conditions

At room temperature, i. e. 18–28 °C. Keep bottle closed well, protect from light.

**Durability:** see date of expiry.

NanoVarnish is not to be used when expiry date has passed.

#### Contains

methyl methacrylate; 2-Propenoic acid, reaction product with Pentaerythrite; Diphenyl(2,4,6-trimethylbenzoyl) phosphinoxide, Acrylic resin.

#### Hazard statements

H225 Highly flammable liquid and vapour. H315 Causes skin irritation. H318 Causes serious eye damage. H317 May cause an allergic skin reaction. H361f Suspected of damaging fertility. H335 May cause respiratory irritation. H412 Harmful to aquatic life with long lasting effects.

#### Precautionary statements

P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water / shower. P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P310 Immediately call a POISON CENTER / doctor. P405 Store locked up. P501 Dispose of contents / container in accordance with local / regional / national / international regulations.

#### Delivery form

**D38900**

#### NanoVarnish

20 ml bottle with drop insert

This product was developed for dental use and has to be applied according to the instructions. In the case of damages resulting from ulterior or improper use, the manufacturer assumes no liability. Furthermore the user is bound to check the product according to suitability and possible use in self dependence before using it subject to its purpose, especially in case the purpose is not mentioned in the instructions.

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#### Composition

metacrilato di metile; 2-Propenoic acid, reaction product with Pentaerythrite; Diphenyl(2,4,6-trimethylbenzoyl) phosphinoxide, Résine acrylique.