cleaning. Although it has effect to apply a little alcohol on the cloth, excessive amount, by seeping inside the Angle Piece, could take off the lubricate oil on ball-bearing. In that case, lubricate again.

3. Inside Cleaning:

Take the same steps with the process of spraying lubricate oil. Operate trial run at low speed every after lubrication. In order to drive dirt out of nose, repeat the process at least twice. Do the cleaning once a day.

4. Cleaning Internal Gears:

Grasp the knurled part of the head cap with the wrench provided and turn wrench clockwise. Remove the stainless steel head gear. Examine the oil hole to make sure it is not clogged. Clean all parts. Reassemble by first inserting the head gear. Attach cap and turn cap counter clockwise with your fingers checking for smooth rotation. Using wrench provided, tighten cap turning counter clockwise.

Sterilization

Do not put the Angle Piece through ultrasonic disinfection solution cleaning procedure. Do not use dry heat.

1. Clean with brush and wipe off by cloth. Refer to the section, Maintenance, 2. Surface Cleaning.
2. Lubricate with spray. Refer to the section, Maintenance, 1. Lubrication.
3. Pack every instrument part into an adequate and fitting transparent, sterile pack (sterilisation bag), packed and sealed in accordance with EN 868-5. Strictly adhere to the instructions given by the manufacturer of the sterilisation bag and the manufacturer of the sealing machine, following the normative requirements
4. Autoclave at the temperature 134°C (2Mpa) for 15 min.
5. Take the Angle Piece out of the autoclave and dry it at room temperature.

Do not operate autoclave dry cycle at the temperature higher than 134 °C!

Warning

1. Do not drop the Angle Piece nor apply excessive force or impact on the bur connected. Warranty is invalid for the trouble caused by those acts. (Prevent the bur on Angle Piece from hitting others in time of autoclave.)
2. When water gets inside the Angle Piece, lubricate more than three times.

The above mentioned instructions are validated as SUITABLE for the preparation of medical products prior to re-use. The reprocessing employees ensure that the actually performed reprocessing and the utilization of equipment, materials and personnel achieve the desired results. Validation and routine supervision are required for this purpose. Likewise, each deviation from provided instructions performed by the reprocessing employees should be carefully examined and evaluated regarding efficacy and possible negative outcome.

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Caution: The product must be cleaned, disinfected and sterilised before first use as well as after any further use, according to the enclosed instructions.

ORBATEC® PROPHY ANGLE PIECE WST. PS

REF 452 153

Instructions for Use

This product is designed for dental use. Its operation and maintenance should be done under dentist's supervision!

Specification:

Gear Ratio:	1 to 1 of motor rotation (Run motor by forward motion only.)
Applicable Brushes/Cups:	SCREW IN TYPE
Coupling:	ISO standard (E-type)
Sterilization:	134 °C autoclave

1. Coupling with Motor:

In order to couple with the motor, push the motor nose into the sheath of the Angle Piece all the way until the stopper of the motor is properly hooked in the inside groove of the Angle Piece. Verify proper coupling by pulling the CA Hand Piece and the motor.

2. Operation:

Use with running below 5,000 rpm. Run motor by forward rotation only.

Maintenance:

To keep the Angle Piece under good condition, lubricate oil specifically made for dental hand piece should be used. Use the spray as its manufacturer. Lubrication should be done every day.

1. Lubrication

- a. Remove cup or brush
- b. Thoroughly clean the outside of the angle with a disinfectant moistened 2" X 2" gauze sponge.
- c. Be sure to lubricate after autoclave, and in the case that water and/or dust gets inside. Spray lubrication could also make inside cleaning. It should be done twice each time. Spray lubricate oil at the back of sheath of the Angle Piece. (It has the same effect with spray lubrication that air, with the hose attached, is supplied after 2 drop of lubricate oil is placed at the back of sheath.)
Without being held firm by hand, the Angle Piece could be pushed away by air-pressure of spray. Accidental drop of the equipment could cause damage on it.
- d. After oiling, rotate the shaft to circulate the oil with the head. Excess oil will be expelled from the oil hole. Wipe off excess.
- e. Sterilize by autoclaving before putting back into service.

2. Surface Cleaning:

Clean with brush and wipe off by cloth. Lubrication should be done every time after surface