**DE****GEBRAUCHSANWEISUNG****PARASORB® Cone Genta / PARASORB® Fleece Genta HD**

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GA 0066, Stand 2020-05

## PARASORB® Cone Genta / PARASORB® Fleece Genta HD

EN Please read instructions carefully prior to use.

The information given in this package leaflet is updated regularly. Please read instructions carefully prior to use.

Collagen cone/collagen sponge with antibiotic protection for surgical use - haemostyptic, sterile -

### COMPOSITION

1 cone PARASORB® Cone Genta contains:

- Collagen from equine tendons, 22.4 mg
  - Gentamicin sulfate 16 mg, corresponding to 8.8-11.44 mg gentamicin
- 1 sponge PARASORB® Fleece Genta HD sized 2.5 x 2.5 x 0.5 cm contains:
- Collagen from equine tendons, 35 mg
  - Gentamicin sulfate 25 mg, corresponding to 13.82-17.88 mg gentamicin

### APPLICATION FORM

Absorbable implant with antibiotic protection, local haemostyptic agent

### INDICATIONS

- PARASORB® Cone Genta/PARASORB® Fleece Genta HD is inserted for local haemostasis and as a defect filling into clean, clean contaminated and contaminated wound cavities.
- PARASORB® Cone Genta/PARASORB® Fleece Genta HD is used for haemostasis following extractions or maxillary surgical intervention with a high risk of infection.
- PARASORB® Cone Genta/PARASORB® Fleece Genta HD is also used for reducing alveolar ridge atrophy following tooth extractions.

### MODE OF APPLICATION AND DOSAGE

After removing the focus of inflammation PARASORB® Cone Genta/PARASORB® Fleece Genta HD is inserted into the defect cavity. PARASORB® Cone Genta/PARASORB® Fleece Genta HD is loosely placed or if there is a bone defect mixed with tiny bone fragments and put in place. PARASORB® Cone Genta/PARASORB® Fleece Genta HD is applied in a dry state by exerting slight pressure in order to achieve better adhesion. The sponge/cone can be cut to the desired size using sterile scissors/scalpel. Unless otherwise directed, PARASORB® Cone Genta/PARASORB® Fleece Genta HD is inserted directly into the extraction alveolus. PARASORB® Cone Genta is adapted to the extraction alveolus, brings about haemostasis and is absorbed.

### PREGNANCY AND LACTATION

PARASORB® Cone Genta/PARASORB® Fleece Genta HD must not be used during pregnancy and lactation. No data on the safe use of PARASORB® Cone Genta/PARASORB® Fleece Genta HD during pregnancy and lactation are available. Aminoglycoside antibiotics cross the placenta and reach the fetus. There is therefore a risk of intrauterine fetal damage. Since aminoglycoside antibiotics also pass into breast milk, PARASORB® Cone Genta/PARASORB® Fleece Genta HD should also not be used during lactation, or breastfeeding should be ceased.

### EFFECTS ON ABILITY TO DRIVE AND OPERATE MACHINES

None known as yet.

### SIDE EFFECTS

Under some circumstances at the start of tissue resorption by the PARASORB® Cone Genta/PARASORB® Fleece Genta HD, increased secretion may occur. With simultaneous systemic administration of aminoglycoside antibiotics or in patients with limited renal function an increase in nephrotoxicity and ototoxicity may occur. Foreign body reaction e.g. formation of seroma and granuloma, may occur.

### OVERDOSE

The recommended application quantity of PARASORB® Cone Genta/PARASORB® Fleece Genta HD causes a negligible increase in serum gentamicin levels. Regarding dosage over and above this there are no known findings available. With a deviation from this dosage in each case the serum gentamicin levels and kidney function must be monitored. In cases of severe intoxication, a peritoneal dialysis or haemodialysis may be considered.

### INCOMPATIBILITIES

None known as yet.

### OTHER NOTES

Because of the affinity of collagen for bloody surfaces PARASORB® Cone Genta/PARASORB® Fleece Genta HD should only be applied using dry gloves and instruments. The product is intended for single use only. PARASORB® Cone Genta/PARASORB® Fleece Genta HD must not be reused, there is a risk of infection for the patient and/or user and the functionality of the product cannot be guaranteed. Remaining material must be disposed of according to local procedures and guidelines and must not be used again, as traceability and sterility are no longer guaranteed.

Haemostasis takes place by a physical method. Caused by the surface properties of collagen the contact between collagen and blood leads to adhesion of platelets and triggers platelet aggregation. The collagen matrix ensures additional strengthening of the coagulum and thus inhibits connective tissue healing of the extraction alveolus. The release of the incorporated gentamicin occurs simultaneously with the dissolving of the collagen cone/collagen sponge. This protects the implant from external contamination as well as from contamination by bacteria spread or introduced during the surgical procedure. In cases of diffuse bleeding (e.g. in the sacral cavity) it may prove useful to apply the implant by pressing it onto the affected areas of the wound. Overpacking the wound with collagen may result in the formation of a seroma. PARASORB® Cone Genta/PARASORB® Fleece Genta HD is able, because of its structure, to absorb large quantities of fluids. Through this purely mechanical process, by means of secretion absorption, shed materials such as bacteria and fibrin layers are absorbed. This accelerates the formation of granulation tissue.

The dosage of PARASORB® Cone Genta/PARASORB® Fleece Genta HD is governed by the size of the defect and the weight of the patient. Dependent upon this, the maximum permissible quantity of gentamicin sulfate for patients up to 50 kg body weight is 600 mg gentamicin sulfate or for patients over 50 kg body weight 1000 mg gentamicin sulfate. PARASORB® Cone Genta/PARASORB® Fleece Genta HD is completely absorbable. PARASORB® Cone Genta/PARASORB® Fleece Genta HD may only be used by qualified medical professionals.

### CONTRAINDICATIONS

PARASORB® Cone Genta/PARASORB® Fleece Genta HD must not be used with a known hypersensitivity to collagen and/or gentamicin or other aminoglycoside antibiotics. PARASORB® Cone Genta/PARASORB® Fleece Genta HD should only be used as an adjunct to other procedures for haemostasis in the presence of bleeding requiring control by ligatures or in instances of larger arterial and/or venous haemorrhage requiring suture ligatures.

### WARNINGS

PARASORB® Cone Genta/PARASORB® Fleece Genta HD should only be used under strict indication guidelines in patients with limited renal function also patients with autoimmune diseases or neuromuscular diseases such as Parkinson's or myasthenia gravis. Simultaneous systemic administration of aminoglycoside antibiotics should be avoided or only given if there is strict monitoring of serum-gentamicin levels and renal function.

### NOTES FOR POST-OPERATIVE CARE FOLLOWING DENTAL SURGERY

To be avoided:

- Rinse your mouth out on the day of surgery (may dislodge clot).
- Drink alcohol or hot drinks on the day of surgery (may increase blood pressure or induce hyperaemia).
- Smoke for at least 24 hours (causes vasoconstriction).
- Carry out physical activity on the day of surgery (increasing blood pressure).

Recommended:

- Sip lukewarm drinks carefully (helps avoid clot removal by washing out socket).
- Bite on gauze or clean, damp handkerchief for 10 minutes if area bleeds, preferably sitting upright in a chair.
- Rinse mouth gently with hot salt mouthwash or Chlorhexidine (0.2%) mouthwash after 24 hours. Continue regular mouthwash use over the following week, especially after meals.
- Clean teeth as normal, avoiding area of surgery for 24 hours.

### SHELF LIFE, STERILITY

The use-by date and batch number are printed on the sales packaging and the blister cover. After the use-by date has expired PARASORB® Cone Genta/PARASORB® Fleece Genta HD must no longer be used. The contents of unused, opened or damaged packaging must not be re-sterilised and must therefore be discarded. PARASORB® Cone Genta MK11 and PARASORB® Fleece Genta HD are packaged in double sterilised packaging. The interior package is sterile inside and outside and therefore can be passed on in a sterile state.

### STORAGE

Store between 0°C - 25°C.

Store protected from moisture and direct sunlight. Do not refrigerate. Store out of the sight and reach of children.

### PACK SIZES

PARASORB® Cone Genta, 10 pcs, size Ø max.: 1,2 cm; H: 1,6 cm, REF MK10  
PARASORB® Cone Genta, 5 pcs, size Ø max.: 1,2 cm; H: 1,6 cm, REF MK11  
PARASORB® Fleece Genta HD (6,25 cm<sup>2</sup>), 5 pcs, size 2,5 x 2,5 cm, REF MK25

PRODUCT CONFORMS TO THE ESSENTIAL REQUIREMENTS OF THE COUNCIL DIRECTIVE 93/42/EEC CONCERNING MEDICAL DEVICES.

Cross-allergies between aminoglycosides exist. There are no records about the use of PARASORB® Cone Genta/PARASORB® Fleece Genta HD in patients with immunological diseases or connective tissue diseases available. Although there has been no evidence found that the use of animal collagen leads to a deterioration in disease PARASORB® Cone Genta/PARASORB® Fleece Genta HD should only be used under strict indication guidelines in these patients. There are no records about its use in paediatrics available.

### SPECIAL PRECAUTIONS FOR USE

To guarantee full functionality of the product, the usual application techniques for a haemostypticum should be observed. For safe application the following instructions for its application must be observed:

- Moistening or wetting the haemostypticum prior to implantation can lead to impairment in its action both in terms of haemostyptic properties and also with respect to a loss in the inherent protective characteristics of PARASORB® Cone Genta/PARASORB® Fleece Genta HD because of a premature release of the highly watersoluble gentamicin sulphate.
- The desired platelet aggregation is also achieved by lining the wound with PARASORB® Cone Genta/PARASORB® Fleece Genta HD.

### INTERACTIONS WITH OTHER SUBSTANCES

Although after implantation of PARASORB® Cone Genta/PARASORB® Fleece Genta HD only very low plasma levels are reached, attention should be paid to the interactions of gentamicin with other drugs. A simultaneous administration of aminoglycosides with loop diuretics such as furosemide or ethacrynic acid must be avoided since loop diuretics themselves have an ototoxic effect. If i.v. diuretics are given at the same time these can increase the concentration of gentamicin in serum and tissue and boost toxicity. Simultaneous or successive systemic administration or topical application of potentially neurotoxic and/or nephrotoxic substances such as Cisplatin, other aminoglycosides, Streptomycin, Cephaloridine, Viomycin, Polymyxin B or Polymyxin E can increase the toxicity of gentamicin. Simultaneous local use of beta-lactam antibiotics can lead to changed action (inactivation). A tendency to neuromuscular blockade can be increased by the simultaneous administration of muscle relaxants, e.g. d-Tubocurarine, Suxamethonium or Pancuronium, as well as ether. A neuromuscular blockade can best be remedied by calcium salt. If a neuromuscular blockade occurs with simultaneous administration of Suxamethonium in the case of acquired or genetic cholinesterase deficit, artificial respiration must be given and cholinesterase administered.

1

2

3

### PREGNANCY AND LACTATION

PARASORB® Cone Genta/PARASORB® Fleece Genta HD must not be used during pregnancy and lactation. No data on the safe use of PARASORB® Cone Genta/PARASORB® Fleece Genta HD during pregnancy and lactation are available. Aminoglycoside antibiotics cross the placenta and reach the fetus. There is therefore a risk of intrauterine fetal damage. Since aminoglycoside antibiotics also pass into breast milk, PARASORB® Cone Genta/PARASORB® Fleece Genta HD should also not be used during lactation, or breastfeeding should be ceased.

### EFFECTS ON ABILITY TO DRIVE AND OPERATE MACHINES

None known as yet.

### SIDE EFFECTS

Under some circumstances at the start of tissue resorption by the PARASORB® Cone Genta/PARASORB® Fleece Genta HD, increased secretion may occur. With simultaneous systemic administration of aminoglycoside antibiotics or in patients with limited renal function an increase in nephrotoxicity and ototoxicity may occur. Foreign body reaction e.g. formation of seroma and granuloma, may occur.

### OVERDOSE

The recommended application quantity of PARASORB® Cone Genta/PARASORB® Fleece Genta HD causes a negligible increase in serum gentamicin levels. Regarding dosage over and above this there are no known findings available. With a deviation from this dosage in each case the serum gentamicin levels and kidney function must be monitored. In cases of severe intoxication, a peritoneal dialysis or haemodialysis may be considered.

### INCOMPATIBILITIES

None known as yet.

### OTHER NOTES

Because of the affinity of collagen for bloody surfaces PARASORB® Cone Genta/PARASORB® Fleece Genta HD should only be applied using dry gloves and instruments. The product is intended for single use only. PARASORB® Cone Genta/PARASORB® Fleece Genta HD must not be reused, there is a risk of infection for the patient and/or user and the functionality of the product cannot be guaranteed. Remaining material must be disposed of according to local procedures and guidelines and must not be used again, as traceability and sterility are no longer guaranteed.

### NOTES FOR POST-OPERATIVE CARE FOLLOWING DENTAL SURGERY

To be avoided:

- Rinse your mouth out on the day of surgery (may dislodge clot).
- Drink alcohol or hot drinks on the day of surgery (may increase blood pressure or induce hyperaemia).
- Smoke for at least 24 hours (causes vasoconstriction).
- Carry out physical activity on the day of surgery (increasing blood pressure).

Recommended:

- Sip lukewarm drinks carefully (helps avoid clot removal by washing out socket).
- Bite on gauze or clean, damp handkerchief for 10 minutes if area bleeds, preferably sitting upright in a chair.
- Rinse mouth gently with hot salt mouthwash or Chlorhexidine (0.2%) mouthwash after 24 hours. Continue regular mouthwash use over the following week, especially after meals.
- Clean teeth as normal, avoiding area of surgery for 24 hours.

### SHELF LIFE, STERILITY

The use-by date and batch number are printed on the sales packaging and the blister cover. After the use-by date has expired PARASORB® Cone Genta/PARASORB® Fleece Genta HD must no longer be used. The contents of unused, opened or damaged packaging must not be re-sterilised and must therefore be discarded. PARASORB® Cone Genta MK11 and PARASORB® Fleece Genta HD are packaged in double sterilised packaging. The interior package is sterile inside and outside and therefore can be passed on in a sterile state.

### STORAGE

Store between 0°C - 25°C.

Store protected from moisture and direct sunlight. Do not refrigerate. Store out of the sight and reach of children.

### PACK SIZES

PARASORB® Cone Genta, 10 pcs, size Ø max.: 1,2 cm; H: 1,6 cm, REF MK10  
PARASORB® Cone Genta, 5 pcs, size Ø max.: 1,2 cm; H: 1,6 cm, REF MK11  
PARASORB® Fleece Genta HD (6,25 cm<sup>2</sup>), 5 pcs, size 2,5 x 2,5 cm, REF MK25

PRODUCT CONFORMS TO THE ESSENTIAL REQUIREMENTS OF THE COUNCIL DIRECTIVE 93/42/EEC CONCERNING MEDICAL DEVICES.

### DESCRIPTION OF THE SYMBOLS USED ON THE PACKAGE

	= REFERENCE NUMBER
	= BATCH NUMBER
	= USE BY YEAR - MONTH - DAY
	= CONSULT INSTRUCTIONS FOR USE
	= DO NOT RE-USE
	= DO NOT RESTERILIZE
	= DO NOT USE IF PACKAGE IS DAMAGED
	= STERILIZED USING ETHYLENE OXIDE
	= TEMPERATURE LIMIT
	= KEEP AWAY FROM DIRECT SUNLIGHT
	= HIBC CODE
	= KEEP DRY
	= MANUFACTURER

CE 1434 = CE MARKING AND IDENTIFICATION NUMBER OF THE NOTIFIED BODY

MANUFACTURER  
RESORBA Medical GmbH, Am Flachmoor 16  
90475 Nürnberg, Germany, Tel. +49 9128 / 91 15 0

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RESORBA Medical GmbH, Am Fleckmoor 16