
# Provil<sup>®</sup> novo **Dynamix Putty**

**KULZER**  
MITSUI CHEMICALS GROUP

(DE)	Gebrauchsanweisung.....	4
(GB)	Instructions for use .....	10
(FR)	Mode d'emploi .....	16
(ES)	Instrucciones de uso .....	22
(IT)	Istruzioni per l'uso .....	28
(PT)	Instruções de utilização .....	34
(NL)	Gebruiksaanwijzing .....	40
(SE)	Bruksanvisning .....	46
(DK)	Brugervejledning .....	52
(NO)	Bruksanvisning .....	58
(FI)	Käyttöohjeet.....	64
(GR)	Οδηγίες χρήσης .....	70
(CZ)	Návod k použití .....	76
(HU)	Használati utasítás .....	82
(LV)	Lietošanas instrukcija.....	88
(LT)	Naudojimo instrukcija.....	94
(PL)	Instrukcja obsługi .....	100

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### **Elastomeric, silicone-based precision impression material, addition-cured**

The Provil novo range is an optimally matched system and its different types of material can be combined for all impression techniques to suit specific requirements.

### **Notes on use**

Not to be used in conjunction with other silicone materials. Visually inspect for damage before use. Damaged products must not be used. With addition-cured silicone impression materials changes in setting characteristics can occur if the material comes into contact with latex gloves, retraction solutions or freshly placed composite fillings. Provil novo impression materials should be at room temperature 23°C (73°F) for processing.

### **Contraindications**

The use of this product is contraindicated in case of known or suspected allergies against components of this product. Provil novo Putty and Putty soft contain peppermint oil.

### **Side effects**

This product or one of its components may in particular cases cause hypersensitive reactions. If suspected, information on the ingredients can be obtained from the manufacturer. Provil novo contains surfactants for hydrophilation, which in some patients may cause a superficial anesthetic effect. Surplus material on the impression tray or in the patient's mouth should be avoided as this may result in gagging, swallowing or aspiration. To further reduce these risks during impression taking, the patient should be positioned as upright as possible.

Provil novo	EN ISO 4823	Product colour		Mixing ratio	Working time*	Min. time in mouth	Accessories	
		Base	Catalyst				Mixing Tips	Intraoral Tips
<b>Dynamix Putty</b>	Type 0 putty	turquoise	light grey	5:1	1:30	2:30	dynamic yellow 66035999	–

\*at room temperature of 23°C (73°F), 50% rel. humidity. Higher temperatures shorten, lower temperatures prolong these times.

### Model fabrication

The impression can be cast 30 minutes after removal from the mouth – there are no other time constraints. The impression is stable in long-term storage.

### Recommended model material

All standard model materials can be used.

We recommend: Class III plaster, MOLDANO®

Class IV plaster, MOLDASTONE®/MOLDASYNT®

## **Disinfection**

Impressions can be treated with an aqueous sodium hypochlorite solution (5.25%, for up to 10 minutes\*) or with other disinfectant solutions suitable for silicones. When doing so, please observe the disinfection/processing instructions issued by the relevant manufacturers.

*\* Gelson Luis Adabo, DDS, PhD, Elaine Zanarotti, Renata Garcia Fonseca, and Carlos Alberto dos Santos Cruz, DDS, PhD Effect of disinfectant agents on dimensional stability of elastomeric impression materials. J Prosthet Dent, 81, 621-4, 1999*

## **Adhesion**

We recommend using our Universal Adhesive.

## **Storage conditions**

Do not store above 25°C (77°F).

Do not use the material after the expiration date.

Please quote batch number and the article number in all correspondence about the product.

According to the EU Medical Devices Regulation, users / patients are obliged to report serious events with a medical device to the manufacturer and to the competent authority of the country in which they occurred.

For disposal please follow the instructions on the material safety data sheet and national regulations.

## Dynamix System, 380 ml – step by step

Fig. 1  
Remove cap, insert  
cartridge into the  
mixer


Fig. 2  
Attach Dynamix  
Mixing Tip


Fig. 3  
Position Dynamix  
fixation ring


Fig. 4.1  
Attach Dynamix  
fixation ring, slide  
over recesses


Fig. 4.2


Fig. 5  
Lock by giving  
a ¼ turn  
clockwise

### Inserting the Dynamix cartridge

For information on how to handle the mixer being used please refer to the detailed instructions issued by the relevant manufacturer.

1. Hold the tab on the cap of the outlets and bend it upwards through 90° without applying any excessive force.
2. Take the tab between your thumb and index finger whilst you hold it securely with your other hand and remove the cap from the cartridge completely (Fig. 1).
3. Place one of the enclosed Dynamix Mixing Tips over the outlets of the cartridge (Fig. 2).
4. If it proves difficult to attach the Dynamix Mixing Tip, check to make sure that the central internal hexagon of the mixing tip is properly aligned with the hexagon of the drive spindle.

5. When the Dynamix Mixing Tip is in the correct position (Fig. 3), push the Dynamix fixation ring over the Dynamix Mixing Tip until it stops and turn it clockwise (1/4 turn) until it is securely locked in place (Figs. 4.1, 4.2 and 5).
6. When impression material has been mixed, transfer it to a suitable impression tray.
7. When the material has been dispensed, retract the mixer plungers from the cartridge in order to avoid continuous pressure and thus prevent excess material from continuing to emerge from the mixing tip.

### **Notes about safe use of the Dynamix System**

Dynamix cartridges are heavy and can be damaged if they are dropped.

Damaged cartridges must no longer be used.

Only use Dynamix cartridges in mixers with a closed protective flap.

Before using for the first time, adjust the levels of the two material components.

Dispense approx. 2–3 cm material and discard this amount before starting to fill the tray.

When the cartridge is used in different mixers, 2–3 cm of material has to be discarded after each device change. Each time before you use the material check the outlet channels of the cartridge for blockages. At the end of a working day or if the mixer is not going to be used for a lengthy period, remove the cartridge from the device and store it in a horizontal position with the outlets closed off by the mixing tip or cap.

**After use, leave the mixing tip on the cartridge to keep it closed. Immediately before the next use, remove the mixing tip and attach a new one. To do so, detach the Dynamix fixation ring and use it again for the new mixer attached. Mixing Tips and Intraoral Tips are single-use products. To avoid cross contamination Mixing Tips and Intraoral Tips must be disposed of immediately after initial use.**

Only use Dynamix System components for dispensing and mixing the impression material in Dynamix cartridges. Check the system components for physical damage on a regular basis and only use components which are intact.

### Warnings and safety instructions

Avoid contact with skin and eyes. May be irritating to eyes and skin. Polyvinyl siloxanes are chemically resistant. Material which has not set can soil clothing. **Wear protective gloves according to EN 420, category III, with protection against chemical risks according to EN 374, type C.** On contact with skin wash off with water and soap. If in eyes rinse with water for several minutes. Contact a physician in case of persistent complaints. Do not swallow or ingest. Do not leave any residual material in the sulcus or oral cavity, as this can lead to abscess formation. If health problems arise after swallowing impression material, seek medical attention immediately. Intestinal blockage may arise in rare cases.

**For use by dental personnel only.**

Dated: 2019-09

#### Description of applicational / packaging icons

- A** = Only use pastes with same LOT No
- B** = Icon "Use of the special program at the Dynamix® *speed* Mixing Unit possible"
- C** = Mixing ratio
- D** = Mixing time

- E** = Extraoral working time
- F** = Intraoral working time
- G** = Time in mouth
- H** = Static mixing tip
- K** = Intraoral tip

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Kulzer GmbH
Kulzer, LLC