- DE Gebrauchsanweisung
- GB Instructions for use
- FR Mode d'emploi
- ES Instrucciones de uso
- IT Istruzioni per l'uso
- PT Instruções de uso
- NL Gebruiksaanwijzing
- SE Bruksanvisning
- DK Brugervejledning
- NO Bruksanvisning
- FI Käyttöohjeet

Pastós
Pasty
Pâteux
Pastoso
Pasteus
Tröglyflande
Pasta-agtig
Pastaliggins
Pastamainen

Manufacturer:
Kulzer GmbH
Leipziger Straße 2
63450 Hanau (Germany)

+49 800 4372522

Made in Germany

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Stand: 2021-05

GB Instructions for use Pala cre-active®

Medical Device - for use by dental health care professionals only.

Intended purpose
This dental material is suitable for patients requiring dental treatment for the following indications with consideration of the contraindications. For pregnant and nursing women, the treatment risks must be weighed carefully against the benefits taking into consideration the unborn child or infant.

KULZER medical devices ensure the rehabilitation of oral functions such as chewing, speaking and aesthetics. They stabilize the remaining dentition and/or the alveolar ridge restoratively or prosthodontically.

Product description
Pala cre-active is a light-curing material for the customised characterisation of methyl methacrylate-based prosthetic resins. It is paste-like and available in a range of colours (white, polar, pink, red, maroon, black, gum, gingiva shade 200, gingiva R50, gingiva light pink, gingiva pink, gingiva clear).

Composition Pala cre-active colour fluids
Multifunctional methacrylic acid esters. Silicane dioxide (silanated), photo-initiators, stabilisers, inorganic pigments (some fluorescent). Pala cre-active colour fluids are cadmium-free.

Pala cre-active gingiva
Multifunctional methacrylic acid esters. Silicane dioxide (silanated), bead polymer, photo-initiators, stabilisers, inorganic pigments.
Pala cre-active gingiva is cadmium-free.

Indications

- Customisation of (methyl methacrylate-based) total and partial prostheses
- Customisation of implant restorations

Contra-indications
In case of known or suspected allergy to components of the product, its use is contraindicated. Do not use this product in case of known or suspected allergies against (meth)acrylate compounds.

Side effects
Hypersensitivities to the product or its components cannot be excluded in individual cases. In case of suspicion please contact the manufacturer for ingredients information.

Advice
Pala cre-active gingiva can be mixed with Pala cre-active colour fluids to create customised hues.

Handling
Prior to the application of Pala cre-active, the surface of the prosthesis is roughened with a coarse diamond bur or via sandblasting (50 µm) and cleaned. Apply the Signum® connector with a brush, allow to set for 2-3 minutes and then polymerise for 90 seconds. The ready-to-use Pala cre-active mixture can then be applied directly to the cured prosthetic resin with a brush or modelling instrument. The different consistencies and hues make it possible to give the different parts of the gums smooth or rough surfaces in pale red or dark red tones.

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© = registreret varemærke tilhørende Kulzer GmbH
© = Kulzer GmbH:n rekisteröity tavaramerkki