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Warning notices: No special notices

Limitation of reprocessing: Frequent reprocessing has only got a small impact on this instrument. The instrument's shelf life is mainly determined by wear and damage by use.

Place of handling: Clean surface with a disposable cloth or paper

Storage and transport: No special requirements.

Preparations for cleaning: Disassemble product into its components (remove adapter)

Cleaning: Only use suitable cleansers, carefully following the instructions of the respective manufacturer. After cleaning, remove cleanser thoroughly under running water. Do not exceed 93 °C for mechanical cleaning.

Disinfection: Only use suitable disinfecting solutions, carefully following the instructions of the respective manufacturer. Remove disinfecting solution thoroughly under running water. Do not exceed 93 °C for mechanical disinfection.

Maintenance: No special requirements.

Packaging: Standardised packaging material may and should be used for sterilisation processes.

Sterilisation: Vapour sterilisation at 134 °C for 5 minutes, vapour sterilisation at 121 °C for 15 minutes

Control/Functional check: Sight check on damages, wear, deformation.

Storage: No special requirements.

The a.m. instructions have been validated as SUITABLE for preparation of a medical device and its re-use by the medical device manufacturer. It is the reprocessor's responsibility that the actual reprocessing achieves the required results with the equipment and materials applied and personnel involved in the reprocessing site. Normally, validation and routine controls are necessary. Furthermore, the reprocessor should evaluate any deviation from the provided instructions regarding efficiency and possible adverse consequences.

Fabricante: HAGER & WERKEN GmbH & Co. KG
Ackerstraße 1, 47269 Duisburg, Germany
T +49 (203) 99269-0, F +49 (203) 299283
www.hagerwerken.de

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Fabbricante: HAGER & WERKEN GmbH & Co. KG