# PremEco® Line Gießkunststoff

Kaltpolymerisat

DENTAL

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# PremEco® Line casting resin

Cold-curing polymer

**(en)** Instructions for use, please read carefully!

**(USA)** Caution: Federal law restricts this device to sale by or on the order of a licensed dental practitioner.

## Indications

Colour-fast, cold-curing, customisable methylmethacrylate acrylic base specially developed for casting partial and full dentures.

## Composition

PremEco® Line casting resin powder contains

- PMMA (polymer and copolymer)
- barbituric acid catalyst system
- organic colorants
- inorganic pigments

PremEco® Line casting resin liquid contains

- MMA (methylmethacrylate)
- dimethacrylate
- barbituric acid catalyst system

## Dosage and application

- See table for mixing ratios and processing times.
- Use class IV plaster for making casts where possible.
- Before duplication always wet casts sufficiently for 5 to 15 minutes depending on the degree of dryness.
- Use standard alginate isolation to separate plaster from acrylic.
- Dosage aids for PremEco® Line casting resin monomer and polymer can be used for dosage of the acrylic - REF 1090194\*.
- Do not compress the powder when using dosage aids.
- As soon as the acrylic surface in the casting holes is satin smooth and dry, the casting holes can be filled with cold water and the flask can be placed in the polymerisation unit. The flask should be three-quarters immersed in the warm water bath.

## Times and quantities

Mixing ratio*	30 g : 18 ml
Polymerisation temperature	45 °C / 113 °F
Polymerisation time	15 min <sup>1</sup> / 30 min <sup>2</sup>
Pressure required during polymerisation	2 - 2,5 bar

Mixing time, stir thoroughly for 30 sec

Dwell time

Processing time (pourable consistency)

\*\*Based on a room temperature of 21 °C / 69,8 °F

<sup>1</sup> in matrix technique, plaster, putty

<sup>2</sup> in flask duplicating technique, duplicating gel, duplicating sili-

PremEco® Line casting resin is a product of the PremEco® Line production system. For more information see [www.merz-dental.de](http://www.merz-dental.de).

### **Processing the acrylic**

To reduce the residual monomer content, the finished (sanded) but not polished polymerised workpiece should be stored in water at room temperature for about 48 hours before delivery.

### **Trimming**

To prevent inaccurate fitting after polymerisation, avoid excessive heat development during trimming and polishing.

### **Bond**

- Roughen the basal surfaces of the acrylic teeth before bonding.
- Ensure that no impurities such as wax or isolating agent are present.
- Mechanical retention is recommended.
- Adequate mechanical retention is essential for ceramic teeth.

### **Restoration**

Products of PremEco® Line casting resin can be restored and supplemented at any time with all standard MMA-based cold-curing polymers, such as Weropress®.

### **Special notes**

- Close all powder and liquid containers tightly after use.
- Store in a dry and cool area away from direct sunlight.
- Take precautionary measures against static discharges.
- Do not use the material after the expiration date.
- Record the LOT number with every process that requires identification of the material.
- The residual monomer content after polymerisation of PremEco® Line casting resin is < 2,5%.

### **Disposal information**

Disposal of contents/container in accordance with local/regional/national/international regulations.

### **Danger and safety instruction PremEco® Line casting resin Monomer**

H225 Highly flammable liquid and vapour.

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H335 May cause respiratory irritation.

P233 Keep container tightly closed.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P272 Contaminated work clothing should not be allowed out of the workplace.

P262 Do not get in eyes, on skin, or on clothing.

P210 Keep away from heat. No Smoking.

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P333+P313 If skin irritation or rash occurs: Get medical advice/attention.

**Danger and safety instruction PremEco® Line casting resin Polymer**

May cause an allergic skin reaction.

Avoid breathing dust/fume/gas/mist/vapours/spray.

**Classification**

PremEco® Line casting resin conforms to EN ISO 20795-1 type 2 class 2.

**The product properties are based on compliance with the instructions for use.**

For detailed information see the applicable Material Safety Data Sheet.

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