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PROVI K+B

CE 0482

Hersteller:R-dental Dentalerzeugnisse GmbH • Winterhuder Weg 88 • 22085 Hamburg
T 040-30707073-0 • F 040-30707073-73 • I www.r-dental.com

REF PRO2010 (Farbe A2)

REF PRO2012 (Farbe A3,5)

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- REF PRO2010 (Shade A2)**
- REF PRO2012 (Shade A3,5)**

Description and indication

PROVI K+B (Automix 4:1) is a self-curing, temporary crown and bridge material based on multifunctional methacrylic esters for manufacturing short - and long-term temporary crown and bridges, inlays, onlays and veneers (temporaries) directly at the patient (chairside).

Characteristics

The very hard and very elastic PROVI K+B is especially suitable for manufacturing longer bridge spans. The high stability is particularly evident on thin crown edges.

In UV-light PROVI K+B shows fluorescence. PROVI K+B is free of *methyl* methacrylate! The intraoral curing temperature is lower than 104 °F (40 °C).

The polymerized temporary shows an increasing transverse strength, a good abrasion resistance and a low polymerization shrinkage. Further on it is characterized by a good aesthetic, an optimal polishability and therefore an extremely low surface roughness, a high colour stability and a low oxygen inhibition layer. The temporary protects the prepared teeth against external influences and preserves the occlusion.

In the case of a fracture PROVI K+B-temporaries are easy to repair, because cured material can be repaired with fresh material or with light-cure flowable composites.

In nearly every property PROVI K+B is superior to the common materials based on monomer/polymer systems on *methyl* meth-acrylate base.

PROVI K+B is delivered in 4:1 auto-mixing cartridges and can be mixed automatically.

Application

1. Situation impression

Take an impression preferably with alginate INDIKATOR-ALGINAT (REF ALG1000) or R-SI-LINE® PUTTY (REF PUT1010) or R-SI-LINE®-impression materials (long storage stability!) before preparing the teeth or extract any tooth. The interdental areas should be cut off so that the temporary crown and bridge material can bridge all teeth. In molar areas with absent teeth it may be necessary to cut a groove in the impression between the abutments to create a bridge-like connection between the tooth units.

Attention: Block out undercuts in the silicone impression and if necessary cut grooves into the impression.

2. Preparing the auto-mixing cartridge/R-SI-LINE® MIXING PISTOL PN

Introduce the auto-mixing cartridge into the MIXING PISTOL PN (REF MIP1086) and lock it. Before first use screw off the original cap and dispose it. When screwing on a new mixing cannula the cartridge should be tested if both pastes will be extruded equally. Then screw on the mixing cannula R-SI-LINE® MIXING CANNULA PN (REF MIK1096). The mixing pistol is ready for application.

The filled MIXING CANNULA functions as a seal until the next application. Do not re-use the original cartridge cap.

3. Application of PROVI K+B into the impression

Before application a small amount of material (size of a pea) should be pressed out and should be wasted immediately. Then apply the extruded PROVI K+B material into the impression (from incisal to gingival areas). The top of the mixing tip should remain in the extruded material during the application so that the filling should occur from bottom upward to prevent voids.

4. Placing the loaded impression onto the prepared areas

Place the loaded impression onto the prepared areas in the mouth or into the model (in the lab for an extraoral application) within 50 s (working time at 73 °F, 23 °C).

Attention: Teeth with core build-ups made of dual-cure composites must be isolated (e. g. with FANTESTIC® OxyBlock, REF FOS3970), otherwise PROVI K+B could be adhere to the core build-up materials.

5. Take the impression out

About 1-2 min (setting time) after placing into the mouth (after 3-4 min extraoral on the model) the impression with the temporary can be taken out. The temporary still shows an elastic condition.

Attention: The temporary can only be removed without destruction during the elastic state.

6. Removal of the temporary

After removal of the temporary from the impression (or from the tooth preparation) excess material and proximal undercuts will be removed.

The temporary can be post-cured in warm water (113-131 °F, 45-55 °C) or in a polymerisation medical device.

If the last step (hot water) is omitted, the temporary can be worked out after 6 min from the beginning of application with rotary instruments and can be polished to high gloss like a composite.

Attention: Do not breathe the polishing dust; use suitable mouth protective device or aspiration!

Attention: The oxygen inhibition layer ("smear layer") on the surface caused by the acrylic system cannot be fully avoided and should be removed before working out and high gloss polishing. It can easily be removed by alcohol or other suitable solvents.

7. Cementing the temporary

The temporary should be cemented preferably with an eugenol-free temporary cement like the self-curing ZAKK® (REF ZKS2070). Do not use eugenol containing cements. Eugenol containing cements will inhibit the hardening of final cements.

8. Repairing a fracture of a temporary shortly after manufacturing

If a freshly manufactured temporary despite a high mechanical stability breaks, the following procedure is recommended: refix both ends with freshly extruded PROVI K+B from the cartridge. Alternatively the light-cure and flowable repair - and correction material ZAKK® C+B Correct (REF ZKS2097 (Shade A2) bzw. REF ZKS2098 (Shade A3,5)) can be applied.

The oxygen inhibition layer ("smear layer") on the surface caused by the acrylic system cannot be fully avoided and should be removed before working out and high gloss polishing. It can easily be removed by alcohol or other suitable solvents.

9. Repairing a fracture of a worn temporary

Clean and roughen the areas of the fracture with a sand blaster and prepare some mechanical retentions. The fractured areas can be bonded with fresh, automatically mixed material from the cartridge. To enhance the curing the temporary can be place for some minutes in warm water.

Alternatively the light-cure and flowable repair - and correction material ZAKK® C+B Correct (REF ZKS2097 (Shade A2) bzw. REF ZKS2098 (Shade A3,5)) can be applied in layers not exceeding 2 mm. Then the fractured parts can be fit and the layer can be light-cured for 40 s with commercial polymerization devices.

The oxygen inhibition layer ("smear layer") on the surface caused by the acrylic system cannot be fully avoided and should be removed before working out and high gloss polishing. It can easily be removed by alcohol or other suitable solvents.

Special storage and/or handling conditions

PROVI K+B is free of *methyl* methacrylate, but contains other acrylates. With susceptible patients, sensitization to PROVI K+B cannot be excluded. PROVI K+B should not be used any more, if allergic reactions are observed. Do not use for patients with allergic reactions against acrylates.


Avoid contact of the PROVI K+B with soft tissues, skin or eyes. If there is an accidental skin contact, flush immediately with soap and numerous amounts of water. In case of contact with eyes, rinse immediately with plenty of water and consult a physician if necessary.

Keep attention on the material safety data sheet (MSDS >> www.r-dental.com).

Expiry date see label. Do not use after expiry date.

Do not store above 78 °F! Avoid exposure to sunlight. Keep out of the reach of children! The cartridge should be stored horizontally.

Time scale for processing and curing of PROVI K+B

R-SI-LINE, FANTESTIC and ZAKK are registered trademarks of R-dental Dentalerzeugnisse GmbH.

Caution: For dental use only.

Warranty

R-dental Dentalerzeugnisse GmbH warrants this product will be free from defects in material and manufacture. R-dental Dentalerzeugnisse GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. The user is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and R-dental Dentalerzeugnisse GmbH's sole obligation shall be repair or replacement of the product.

Limitation of Liability

Except where prohibited by law, R-dental Dentalerzeugnisse GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

© 2020 R-dental Dentalerzeugnisse GmbH. All rights reserved.