R-PILOT™

de	Gebrauchsanweisung	2 – 3
en	Directions for Use.....	4 – 5
fr	Mode d'emploi.....	6 – 7
it	Manuale d'uso.....	8 – 9
es	Manual de usuario	10 – 11
nl	Gebruiksaanwijzing	12 – 13
cs	Návod k použití	14 – 15
pl	Instrukcja użytkowania.....	16 – 17
ro	Instrucțiuni de utilizare	18 – 19
tr	Kullanım Talimatı	20 – 21
el	Οδηγίες χρήσης.....	22 – 23
pt	Manual de instruções.....	24 – 25
ru	Инструкция по применению	26 – 27
zh	使用说明	28
ko	사용 설명서.....	29 – 30
sr	Uputstvo za upotrebu	31 – 32
sv	Bruksanvisning	33 – 34
da	Brugsanvisning	35 – 36
hu	Használati útmutató.....	37 – 38
id	Petunjuk Penggunaan.....	39 – 40
ja	取扱説明書	41 – 42
uk	Інструкція з використання	43 – 44
ar	امامتسسال ليلد	45 – 46

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For dental use only, Rx only

R-PILOT™ instruments for root canal preparation

1. Composition

The working part of the R-PILOT™ instruments is made of nickel-titanium alloy (M-Wire).

2. Contents

R-PILOT™ blister package with 1, 2 or 6 sterile instruments each

3. Indications for Use

R-PILOT™ instrument is a Glide Path instrument used in reciprocating motion to prepare the root canal system before the shaping with a rotary or a reciprocating instrument.

R-PILOT™ instrument may only be used in reciprocating motion with a designated drive system with original VDW RECIPROC® settings, e.g. VDW.CONNECT Drive® or VDW.GOLD® RECIPROC®. The use of reciprocating drive systems without original VDW RECIPROC® settings can lead to instrument fracture and misuse.

4. Contraindications

Similar to all machine driven root canal instruments, R-PILOT™ should not be used in a root canal with an abrupt apical curvature in the apical region which is managed with pre-curved hand files.

5. ⚠ Warnings

None known.

6. ▽ Precautions

- The instruments may only be used by qualified persons in dentists' offices or clinics.
- Use a rubber dam.
- R-PILOT™ instruments are sterile, single use instruments, designed to be used in no more than one molar. The instrument cannot be re-sterilized.
- Examine the instrument during treatment for signs of wear such as untwisting, bending or damaged cutting edges. In these cases, instruments should be discarded.
- Clean the flutes from debris in an interim stand every three pecks and irrigate the canal according to the appropriate irrigation protocol. (INTERIM-STAND [REF] 495)

- Only light pressure should be applied on R-PILOT™ instruments. If resistance is encountered or if the instrument does not advance, clean the instrument, irrigate the canal according to the appropriate irrigation protocol and check the canal for possible blockage with a C-PILOT® file ISO size 08 at minimum or similar hand file. Then use the R-PILOT™ instrument briefly with a lateral brushing motion.

7. Adverse Reactions


None known.

8. Step-by-Step Instructions

The application of R-PILOT™ is recommended for the following cases

- **Initially** to render the existing natural canal path smooth prior to the use of shaping instruments
- **During root canal preparation with RECIPROC®** whenever advancement with RECIPROC® becomes difficult (e.g. RECIPROC® R25 or RECIPROC® *blue* R25)

Symbols used

	Explanation
	Manufacturer
	Packaging unit
	Catalogue number
	Batch code
	Use-by date
	Sterilized using irradiation
	Do not use if packaging is damaged
	Consult instructions for use
	Do not re-use
	CE marking
	Warnings
	Precautions
	RECIPROC®
	M-Wire

Prior to the use of shaping instruments

1. Establish a straight-line coronal and radicular access to the root canal entrance.
2. Use a C-PILOT® file ISO size 08 at minimum or similar hand file to ensure patency and determine working length with the help of a preoperative radiograph [mandatory] and/or an apex locator.
3. Afterwards use R-PILOT™ as described in sections a-f until working length has been reached.
 - a. Place irrigant in the access cavity of the root canal.
 - b. Introduce R-PILOT™ free without binding into the canal. Press the ON/OFF button or foot pedal of the motor when the instrument is at the root canal orifice.
 - c. Move the instrument in a slow in-and-out pecking motion. The amplitude of the in-and-out movements should not exceed 3 mm. Only very light pressure should be applied. It is also possible to brush laterally. The instrument will advance easily in the canal. One in-and-out movement = one peck.
 - d. Clean the instrument in the interim stand after 3 pecks.
 - e. Irrigate the canal.
 - f. Repeat steps c-e until full working length has been reached.
4. Continue your shaping with the VDW file system of your choice.

During root canal preparation with RECIPROC®

In some canals the RECIPROC® instrument may stop from advancing in the canal or its advancement may become difficult.

1. Use a C-PILOT® file ISO size 08 at minimum or similar hand file to determine working length with the help of a preoperative radiograph [mandatory] and/or an apex locator.
2. Afterwards use R-PILOT™ as described in section 8. 3a-f until working length has been reached.

9. Disinfection, cleaning and sterilization and disposal

R-PILOT™ is intended for single-use only. Disinfection, cleaning and sterilization procedures are not recommended. Re-use of R-PILOT™ can increase the risk of cross contamination and breakage. The country-specific and hospital-specific disposal requirements have to be followed.

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