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R-SI-LINE® PUTTYSOFT

Hersteller:

R-dental Dentalerzeugnisse GmbH • Winterhuder Weg 88 • 22085 Hamburg
T 040-30707073-0 • F 040-30707073-73 • I www.r-dental.com

REF PUS1020 (R-SI-LINE® PUTTYSOFT)

REF PUS1020FS (R-SI-LINE® PUTTYSOFT FS (fast set))

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REF PUS1020 (R-SI-LINE® PUTTYSOFT)

REF PUS1020FS (R-SI-LINE® PUTTYSOFT FS (fast set))

REF PUS1020SH (R-SI-LINE® PUTTYSOFT SH (super hydro))

REF PUS1020SHFS (R-SI-LINE® PUTTYSOFT SH FS (super hydro, fast set))

Description and indication

R-SI-LINE® PUTTYSOFT (- FS, - SH, - SH FS) is used as a precision base impression material for the two-phase-impression-technique (one step) in conjunction with R-SI-LINE®-correction impression materials.

Characteristics

R-SI-LINE® PUTTYSOFT (- FS) is a kneadable addition curing impression material on vinylpolysiloxane base. It is hydrophilic and therefore compatible with moist surfaces.

R-SI-LINE® PUTTYSOFT SH (- SH FS) is super hydrophilic and preferred used for impressions in conjunction with the super hydrophilic R-SI-LINE®-correction impression materials (super hydro).

R-SI-LINE® PUTTYSOFT (- FS, - SH, - SH FS) is easy to mix and not sticky.

R-SI-LINE® PUTTYSOFT (- FS, - SH, - SH FS) meets the requirements of ISO 4823 Type 0 (very high consistency - putty).

Recommended procedures

1. Preparation of the tray

It may be used perforated or non-perforated metal or plastic trays. In order to improve the adhesion of R-SI-LINE® PUTTYSOFT (- FS, - SH, - SH FS) to the tray, the R-SI-LINE® SILICONE-ADHESIVE (REF SIA1070) should be applied as a thin even film according to the instructions for use.

2. Dosage

Use equal parts (volume spoon or weight) of base paste and catalyst paste (dark spoon for base paste, white spoon for catalyst paste, do not interchange).

Then knead both pastes at least 30 s until the colour is homogeneous.

Attention: Use Base und Catalyst only with same LOT-No. Do not interchange lids of pots or spoons of base and catalyst!

3. Impression taking

The ready kneaded material is then loaded into the prepared tray. Insert the loaded tray into the mouth within 1 min 30 s (the fast set materials - PUTTYSOFT FS, - SH FS within 60 s) after start of mixing together with the R-SI-LINE®-correction material (working time).

Apply the parallel mixed correction impression material directly on the top of the base material on the tray and around and onto the prepared areas.

Maintain these two-phase impressions with normal set correction impression materials and with PUTTYSOFT (- SH) at least 3 min, and leave the two-phase impressions with the fast set correction impression materials and with PUTTYSOFT FS, - SH FS at least 2 min, in the mouth.

Times mentioned for 74 °F (23 °C). Higher temperatures shorten, lower temperatures prolong these times.

4. Casting and disinfection

The impression should be stored at normal room temperature (not above 77 °F, 25 °C) and can be poured 30 min at the earliest after removal from the mouth. There are no other time limitations. No separating materials are needable.

The impressions can be copper - and silver plated and disinfected with standard disinfection solutions.

Special storage and/or handling conditions

Do not interchange lids of pots or spoons of base and catalyst!

Use Base und Catalyst only with same LOT-No.

Do not use Latex gloves, these may interfere with the setting reaction of the A-silicone impression material. Recommended are gloves based on polyethylene (PE).

Pay attention to the material safety data sheets (MSDS >> www.r-dental.com).

Shelf life see label. Do not use after expiry date.

Do not store above 77 °F (25 °C)! Avoid exposure to sunlight. Keep out of the reach of children!

Do not use R-SI-LINE® PUTTYSOFT (- FS, - SH, - SH FS) with standard condensation reaction type silicone impression materials.

Technical data

Mixing time: at least 30 s

Working time incl. mixing time (- PUTTYSOFT, - SH): < 1 min 30 s

Working time incl. mixing time (- PUTTYSOFT FS, - SH FS): < 60 s

Recommended time in mouth (- PUTTYSOFT, - SH): ≥ 3 min

Recommended time in mouth (- PUTTYSOFT FS, - SH FS): ≥ 2 min

Time (Start mixing - withdrawal of mouth) (- PUTTYSOFT, - SH): ≥ 4 min 30 s

Time (Start mixing - withdrawal of mouth) (- PUTTYSOFT FS, - SH FS): ≥ 3 min

Dimensional change after 24 hours: < 0.1 %

Elastic recovery: > 98.6 %

Strain in compression: about 2.6 %

Shore A-hardness: about 67

R-SI-LINE is a registered trademark of R-dental Dentalerzeugnisse GmbH.

Caution: For dental use only. US Federal Law restricts this device to sale by or on the order of a dentist or licensed professional, to use or to order the use of this device.

Warranty

R-dental Dentalerzeugnisse GmbH warrants this product will be free from defects in material and manufacture. R-dental Dentalerzeugnisse GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. The user is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and R-dental Dentalerzeugnisse GmbH's sole obligation shall be repair or replacement of the product.

Limitation of Liability

Except where prohibited by law, R-dental Dentalerzeugnisse GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

© 2020 R-dental Dentalerzeugnisse GmbH. All rights reserved.

Bite to Perfection

R-SI-LINE® PUTTYSOFT

Fabricant:

R-dental Dentalerzeugnisse GmbH • Winterhuder Weg 88 • 22085 Hamburg • DE
Téléphone +49 40 - 30707073-0 - Fax +49 40 - 30707073-73 - www.r-dental.com

REF PUS1020 (R-SI-LINE® PUTTYSOFT)

REF PUS1020FS (R-SI-LINE® PUTTYSOFT FS (fast set))

REF PUS1020SH (R-SI-LINE® PUTTYSOFT SH (super hydro))

REF PUS1020SHFS (R-SI-LINE® PUTTYSOFT SH FS (super hydro, fast set))

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R-SI-LINE est une marque déposée de R-dental Dentalerzeugnisse GmbH.
© 2020 R-dental Dentalerzeugnisse GmbH. Tous droits réservés.