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
## **1. PRECAUTIONS**

### **1.1 Safety precautions**

**BEFORE starting to set up the unit  
PLEASE NOTE:**

**Work on electrical appliances may ONLY be carried out by qualified technicians (in accordance with EN 50110-1) and/or by our customer service staff.**

**Use only original parts when carrying out repairs, otherwise the warranty will be rendered invalid.**

**Please observe the corresponding safety and inspection regulations, both when combining electrical appliances and when repairing their electrical components:**

- **VDE 0100, IEC 364, VDE 0105, 0107, 0751, IEC 601-1 UVV VBG 1, 4, 103**
- **Accident-prevention regulations UVV-VBG 1, 4, 5, 103.**
- **Directive regulations regarding medical products- Directive 93/42/EEC**
- **The „MPG“ (German Law on Medical Products), along with corresponding national legislation.**

### **1.2 Electrical Safety Inspection**

In accordance with IEC 601-1 and VDE 0751 an electrical safety inspection is required whenever an electrical appliance is installed, put into operation, subject to maintenance work, serviced or repaired. Written evidence should be kept of the test results.

### **1.3 Warning signs and symbols**

The following indications and/or signs are used in the operating instructions to mark especially important points:

**Information and/or instructions or prohibitions regarding personal safety or extensive material damage.**

specific data relating to the efficient operation of the machine and other information

### **1.4 Intended use**

The DÜRR Regenerating Unit is designed to be operated exclusively in conjunction with the DÜRR XR25 X-ray film developer. Any use other than or in excess of specifications does not constitute intended use. The manufacturer accepts no liability for damage or injury resulting from such use. All such use is entirely at the owner's risk.

The definition of „intended use“ also covers proper use according to the operating instructions and observance of the installation, operating and maintenance conditions.

## 2. DELIVERY CONTENTS

### RU25 Regenerating Unit

230V, 50Hz/60Hz .....	1700-840-00
110V, 50Hz/60Hz .....	1700-850-00

### 2.1 Accessories

2 Locking rings .....	9000-311-22
1 Operating instructions .....	9000-600-23/01

## 3. POSSIBLE INSTALLATION CONFIGURATIONS FOR THE REGENERATING UNIT

- **Installation of the RU25 Regenerating Unit alongside the XR25 Developer**

This configuration should be selected whenever sufficient space is available.

- **Installation of the RU25 Regenerating Unit in a cabinet below main unit:**

Refer to instructions on regular check of conditions of chemicals.

- **Installation of the RU25 Regenerating Unit above the XR25 is not permitted,**

if installation were carried out in this position, chemicals might run into the developer of their own accord.

Please ensure when installing that there are no kinks in the tube connecting the RU25 Regenerating Unit to the XR25 Developer. The distance between the Regenerating Unit and Developer should not exceed 2 metres.

Ensure when installing that sufficient access room is left for the Developer and containers.

## 4. TECHNICAL DATA

### Dimensions (cm)

Height .....	37
Width .....	36
Depth .....	32

### Weight when empty (kg) .....

4

### Voltage (V)

1700-840-00 .....	230
1700-850-00 .....	110

### Frequency (Hz)

1700-840-00 .....	50 / 60
1700-850-00 .....	50 / 60


### Amperage (A)

1700-840-00 .....	0.7 / 0.5
1700-850-00 .....	0.9 / 0.7

### Power consumption (W)

1700-840-00 .....	75 - 100
1700-850-00 .....	40 - 65

## 5. FUNCTION DIAGRAM

## 6. FUNCTION DESCRIPTION

Each time the XR25 Developer is connected at the main switch, the regenerating pumps in the XR25 start up automatically and run for 30 seconds (factory-adjusted).

The regenerating pump running time setting is non-stepped and made via a potentiometer with a time function element.

In 30 seconds, the pumps carry approx. 190 ml of fresh chemicals from the 5-litre regenerating containers to the developing and fixing tank of the XR25 Developer.

If the Developer is left switched to standby, the regenerating pumps automatically run for 30 seconds every 24 hours (factory-adjustment).

This timer period remains programmed until the power supply is switched off. When the unit is reconnected, the regenerating pumps run for 30 seconds and the timer begins to count again from zero.

If the **unit is normally switched off during the lunch break**, we recommend adjusting the **regenerating cycle period to 15 seconds** (see section 10.2).


The non-return valves installed on top of the pump are to prevent chemicals flowing back.

**A** XR25 Developer

**B** RU25 Regenerating Unit

**B1** Time-switch built into the RU25 Regenerating Unit

**B2** Non-return valves to prevent back-flow of chemicals



## 7. PREPARATION TASKS

- Drain chemicals from tanks of XR25 Developer.
- Rinse out tanks with warm water.
- Close water shutoff tap.
- Pull out mains plug.
- Remove cover panels from XR25 unit.


## 8. CONNECTION TO THE XR25 DEVELOPER

D = Developing tank  
F = Fixing tank


- Push the connecting tubes for developing solution (black) and fixing solution (red) through the holes (1) at the back of the XR25 and fit them onto the connecting nozzles (3) of the tanks.

- Unscrew the union nuts (2) from the connection nozzles (3) and pull out the black drain plug.
- Push the clamping nut and grey clamping ring (4) onto the regenerating tube. Fit the connecting tubes for developing solution (black) and fixing solution (red) to the tanks by pushing them as far as they will go onto the corresponding connection nozzles and securing them with union nuts.


- Remove side cover panels from the XR25 unit.

## 9. ELECTRICAL CONNECTION

- Fit the supplied electrical connection plug (6) into the opening provided at the back of the XR25 and secure it using 4 B2.9 x 9.5 screws.

- Run the wires from the electrical connection plug (6) along the wire conduits (5) to the terminal strip (7) under the drive motor, on the right-hand side of the XR25. Connect as follows:
- Attach PE wire, (earth - green/yellow), to the free terminal 
- Attach N wire (blue) to the free terminal marked **N**.
- Attach L1 wire (black) to the free terminal marked **1**.

## 10. PUTTING INTO OPERATION

- Fill the tanks of the XR25 Developer with water up to the maximum level mark.
- **Check flow rate of regenerating pumps:** Fill two measuring jars with water (approx. ½ litre), place them in the Regenerating Unit and insert the suction tubes in the measuring jars.

When inspecting the **Regenerating Unit** the pump flow rates should be checked. The developer and fixing solution pumps should pump **identical** quantities into tanks E and F. If there is any difference between the flow rates of the pumps, the piston stroke can be adjusted accordingly, see section 10.1.

- **Checking the regeneration period:** Turn main switch of XR25 Developer to ON. The regenerating pump will run for approx. 30 seconds (factory-adjusted).


The regeneration period can be individually adjusted, e.g. if the main switch of the establishment is disconnected during the lunch break. See section 10.2.

**Never run the unit without an D/E-Roller set fitted.** If the unit is operated without an D/E-Roller set the liquids cannot be diverted to the tanks and spray out upwards.

- **Danger of eye injuries !**
- **Danger of chemicals mixing**

- Check the unit for water-tightness.
- Once the examination is over, drain the water from the regenerating containers and the tanks of the unit and then add chemicals, referring to the instructions for the XR25 (9000-600-03/..).

## 10.1 Adjusting the piston stroke of the regenerating pump

If the feed rates for developer and fixing solution are not identical, they can be made to match by adjusting the stroke of both pistons

- Loosen Allen screw (9).
- Adjust piston stroke by turning the hexagonal bolt (8) with a spanner or socket wrench: Indicator towards 100%=increased feed rate  
Indicator towards 0%=decreased feed rate
- Retighten the Allen screw.

## 10.2 Adjusting the regenerating period $t_1$

$t_1$  Regenerating period (pump running time)

Regenerating period adjustment range pre-set at factory to 30 seconds.

$t_2$  Pump downtime

Pump downtime adjustment range, pre-set at factory to 24 hours. **Do not alter !**

### 10 Time-switch

Each time the main switch of the XR25 unit is turned on the time switch, in its present position, transmits a regenerator pulse. **Do not alter.**

11 Regenerating period adjuster – **coarse** adjustment

12 Regenerating period adjuster – **fine** adjustment

Each time the XR25 Developer is connected at the main switch, the regenerating pumps start automatically and run for **30 seconds** (factory-adjusted).


If the main switch of the establishment is disconnected during the lunch break, we recommend adjusting the regeneration period to **15 seconds**. This is due to the fact that when the unit is switched on every morning and afternoon, about 100 ml each of developing and fixing solution are pumped into the tanks of the XR25 unit.

**Example:** Desired regenerating period: 15 sec.

- Set adjuster **11** to 1min. (60 sec.)
- Set adjuster **12** to 2.5 (=0.25)
- 0.25 x 60 sec. (1min) = 15 sec.

**or**

- Desired regenerating period: **30 sec.**
- Set adjuster **11** to 1min. (60 sec.)
- Set adjuster **12** to 5 (=0.5)
- 0.5 x 60 sec. (1min) = 30 sec.

## 11. MAINTENANCE

Drain chemicals from Developer every three months. Fill the **unit and the regenerating container** with warm water and shut the lid.


**Never run the unit without an D/E-Roller set fitted.** If the unit is operated without an D/E-Roller set the liquids cannot be diverted to the tanks and spray out upwards.

- **Danger of eye injuries !**
- **Danger of chemicals mixing**

- Switch on the XR25 Developer and wait for the 30-second regenerating period to end.
- Switch the XR25 off again.  
Repeat this procedure **at least twice** to ensure that the regenerating unit (regenerating pumps, tubes) is thoroughly rinsed out.
- Switch on the film conveyor of the XR25 and run it for a working cycle. This will allow the circulation pumps to be rinsed out.
- Drain water from the tanks of the XR25 unit and from the regenerating containers.
- Remove D/E-Roller set and tub, and sponge clean with warm water.
- Refill with chemicals (referring to operating instructions XR25, 9000-600-03), refit tub, Roller set etc. and shut lid.
- Switch the XR25 Developer off once it returns to standby status.

# 12. CIRCUIT DIAGRAM

- C1** Anti-interference capacitor
- K1** Time lag relay
- M1** Pump motor
- X1** Terminal strip


# 13. TROUBLESHOOTING TIPS

<b>Malfunction</b>	<b>Possible Cause</b>	<b>Correction</b>
<b>1. Regeneration container is empty too quickly</b>	a) System or main surgery power supply switch operated too frequently.	a) Only turn off the system during the midday and leave on stand-by at other times.
	b) Regeneration time incorrectly set	b) Reset regeneration time Refer to section 10.2
<b>2. Films too grey due to insufficient fixing</b>	a) Regeneration tank empty	a) Refill the regeneration tank and operate the main switch with regeneration switched to on. If there is no improvement, replace the fixer and developer in the XR25.
	b) Pump timing too short.	b) Alter pump time by extending the regeneration time. Refer to section 10.2
	c) Either regeneration tank or the hoses wrongly connected	c) Drain chemicals from the unit. Clean thoroughly. Connect tank correctly and replace chemicals.
	d) Non-return valve on regeneration pump blocked or incorrectly mounted.	d) Clean non-return valve or replace if necessary.
		 <b>Be careful to position correctly.</b>
	e) Kink in regeneration hose (to developer)	e) Straighten hose.
	f) Shelf-life of film exceeded	f) Replace film. Note shelf-life of film.
	g) Shelf-life of chemicals exceeded	g) Replace chemicals. Note shelf-life of chemicals.
<b>3. Regeneration pump out of order</b>	a) Cable connection to regeneration loose or disconnected.	a) Check cable connections.
	b) Regeneration pump blocked	b) Check motor for mechanical faults (eg turn fanwheel)
	c) Regeneration time set to 0 secs.	c) Resert regeneration time. Refer to section 10.2

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# REGENERATING UNIT RU25 8/98

## MODEL 1700-840-00, 1700-850-00

When ordering spare parts the Unit model and, if present, the serial no. must be stated.

This information can be found on the nameplate of the Unit.

Pos. no.	Description	Order number
1	Connection cable .....	1700-844-00
2	Protective rubber sleeve .....	9000-123-14
3	Clamping nut .....	9000-312-39
4	Hose liner 90° .....	9000-315-57
5	Non return valve set .....	1700-825-00
6	Pump cover compl. ....	1700-826-00
7-8	Motor with bearing bush	
	230V .....	1700-813-00
	110V .....	1700-814-00
9	Hose clamp .....	9000-160-73
10	Spacer .....	1700-821-03
11	Captive nut .....	9000-429-10
12	Mounting plate .....	1700-821-01
13	Traction relief clamp .....	9000-137-04
14	Hose, red .....	1700-100-33
	Metre(s) .....	9000-318-54
15	Hose, black .....	1700-100-32
	Metre(s) .....	9000-318-53
16	Terminal block .....	9000-125-39
17	Cable band .....	9000-138-08
18	Hose, black .....	1700-821-04
	Metre(s) .....	9000-318-53
19	Hose, red .....	1700-821-05
	Metre(s) .....	9000-318-54
20	Collection trough .....	1700-820-01
21	Suction pipe .....	1700-822-00
22	Mixing vessel "Developer" .....	1330-015-01
23	Suction pipe .....	1700-823-00
24	Mixing vessel "Fixer" .....	1330-016-01
25	Cover .....	1700-821-02
26	Time lag relay, 230V AC, 50/60 Hz .....	9000-167-51
27	Condensator .....	1700-847-00

### Spare part sets

3-8	Double bellows pump compl.	
	230V .....	9000-341+06
	110V .....	9000-341+08

### Comments

+ = Exchange part