Conditions de stockage

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Stand der Information: 12.2018

Intended use

Bite registration material on vinyl polysiloxane basis.

Indications for use: Regofix is an addition-vulcanising silicone for bite registrations with high final hardness.

Material description

Regofix opaque and transparent are materials on the basis of A-silicones (vinyl-polysiloxanes) for bite registrations and as transparent matrices in extremely thixotropic consistency and high final hardness. The practical double cartridge system and the thixotropic characteristics ensure fast and clean application of the material onto the teeth rows.

Indication

- Bite registration
- Registration key in the field of intraoral bearing devices

Technical data according to DIN EN ISO 4823

Regofix	opak	transparent
Konsistenz	type 1	type 2
Processing time 23°C	30 ± 5 sec.	30 ± 5 sec.
Demoulding time 37°C	≥ 90 sec.	≥ 90 sec.
Final hardness	90 ± 5 Shore A	70 ± 2 Shore A
Recovery from deformation	> 99.8%	> 99.6%
Linear dimensional change	< 0.2%	< 0.2%

The mixing and processing times refer to a room temperature of 23°C and a relative air humidity of 50%. Lower temperatures prolong, higher temperatures shorten these times.


Processing

Regofix opak and **Regofix transparent** are offered exclusively in our double cartridge system. By help of an injector and mixing tips an ideal mixing of the components without loss of time is guaranteed.

Application in the dental surgery:

Apply **Regofix** by means of an injector directly from the double cartridge onto the teeth row in the desired quantity. This working step is to be carried out without interruption, as the mouth temperature leads to an accelerated curing. A controlled application is guaranteed, as **Regofix**

remains on the teeth rows due to the high thixotropy and does not flow away. Ensure that the bite is brought into the maximum intercuspidation within the given processing time and keep it closed for 90 sec. The fast polymerisation

evokes a slight, physiologically harmless temperature increase. After curing the registrate may be removed from the patient's mouth easily and quickly.

Application in the dental laboratory:

Regofix transparent can also be used as transparent form material during the hardening of composites. For this apply Regofix directly onto the wax modellation and dimension like a form.

Important: In order to guarantee a preferably straight forwarding of the light through the form, the external surface should be smoothed.

The hardening time at room temperature is approx. 3–4 min. In order to use the physical characteristics of Regofix transparent in the best possible way a waiting period of approx. 10 min. should be kept between the hardening and the further processing. An insulation of Regofix transparent is not necessary.

Processing 50 ml cartridge

When using double cartridges, clamp these into the injector and apply the mixing tips; the mixing is done automatically inside the mixing tip. It is recommended to dispose the initially mixed (pea-sized) amount of material.

Disinfection

After removal from the mouth rinse the bite registration under running water. Disinfection for 15 min. in a 2 % glutaraldehyde or another appropriate solution.

Important advice

- Remove residues of retraction agent with water spray before taking the bite registration.
- The vulcanisation may be disturbed by the use of latex gloves, hand creams, cleaning agents etc.
- Remove all residues of registration material from the patient's mouth.

- The product is exclusively designed for the field of application described above and only to be used by qualified and instructed personnel.
- When using double cartridges, the mixing canula should remain on the cartridge as sealing.

Storage conditions

At room temperature, i. e. 18°C–28°C.

Durability: see date of expiry.

Regofix opak and **Regofix transparent** are not to be used when expiry date has passed.

Keep away from children!

Only for dental use!

Forms of delivery

Regofix opak	
REF D504874	2 x 50 ml double cartridges + 12 mixing tips
REF D504875	8 x 50 ml double cartridges

Regofix transparent

REF D504876	2 x 50 ml double cartridges + 12 mixing tips
REF D504877	8 x 50 ml double cartridges

Accessories

REF D4970X2	Mixing tips, 40 pieces
REF D49701X2	Mixing tips, 100 pieces
REF D49603	Application tips, 40 pieces
REF D5007	DS Injector 1:1

This product was developed for dental use and has to be applied according to the instructions. In the case of damages resulting from ulterior or improper use, the manufacturer assumes no liability. Furthermore the user is bound to check the product according to suitability and possible use in self dependence before using it subject to its purpose, especially in case the purpose is not mentioned in the instructions.

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