SOPIRA® Ject

SOPIRA® Ject vario

-  **Gebrauchsanweisung**
-  **Instructions for use**
-  **Mode d'emploi**
-  **Instrucciones de uso**
-  **Istruzioni per l'uso**
-  **Bruksanvisning**

99001436/11

SOPIRA® Ject und SOPIRA® Ject vario

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SOPIRA® Ject and SOPIRA® Ject vario

Instructions for use

Please note:
Due to US Medical Device Legislation and specific requirements for devices marketed in the USA the wording of the user instructions for the US may slightly differ from the standard English version.

Special observance
Please read the following information and application notes carefully. The product may only be used by experienced medical staff in accordance with these instructions. Kulzer does not give any recommendation regarding a method of treatment. The treating medical specialist is responsible for the manner of application and for the selection of the patient. The non-observance or contravention of the instructions will lead to the expiry of the guarantee and to endangering the safety of the patient. If used in combination with further products, the instructions for use and the compatibility statements of those products shall always be taken into account. Please check the product for completeness and intactness before application. If substantiated cause for doubts regarding the completeness and the intactness are on hand, then do not use the product.

Product specification/compatibility
The SOPIRA Ject is available in two various designs. The SOPIRA Ject with slip-on mechanism to open the cartridge holder (the cartridge holder is firmly connected to the syringe), and the SOPIRA Ject vario with slip-on device for separate cartridge holders (the cartridge holders can be replaced). Both models offer self aspiration as well as manual aspiration by means of hooklets.


Indications
Infiltration- and nerve block anaesthesia in dentistry

Contraindications
The contraindications to be observed are those of the applied anaesthetic.

Complications
No complications are known to occur under professional, proper application.

Warning indications
The aspiration hooklets (C) have to be retracted position when inserting or removing the cartridges. They may otherwise be bent. Extend hooklets only to connect cartridge/plunger, please retract again when separating the connection! Retract the hooklets by turning the knurled screw (H) to the left hand side until the limiting stop becomes effective!

Special warning indications for use with other compatible products:
– It is essential to clean and sterilize non-sterile medical devices prior to first use!
– Treat a contaminated instrument immediately after each use (see pre-cleaning manual).
– The instrument must be visually and functionally checked before each use.
– Sort out damaged and defective instruments and replace them.

- A Screw thread for cannula.
- B Cartridge holder with viewing windows.
- C Three pivotable (swing-in and -out) aspiration hooklets at the end of the plunger provide for stable anchoring in the hollow stopper of the cartridge, and they ensure for smooth guidance when injecting and aspirating.
- D Joint with swing-out mechanism (for opening the cartridge holder).
- E Slip-on device for mounting the cartridge holder.
- F Finger rest.
- G Plunger.
- H Knurled screw to extend and retract the hooklets. The hooklets are extended by turning the knurled screw (H) to the right-hand side (clockwise), and turning left (counter clockwise) will retracted them. The screw – mechanism is equipped with a locking device. It is engaged as soon as the hooklets have been completely extended by turning to the right-hand side. An unintentional separation of the connection is impossible. The effect of the locking device may be noticed as a slight initial resistance when turning the knurled screw to the left-hand side to separate the connection.
- J Thumb plate ring.

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

Sequence of use

1. Insert the cartridge
a) SOPIRA Ject:
Pull the plunger (G) back completely. A spring resistance has to be overcome towards the end. Swing cartridge holder (B) out at a right angle. The holder is in open position. Insert the cartridge. Swing the holder back in alignment with the longitudinal axis. Thereby, completely retract the plunger backwards against the resistance of the spring. Release the plunger against the spring resistance. The cartridge will slip into the correct functioning position, centred in the cartridge holder by a guide ring.
b) SOPIRA Ject vario:
Pull the plunger (G) back completely. A spring resistance has to be overcome towards the end. Remove the cartridge holder (B) from the syringe by pulling it out of its rest laterally. Insert the cartridge into the cartridge holder. Slip the cartridge holder with the inserted cartridge into the rest of the syringe body. Thereby, completely retract the plunger backwards against the resistance of the spring. Release the plunger against the spring resistance. The cartridge will slip into the correct functioning position, centred in the cartridge holder by a guide ring.
2. The plunger (G) can be advanced forward several millimetres against the inserted cartridge. The front end of the hooklets (C) will be fitted tightly into the recessed void of the hollow stopper. The connection is created by turning the knurled screw (H) to the right-hand side (= extension of the hooklets). The locking device to safeguard the connection is engaged as soon as the hooklets have been completely extended.
3. The cannula is tightly screwed onto the cartridge holder.
4. Aspiration
The aspiration device is essentially based on the creation of a stable connection between the hollow stopper of the cartridge and the plunger of the syringe, which will bear tractive force and force of pressure. This jam-free fastened connection permits the performance of aspiration and injection without problems. General complications relating to local anaesthesia during dental treatment are often the result of an unintentional intravenous or intra-arterial injection. Failures of onset may also occur in case of an inadvertent intravascular injection, since the local anaesthetic is carried away quickly from the location of the application.
The aspiration check, the so-called „Jege arts“, will largely provide protection from injecting into a blood vessel.
5. The removal of the cartridge is performed in reverse order.

Treatment

General information

-  When working on contaminated instruments, observe the guidelines of the employers' liability insurance association and equivalent organizations relating to personal protection. Wear suitable protective equipment and ensure adequate vaccination protection.
-  **Risk of infection:** Incorrect treatment of the instruments may present a risk of infection to the patient, users and third parties as well as a risk of performance limitations of the instrument.
-  **In case of suspicion of or known Creutzfeldt-Jakob disease,** Kulzer recommends disposing of the instrument after a single use in accordance with the country-specific requirements.
-  At all times, observe the procedures, equipment and devices validated by the user/operator/central sterilisation and check them for compatibility with the information given here.
-  When preparing and using the solutions, the information given by the chemical manufacturer regarding concentration and exposure time must be observed. Failure to do so may result in damage to the instrument. Use only approved chemicals for treatment!

| Chemical | Manufacturer | Category | pH value |
|--------------------------|---------------|---|---------------|
| Deconex 23 | Borer Chemie | Neutral enzymatic cleaner | approx. 8.8* |
| Mucapur-AF | Merz Consumer | Alkaline cleaner | approx. 11.5* |
| Needsher mediclean forte | Dr. Weigert | Alkaline cleaner | 10.4–10.8* |
| Needsher z | Dr. Weigert | Acid neutralizing and cleaning agent/citric acid for neutralization | 3.0–2.6* |

* Data according to manufacturer data sheet
Further information on instrument preparation can be found at: www.a-k-i.org

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Preparation at the place of use

Clean a contaminated instrument immediately after each use. To prevent dryness, coarse contaminants, corrosive solutions and drugs must be removed immediately after use of the medical device e.g. removed by wiping and rinsing.

Transport

Use suitable transport containers to transport to treatment in order to exclude hazards or contamination from third parties.

Manual pre-cleaning in the treatment phase
Treat a contaminated instrument immediately after each use. Disassemble the instrument (instruments intended for disassembly) down to its individual parts (see chapter Disassembly Instructions).

Pre-cleaning the surfaces

Remove any visible contamination or coarse dirt from the surface of the instrument using a brush or sponge under running cold water.

Pre-cleaning the cavities/lumens

Clean the instrument's working channels, lumens and cavities under running cold water with the help of a suitable brush, water pressure gun and, if necessary, a flushing attachment.

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