Instructions for use
StarLoy[®] C

Gebrauchsanweisung
Mode d'emploi
Istruzioni per l'uso
Instrucciones para el uso

 Dentsply
Sirona

StarLoy® C

en	Instructions for use	04 – 09
de	Gebrauchsanweisung	10 – 15
fr	Mode d'emploi	16 – 21
it	Istruzioni per l'uso	22 – 27
es	Instrucciones para el uso	28 – 33

CoCr Dental Casting Alloy for the Ceramic Bonding Technique

StarLoy C is a highly corrosion-resistant non-noble ceramic-bonding alloy based on CoCr. It can be veneered with all dental ceramics covering the CTE range of $14.3 \mu\text{m}/\text{m}\cdot\text{K}$ (25 – 600 °C).

StarLoy C meets the requirements of the EN ISO 9693 Standard.

Material data:

Composition:

Co	59.4 %
Cr	24.5 %
W	10.0 %
Nb	2.0 %
V	2.0 %
Si	1.0 %
Mo	1.0 %
Fe	0.1 %

Technical Data:

Type:	Non-precious ceramic bonding alloy, type 5	
Vickers hardness:	330 HV	
0.2% yield strength:	630 MPa	
Tensile strength:	790 MPa	
Modulus of elasticity:	200 GPa	
Elongation at fracture:	3%	
Density:	8.8 g/cm ³	
Solidus point:	1270 °C	
Liquidus point:	1370 °C	
WAK value:	(25 – 500 °C)	14.0 µm/m·K
	(25 – 600 °C)	14.3 µm/m·K

Contraindications: Do not use in known insensitivity to one or several metals contained in the alloy. For dental use only.

Keep under lock and out of the reach of children.

Warning! StarLoy C contains cobalt and chromium. There are maximum allowable concentrations (MAC values) for the workplace for certain compounds containing these metals. All finishing with rotary instruments should be undertaken with a vacuum system in place. In addition, we recommend wearing a mouth guard and face shield.

This alloy should not be delivered to patients or processed by persons with a known hypersensitivity to one of its components.

Adverse effects: Reactions that may occur are allergic reactions against metals contained in the alloy or paraesthesia caused by electrochemical reactions. Systemic side effects caused by metals contained in the alloy have been reported in isolated cases.

Interactions: Avoid occlusal and proximal contact with different alloy types.

Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.

Processing:

Modelling: Before modelling, the dies should be coated with die spacer. The crown should be modelled to reduced anatomical shape. Care should be taken in the modelling process to ensure that no sharp angles or edges are created that might later result in cracks in the ceramic veneer.

When modelling and finishing the frameworks, make sure to maintain a minimum wall thickness of 0.3 mm for single crowns and 0.5 mm for abutment crowns.

Sprueing: If single-crown sprueing is preferred to the beam-casting method, the sprues should have a diameter of 3 – 4 mm. When using the beam-casting method the following sprue dimensions should be observed.

Sprue from the cone	3 – 4 mm diameter
Transverse sprue (beam)	4 – 5 mm diameter
Connection to crown	3 – 3.5 mm in diameter and 5 – 8 mm in length

Investing: Any phosphate-bonded, graphite-free investment compound can be used that is suitable for alloys that melt at high temperatures. When using metallic muffle rings, an adequate – if necessary multi-layered – muffle lining is required.

Follow the Instructions for Use for the investment compound employed.

Wax elimination/preheating: At 300 °C and 950 °C depending on flask size:

Flask size	1 x 20 min
	3 x 30 min
	6 x 45 min
	9 x 60 min

Required metal quantity: The quantity of metal required for casting is calculated according to the following rule of thumb: “Wax weight of the mould multiplied by the alloy density (8.8 g/cm^3) + approx. 10 g”.

Melting procedure: StarLoy C can be melted using the high-frequency procedure, induction casting or with an open flame.

High-frequency procedure: The casting is initiated 3 – 5 seconds after complete melting.

Induction casting (vacuum-die casting): The casting is initiated 3 – 5 seconds after the oxide film is broken.

Open-flame casting: To melt StarLoy C, the burner must be used such that the heat required for melting is uniformly distributed. Use the maximum burner setting.

Follow the Instructions for Use for the burner used.

The alloy will form an oxide skin on melting. It must be further heated until the metal under the oxide skin moves visibly due to the pressure of the flame. The casting is initiated after 3 – 5 seconds.

Caution! Do not use any melting /flux powder.

Crucible materials: Use only ceramic melting crucibles to cast StarLoy C. Before casting, the crucible (without metal) should be preheated in an oven.

Finishing: Carbide cutters are suitable for finishing. Diamond tools should not be used. Always ensure that the tools are sharp. Always use high cutting speeds and low pressure.

Ceramic veneering: After removal from the mould, the surface of the StarLoy C is roughened by intensive sandblasting with Al_2O_3 (250 μm) at a pressure of 3 – 4 bar. Ceramic bonding is performed in accordance with the Instructions for Use for the ceramic veneering material.

Comply with the recommendations of the ceramic manufacturer for NP alloys.

Oxide firing: There is no fundamental need for oxidizing.

Polishing: To reduce oxidation on the uncovered alloy fractions, it is advisable to pre-polish them before ceramic firing. After firing, the alloy surface is first polished to a satin lustre using a rubber stone or wheel and the polished to a high gloss using polishing or diamond paste. Polishing the metal surface will increase its corrosion resistance. It reduces the release of alloy components and further improves the biocompatibility of the alloy. Residual polish can be removed with a steam jet cleaner or in an ultrasound bath.

Soldering:

Before firing:

Degudent solder U1W (1120 °C)

Flux:

Oxynon

Reusing old material: Neither sprues or cones nor failed castings may be reused, as this may result in oxides being embedded in the casting object, compromising its quality.

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DeguDent GmbH