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GB

Instructions for use

## Composite materials for provisional fillings and temporary sealings, light curing

Fig. 1

Fig. 2


Fig. 3

Fig. 4

Fig. 5

### 1. Range of applications

Provisional fillings for one to multisurface cavities. Temporary restorations in the inlay and onlay technique without the use of a temporary cement. Sealing accesses to implant screws. Cavity lining of prefabricated transparent plastic temporary restorations.

### 2. Preparation

Prepare, place a cavity liner and take the impression in the usual manner. Rinse and dry the cavity; place a matrix if required.

### 3. Applying smartprotect® Desensitizer

Clean and dry the dentine surfaces. Use cotton wool rolls to keep the working area dry. Apply smartprotect® with a smartbrush or application brush to the dentine and massage in for 10 sec. (Fig. 1). Carefully blow off any excess and dry with the air jet (Fig. 2). Do not overdry the dentine.

### 4. Applying tempofill®2

Insert tempofill®2 into the cavity with a suitable instrument, e.g. ball-end plunger, spatula etc. (Fig. 3). When filling the cavity, contour with a slight excess of material. With inlays and onlays keep the approximal contact points free of tempofill®2 by using dental floss or a matrix before curing. Close the tempofill®2 syringe immediately after use. Store in a cool place if it is intended to be used for a longer period.

### 5. Curing and preparation

Light cure with a halogen lamp for 20 - 40 sec., depending on the layer thickness. Cure layers up to 4 mm for 20 sec. and layers of 5 mm for 40 sec. (Fig. 4). Then check the occlusion. Any excess can be cut away with a scalpel. Silicone rubber finishers and tungsten carbide finishers are suitable for trimming.

### 6. Removing the temporary

To remove the temporary, insert a suitable instrument (e.g. probe) into the temporary and remove the temporary. Then clean the cavity and try in the restoration (Fig. 5).

### Important working hints

- With large cavities a composite (e.g. smartrepair®) can be used to increase stability in the occlusal area of the temporary restoration. Polymerise the tempofill®2 separately or together with the filling material.
- tempofill®2 can be more easily contoured by coating the instrument lightly with smartprotect® Desensitizer.
- If there is minimum or inadequate retention, the temporary can be retained using a eugenol-free temporary cement.
- tempofill®2 may bond with denture acrylics or cavity liners (e.g. light curing glass ionomer cement). In this case use a suitable separating medium (e.g. glycerine gel).
- Avoid cavity lining materials containing eugenol, as these can detrimentally affect curing.
- Avoid uncured tempofill®2 coming into contact with the skin, mucosal tissue and eyes. Wear suitable eye protection and gloves! In rare cases repeated contact can cause sensitisation.
- Biocompatibility is only guaranteed with complete polymerisation.
- Please follow the instructions of the safety data sheet!

Only for dental use by qualified personnel.

### Indications for use:

- provisional fillings of cavities, temporary inlays and onlays
- temporary sealings in the coronal area in the field of endodontia
- sealing of spiral canals in implant abutments
- for lining / filling of temporary crown forms

### Technical data:

- **Weight:** 3 g (syringe)
- **Mixing ratio:** one component material
- **Product colour:** light opaque
- **Setting time:** light curing 20 - 40 sec.
- **Wavelength:** 430 nm - 480 nm
- **Application:** At 23 °C ± 2 °C / 73 °F ± 3.6 °F
- **Storage:**

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