KULZER
MITSUI CHEMICALS GROUP

Translux® 2Wave

- DE** Gebrauchsanweisung
- GB** Instructions for use
- FR** Mode d'emploi
- ES** Instrucciones de uso
- IT** Istruzioni per l'uso
- PT** Instruções de uso
- NL** Gebruiksaanwijzing
- SE** Bruksanvisning
- DK** Brugervejledning
- NO** Bruksanvisning
- FI** Käyttöohjeet
- GR** Οδηγίες χρήσης
- RU** Инструкция по применению

DE

GB

FR

ES

IT

PT

NL

SE


DK

NO


FI

GR


RU


4.2.1 / 4.2.2




4.2.3

6.1

5.4

5.5


6.3


5.3


4.2.4

8.5.2


8.6.1

8.6.2

Made in Italy

Manufacturer:
Kulzer GmbH
Leipziger Straße 2
63450 Hanau (Germany)

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Stand: 2017-08

Table of contents

1	Application	3
1.1	In general	3
1.2	EC – Declaration of Conformity	3
1.3	Manufacturer's declaration	3
2	Instructions on using the unit safely	3
2.1	Explanation of symbols and accompanying word(s)	3
2.2	Transport damage – Unpacking and inspection	4
2.3	Owner's obligations	4
2.4	Machine journal	4
3	Intended use	4
3.1	Description of the device	5
3.2	Safety requirements	5
3.2.1	General safety instructions on operating the device	5
3.2.2	Safety Precautions	6
3.2.3	BATTERY – Safety precautions	6
3.2.4	BATTERY – Information	7
3.2.5	LIABILITY	8
4	Unit description	8
4.1	Delivery contents	8
4.2	Description of the control and operating elements	8
4.2.1	Hand piece	8
4.2.2	Allignment key	8
4.2.3	Hand piece display	8
4.2.4	Charging unit	10
4.2.5	Acoustical signals and indications on the display – hand piece	10
4.2.6	LED information concerning the radiation emitted	11
5	Installation and putting into initial operation	12
5.1	Safety requirements during Installation	12
5.2	Connection of the charging unit to the mains power outlet	12
5.3	Inserting the battery into the hand piece	13
5.4	Inserting the fibre optic into the hand piece	14
5.5	Battery charging	14
6	Operation	15
6.1	Selecting the curing time	16
6.2	Activating and deactivating the exposure cycle	16
6.3	Measurement of the Light Intensity	16
6.4	Sleep Mode	17
6.5	Forced Sleep Mode	17
6.6	Additional information on the unit operations	17
7	Diagnostic signals	17
7.1	Faulty battery	17
7.2	Faulty LED	17
7.3	Overheat protection	18
7.4	BUILT-IN Temperature overload sensor defective	18
7.5	Battery-low signal	18


GB

- 1 -


NOTE

Please use the book icon to find the chapter images or further information on the fold-out cover.

Please
fold-out front
cover

1x/2x

Please
fold-out rear
cover

8	Cleaning, disinfection and sterilisation	18
8.1	Disassembly	18
8.2	LIGHT GUIDE – Cleaning and disinfection	18
8.2.1	LIGHT GUIDE – Inspection before the sterilisation	19
8.3	PROTECTIVE LIGHT CONE – Cleaning and disinfection	19
8.3.1	PROTECTIVE LIGHT CONE – Inspection before the sterilisation	19
8.4	LIGHT GUIDE and PROTECTIVE LIGHT CONE – Packaging before the sterilisation	19
8.4.1	LIGHT GUIDE and PROTECTIVE LIGHT CONE – Sterilisation	19
8.4.1.1	Sterilisation method	19
8.4.1.2	Sterilisation parameters	20
8.5	HAND PIECE – Cleaning and disinfection	20
8.5.1	Cleaning Mode	20
8.5.2	Cleaning and disinfection	20
8.6	CHARGING UNIT – Cleaning Casing	21
8.6.1	CHARGING UNIT – Cleaning charging contacts	21
8.6.2	CHARGING UNIT – Replaceable charging contacts	22
9	Shelf storage	22
10	Disposal	22
11	Troubleshooting	23
12	Technical specifications	25
12.1	Electromagnetic compatibility EN 60601-1-2	26
12.2	Charging unit identification plate	32
12.3	Hand piece serial number	32
13	Warranty	32
14	Service	32
15	Document history	32


GB

-2-


NOTE

Please use the book icon to find the chapter images or further information on the fold-out cover.

Please
fold-out front
cover

1x/2x

Please
fold-out rear
cover

1 Application

1.1 In general

Translux 2Wave is a registered trademark of Kulzer GmbH.

These operating instructions apply to:

Order no.	Type Features	Date
66055013	Translux 2Wave – LED Light curing unit	2017-07 / 99000873/11

1.2 EC – Declaration of Conformity

We, Kulzer GmbH, Leipziger Straße 2, 63450 Hanau, hereby declare that the design, structure and version introduced to the market of the unit described below correspond to the relevant safety and health requirements of the EU Directive. This declaration shall become invalid in the event of an alteration made to the unit without our prior authorisation.

1.3 Manufacturer's declaration

We hereby confirm that the unit Translux 2Wave complies with the directive 93/42 EC and the standards IEC 60601-1 and IEC 60601-1-2.

We hereby confirm that system of Quality Assurance is guaranteed in conformity with the EN ISO 13485.

Hanau, 2017-07 – Kulzer GmbH, Leipziger Straße 2, 63450 Hanau (Germany).


GB

- 3 -

2 Instructions on using the unit safely

Please read this manual and follow its instructions carefully. The words **WARNING**, **CAUTION** and **NOTE** carry special meanings and must be carefully reviewed (see the following paragraph). Always keep this manual within easy reach.

2.1 Explanation of symbols and accompanying word(s)

Symbols	Accompanying word(s)	Explanation
	WARNING	<i>Implies that death or serious injury could occur if the advice is not followed.</i>
	CAUTION	<i>Implies that minor injury or device damage could occur if the advice is not followed.</i>
	NOTE	<i>Implies advice that is not related to harm.</i>
	WARNING OF EXPLOSION	<i>WARNING of explosion – observe information and warnings in the chapters and sections marked with the following pictorial symbols.</i>
	DISCONNECT FROM MAINS	<i>WARNING of electric shock. Before carrying out the cleaning and disinfection procedures disconnect the charging base from the electrical power source!</i>
	NEMKO	<i>Nemko mark – UL / CSA conformity.</i>
	WEEE / ElektroG	<i>Disposal of old equipment according to WEEE (European Directive) or Electrical and Electronic Equipment Act (German Law ElektroG).</i>
		<i>Applied part: Type B applied part-in conformity with technical standards IEC 60601-1.</i>
		<i>Double insulation: Class II device.</i>
		<i>Power supply: Alternate current.</i>
		<i>Direct current.</i>
		<i>Switching power supply plug polarity – Centre positive.</i>
		<i>Use in closed spaces only.</i>

Symbols	Accompanying word(s)	Explanation
		Device manufactured in conformity with directive 93/42/EEC including technical standards IEC 60601-1 and IEC 60601-1-2.
		Symbol ISO 7010-W001 Generic warning symbol
		Consult instructions for use.
		CAUTION, consult accompanying documents.
		Operating switch: "STAND-BY".
		"Time / Program" selecting switch.
		Keep away from rain.
		Storage temperature limitation (–20°C (–4°F) to 40°C (104°F)).
		Storage humidity limitation (45% – 85% RH).
		Storage atmospheric pressure limitation (500 – 1060 hPa).

2.2 Transport damage – Unpacking and inspection

The unit is sensitive to shock, because it contains electrical components. Particular care must therefore be taken during transportation as well as storage. The equipment shipped by Kulzer was checked thoroughly prior to shipment. The unit is correctly protected and packed when delivered.

Check the unit for damage after receiving it. If it is damaged, report this to the transportation company within 24 hours of delivery. Under no circumstances, install or work with a damaged unit and / or accessories.

2.3 Owner's obligations

In addition to complying with the statutory regulations specified, the owner must ensure that statutory requirements are observed and implemented in the workplace, i.e. he must train his personnel and comply with industrial safety legislation and all further valid regulations or laws in force.

For work on and with the unit, written instructions in an understandable form and in the relevant language of the respective employees must be provided by the owner using the operating instructions and on the basis of the work to be carried out.

2.4 Machine journal

We recommend that you keep a "machine journal" and document all tests and major work (e.g. servicing, changes) in it.

3 Intended use

The Translux 2Wave is a LED dental curing light unit intended for use in the oral cavity for the polymerisation of light curing dental materials such as adhesives and filling materials that are activated in the 385 – 510 nm wavelength range.

NOTE

Although most light curing dental materials are activated within this wavelength range, in case of uncertainty, verify the manufacturer's specifications or contact the manufacturer.

WARNING – Qualified and specialised personnel.

The device must be used exclusively by specialized personnel with appropriate medical background; no training activities for the use of the device are foreseen by Kulzer. The use of the device does not cause side effects if it is used correctly.

3.1 Description of the device

The Translux 2Wave uses a source of light consisting of a very-efficient duo-chromatic LED diode with a wavelength comprised between 385 and 510 nm.

The entire light radiated by the Translux 2Wave is used to activate the camphorquinone photoinitiators as well as Lucirin TPO and PPT. This produces excellent polymerisation results with low heat emission.

The Translux 2Wave is shipped with a 360° rotating light guide with 8 mm diameter.

The Translux 2Wave consists of a charging unit and a hand piece, powered by a removable rechargeable Lithium-ion battery.

The hand piece is equipped with a “**Sleep Mode**” to minimise the unit’s energy consumption.

The hand piece switches to “**Sleep Mode**” if it is not used for a period of about 5 minutes.

The hand piece features a display with symbols / icons for different modes, curing time selection and the battery status.

The Translux 2Wave can operate with four different exposure cycles:

- **Exposure with slowly increasing output: “Soft-Start”** (cycle duration 20 seconds).
The “**Soft-Start**” means increase of light intensity from 50% up to 100% within 2 seconds.
- **Exposure with constant output:** cycle duration 20 seconds, 10 seconds or 5 seconds.

NOTE – Observe the exposure times of the materials specified by the manufacturer.

For materials with specific exposure times (e.g. 30 / 40 seconds for dark composites) repeat the exposure several times if necessary.

3.2 Safety requirements

The effectiveness of the safety instructions with regard to personnel protection, handling the unit and handling the processed product is to a large degree dependent on the behaviour of the staff working with the equipment.

WARNING

Prior to initial operation, carefully read through these Operating Instructions and observe the information in order to avoid errors and subsequent damage, particularly damage to health. In addition to the instructions in this operating manual, comply with the national laws, regulations, and directives in your country when setting up and operating the unit.

3.2.1 General safety instructions on operating the device

WARNING – INTENDED USE

Use the unit only for the purpose for which it is intended (see paragraph 3 Intended use). Any other uses are strongly contra-indicated. Any non-compliance of this regulation can result in serious injury of the patients or operators of the unit as well as damage or faults to the unit. Kulzer is not liable for consequential damages due to the non-observance of the aforementioned provision!

WARNING

The user is responsible for testing the Translux 2Wave for its use and suitability for the intended purposes.

***Never** point the light towards the eyes! Direct or indirect exposure of light into the eyes must always be prevented by wearing protective goggles that filter out blue and ultraviolet light. When curing, position the light guide tip directly over the material that needs to be cured and avoid extraneous exposure to soft tissues. Eyes, gums, soft tissue and other body parts must not be exposed to the light beam (if necessary, these parts must be covered with protection). Use a rubber dam to protect the soft tissues, i.e. gingiva during treatment. Exposure must be restricted to the area of the oral cavity in which clinical treatment is executed.*

WARNING – CONTRA-INDICATIONS

***Do not** use the Translux 2Wave on patients fitted with a pacemaker or other implantable electronic devices. This also applies to the operator.*

The unit must not be used by children, the blind or the deaf and also not by / on persons who are prone to possible epileptic seizures. The above pre-named persons may not recognise the risks of use and operation. Due to safety reasons for the operator and the patients with potential risk of epileptic seizures, use and operation of the device is prohibited.

WARNING

No modification of this equipment is allowed.

The owner / user must ensure that:

- The unit is not used on patients with a positive anamnesis after light testing, for example with urticaria solaris and / or porphyria as well as on patients who have been treated with photosensitising medication.
- The unit is not used on patients that have had surgery for eye cataracts, which are particularly sensitive to light, unless they use suitable protective measures, such as the use of eye protectors / spectacles that filter out blue and ultraviolet light.
- Patients whose anamnesis includes retinal diseases must have treatment with Translux 2Wave authorised by their ophthalmologist.
- A specialist doctor is consulted in all potential risk cases.

WARNING – Temperature of applied parts

Max. temperature of the applied parts (light guide and protective light cone): 41 °C (106 °F).

Max. temperature of the part considered as applied part (metal front cone) as it can come into contact with the patient: 46 °C (115 °F).

Kulzer declares that if the device is used in conformity with the instructions contained in this manual, no tissue heating will occur.

WARNING – Control of infections.

The light guide and the protective light cone must be cleaned, disinfected and sterilised (steam) before every use.

*The light guide and the protective light cone are delivered in **NON STERILE** condition and must be sterilised before being used for the first time. See paragraph 8.4.1.2 Sterilisation parameters.*

GB

-6-

3.2.2 Safety Precautions

WARNING – RISK OF EXPLOSION

Do not install the device where there is the risk of explosions.

The device must not be used where an inflammable atmosphere is present (anaesthetic mixtures, oxygen, nitrous oxide, etc.). Use the device in a well ventilated area.

WARNING

Before using the unit, check the mains power cable and plug for damage. If they are damaged, do not connect the unit to the mains.

Use only original spare parts or accessories from Kulzer GmbH. Kulzer provides a light guide suitable for use with the Translux 2Wave. Other light guides must not be used. Kulzer GmbH will not accept any liability for damage resulting from the use of non original Kulzer spare parts or accessories.

3.2.3 BATTERY – Safety precautions

Kulzer GmbH assumes no liability for battery problems that occur when the Safety Precautions listed in the paragraph below are not carefully adhered to.

WARNING

Misusing the battery may cause the battery to get hot, rupture, or ignite and cause serious injury. Be sure to follow the safety rules listed in the below paragraph.

Use original Kulzer batteries only!

The use of non Kulzer batteries or non-rechargeable batteries and / or primary batteries is a potential danger and can damage the device.

Keep the battery out of the reach of children!

Never open, puncture or crush the battery – it contains poisonous substances.

If the battery shows corrosion, emits an unusual smell, or loss of fluids, remove the battery from the hand piece immediately.

Do not expose the battery to water or salt water, or allow the battery to get wet. Do not store the battery in high humidity or in a place which may expose the battery to rain.

In the event that the battery leaks and the fluid gets into one's eye, do not rub the eye. Rinse well with water and immediately seek medical care. If left untreated the battery fluid could cause damage to eye.

WARNING – RISK OF EXPLOSION

Use only the Translux 2Wave charging unit and the input power supply, which are provided with the unit, to charge the hand piece's battery. Never attempt to charge the Translux 2Wave battery by using any other charging unit. The use of any other charging unit can result in the damage of the battery, risk of explosion and fire!

Never use the Translux 2Wave recharging unit to charge other types of battery or other units with a rechargeable battery! Use only original batteries from Kulzer.

Do not pierce the battery with any sharp objects, strike the battery with a hammer / tools, step on the battery, or otherwise subject it to strong impacts or shocks.

Do not place the battery in fire or apply heat to the battery. Never short circuit the end elements of the battery with any metal objects due to the risk of burns, fire and explosion. Never carry or store the battery together with necklaces, hairpins, or other metal objects.

WARNING – Do not disassemble or modify the battery!

Safety mechanisms are incorporated in the battery which if damaged may cause the battery to generate heat, rupture, explode or ignite.

WARNING – Do not place the battery on or near fires, or other high-temperature locations.

Do not place the battery in direct sunshine. Doing so may cause the battery to generate heat, rupture or ignite. Using the battery in this manner may also result in a loss of performance and shortened life expectancy.

WARNING – In case of fire, DO NOT THROW WATER ON A BURNING BATTERY!

A class C fire extinguisher (according to European Regulation EN 3) must be used.

WARNING – Damaged or exhausted battery disposal.

Dispose a damaged or exhausted battery after taking the necessary steps to prevent external short circuit by the following method:

After having insulated the terminals of the battery with insulating tape, dispose off under the law or in accordance with the local council.

3.2.4 BATTERY – Information

- **Time to charge empty or new battery: approx. 3 hours.**

NOTE – NEW Battery: first charging

For the first charge, this process will take about **3 hours**. New batteries or batteries stored for a long time may take more time to charge. The battery will reach its full capacity after a few full charge / discharge cycles.

NOTE

Place the hand piece into the charging unit after each treatment, or when not in use.

- **Operating condition:**

Operating temperature: 10°C (50°F) up to 35°C (95°F).

Relative humidity: 45% to 85% RH.

Atmospheric pressure: 800 hPa to 1060 hPa.

CAUTION

The temperature and relative humidity ranges over which the battery can be charged is 10°C (50°F) to 35°C (95°F), and 45% to 85% RH. Charging the battery at temperatures outside of this range may cause the battery to become hot or damaged. Charging the battery outside this temperature range may also harm the performance of the battery or reduce the battery's life expectancy.

- **Battery transportation and storage conditions:**

Recommended temperature and relative humidity ranges:

–20°C (–4°F) to 40°C (104°F) and 45% to 85% RH.

Always store the battery charged and not for longer than 5 months.

Ambient air pressure: 500 hPa to 1060 hPa.

NOTE

Always fully charge the battery before an extended period of inactivity of the device and at least once every 5 months.

In case of extended period of inactivity replace the battery from the hand piece.

- **Battery Typical Life Characteristics:**

The Translux 2Wave battery has a typical life expectancy of 300 – 400 charging cycles depending on the application and environmental conditions in which it is used.

NOTE

The battery is a chemical product utilising a chemical reaction, so the performance will deteriorate not only through use but also with the passage of time even if it left unused.

3.2.5 LIABILITY

Work on electronic equipment in the unit must be carried out by **Kulzer, Kulzer service partners or qualified, specialist personnel only** and when the equipment is in a safe status (de-energised).

Only use authorised original spare parts and accessories. Avoid other parts as they hide unknown risks.

The operability and safety of the unit is only guaranteed if the necessary tests, maintenance and repair work is carried out by **Kulzer, Kulzer service partners or trained, specialist personnel**.

For possible damage coming from a fault / malfunction of the unit due to improper repair, which was not carried out by **Kulzer service partners or by personnel trained by us** or in cases where original spare / accessory parts were not used during a part replacement, Kulzer GmbH is **not liable**.

4 Unit description

4.1 Delivery contents


- 1 Charging unit
- 2 Hand piece
- 3 Light guide, $\varnothing = 8 \text{ mm}$
- 4 Protective light cone
- 5 Rechargeable Lithium-ion battery pack
- 6 Universal input power supply (100 V – 240 V (\sim)) including international mains socket adapters:
6a = Europe type 1
6b = UK type 2
6c = US / JP type 3

4.2 Description of the control and operating elements

4.2.1 Hand piece

- 1 "STAND-BY" button
- 2 Display
- 3 "Time / Program" button for selecting the curing program and the exposure times
- 4 Charging unit contacts
- 5 Alignment key


4.2.2 Alignment key

The alignment key (5), located at the back bottom part of the hand piece's case, allows the correct insertion of the hand piece in the charging unit holder.

When inserting the hand piece into the charging unit, the alignment key must be aligned with the guiding slot positioned inside the charging unit.

A beep confirms that the contacts are in the correct charging position!

4.2.3 Hand piece display

The information provided by the hand piece display are described here below:

① Battery status

The battery icon indicates the charging state and eventual anomalies of the battery. The black bars in the battery icon will light up progressively during the charging phase. When the battery is fully charged, all four internal bars light up.

NOTE

If the battery is completely empty, the charging operation starts with a "prequalification state" during which the Translux 2Wave micro-processor checks the functional charging parameters of the battery. During the prequalification state the first black bar inside the battery icon blinks. If the functional charging parameters are correct, the regular charging phase starts and all the black bars in the battery icon will light up progressively.


NOTE – Faulty Battery

In case of a defective battery, the battery icon, with no internal black bars, will continuously blink.

A "beep" sound is emitted every second for a period of 30 seconds.

After 30 seconds only the battery icon frame remains blinking.

IMPORTANT

This error condition is detected and shown only when the hand piece is placed in the charging unit.

NOTE – Battery-low signal


When after frequent use the battery charge drops to the minimum level, the Translux 2Wave microprocessor allows another few exposure cycles without re-charging the battery. This condition is signalled at the end of each cycle by “2 beeps”.

At the end of the remaining battery charge, the unit will have no further exposure cycle and this condition is shown on the display by the battery icon with no internal black bars.

2 Icon “Soft-Start” program

The “Soft-Start” program is selected, by the “Time / Program” button, two round spots simultaneously light up next to the icons and :


The “Soft-Start” means increase of light intensity from 50% up to 100% within 2 seconds.

3 Icons showing the curing programs: 20, 10 or 5 seconds

These programs are selected by means of the “Time / Program” button .

A round spot lights up next to the curing programs selected.

GB

- 9 -

4 “LED” light source indication

This indication blinks at the top left of the display when the Translux 2Wave LED light source is defective.

5 “T” indication – Overheating protection

The hand piece has become overheated in the course of use and the temperature overload sensor has activated.

The activation of the overheat protection temporarily prevents the use of the unit for a few minutes.

In this mode the letter “T” comes on, steady, at the top right of the display and three beeps signals are emitted once.

blinking

steady

NOTE

The hand piece will automatically exit from the cooling state upon reaching operational temperature and the letter “T” (steady) on the display switches off.

6 “T” indication – Temperature overload sensor defect

In case the built-in temperature overload sensor is defective the letter “T” comes on (blinking) at the top right of the display.

NOTE

If the temperature overload sensor fault occurs during an exposure cycle, the software allows completing the exposure. No further exposure cycle is allowed by pressing the “STAND-BY” button and “4 beeps” are emitted.

blinking

4.2.4 Charging unit

Top view elements

- 1 Charging unit cover
- 7 Hand piece holder with alignment key and replaceable charging unit contacts
- 8 Light intensity sensor
- 9  – GREEN colour LED
- 10  – 2-colour LED (GREEN / YELLOW)

Function: It indicates the charging unit is powered up.


Function: Light intensity information.

GREEN: Indicates that the light intensity, measured with the built-in light meter is suitable for effective treatment.

YELLOW: Indicates that the light intensity is insufficient.

Bottom view charging unit

- 11 Power supply socket
- 12 Replaceable charging contact insert

GB

- 10 - 4.2.5 Acoustical signals and indications on the display – hand piece

Function / Error Condition	Hand piece Button operation	Acoustic signal/s emitted	Display indication
“SOFT-START” PROGRAM: 20 seconds and  .	<ul style="list-style-type: none"> Press the “TIME / PROGRAM” button  to select the “SOFT-START” PROGRAM. Press the “STAND-BY” button  briefly to start the exposure. 	<ul style="list-style-type: none"> 1 beep when pressing the “STAND-BY” button. 1 beep after 10 seconds of exposure time. 1 beep at the end of exposure time. 	
POLIMERISATION PROGRAMS: – 20 seconds – 10 seconds – 5 seconds.	<ul style="list-style-type: none"> Press the “TIME / PROGRAM” button  to select the exposure time on the display: 20 s – 10 s – 5 s. Press the “STAND-BY” button  briefly to start the exposure. 	<ul style="list-style-type: none"> 1 beep when pressing the “STAND-BY” button. 1 beep after 10 seconds of exposure time (only 20 s Mode). 1 beep at the end of exposure time. 	A round spot appears side the exposure time selected.
INTERRUPTION OF THE EXPOSURE CYCLE	<ul style="list-style-type: none"> The exposure cycle can be interrupted at any time and program in use by pressing the “STAND-BY” button . 	<ul style="list-style-type: none"> 1 beep when pressing the “STAND-BY” button. 	
BATTERY-LOW SIGNAL The residual charge level of the battery is sufficient for a few exposure cycles.		<ul style="list-style-type: none"> 2 beeps at the end of the exposure cycle. 	
DEFECTIVE BATTERY SIGNAL		Only with hand piece placed into the charging unit. <ul style="list-style-type: none"> 1 beep every second for a period of 30 seconds. 	 blinking
LED Light source defective		No acoustical signal emitted.	 blinking
OVERHEATING CONDITION		<ul style="list-style-type: none"> 3 beeps, the light emission is stopped. 	 steady

Function / Error Condition	Hand piece Button operation	Acoustic signal/s emitted	Display indication
BUILD-IN Temperature overload sensor defective		<ul style="list-style-type: none"> • 4 beeps at the start of an exposure cycle. 	 <p>blinking</p>
Charging circuit of the hand piece defective Note: This condition is detected only with hand piece placed inside the charging unit.		<ul style="list-style-type: none"> • 1 beep per second for the length of 30 seconds. 	 <p>blinking</p>
“CLEANING MODE”	<ul style="list-style-type: none"> • With hand piece in the “Sleep Mode” condition, PRESS and KEEP PRESSED SIMULTANEOUSLY the “TIME / PROGRAM” and the “STAND-BY” buttons   for 5 seconds. 	<ul style="list-style-type: none"> • 1 beep when the “TIME / PROGRAM” and the “STAND-BY” buttons are simultaneously pressed. • 1 beep at the beginning of the “Cleaning time”. • 1 beep at the end of the “Cleaning time”. 	Display is switched off.

GB

4.2.6 LED information concerning the radiation emitted

Risk Group 2
CAUTION. Possibly hazardous optical radiation emitted from this product. Do not stare at operating lamp. May be harmful to the eye
Product tested against IEC62471

WARNING

Photobiological safety of lamps and lamp systems IEC 62471. According to the Standard IEC 62471, the device results in risk class 2 (moderate risk) concerning a retinal risk from blue light or thermal retinal risk. The above warning indications are applied to the device package.

5 Installation and putting into initial operation

The following sections of this paragraph contain the information and instructions which must be followed to ensure successful and faultless operation of the device. Before operating the Translux 2Wave follow the instructions in the sequence described below.

The medical system is in compliance with the safety standards only if it is installed by following the instructions described here below.

The device must be installed in a suitable place that is handy for its use. Place the charging unit on a sturdy, dry, flat and horizontal surface.

IMPORTANT

The hand piece is intended for use within the patient environment whereas the charger unit and the Switching Power Supply adapter are not intended for use within the patient environment.

NOTE

The patient environment is defined as an area 1.5 m from the patient (as per IEC 60601-1, third edition and IEC 60601-1-1).

WARNING

Always inspect the unit and its components for any damage before installation. If damage is apparent, do not proceed with the installation of the device.

The operator must not come simultaneously into contact with the elements outside the patient environment (charger units and switching power supply adapter) and the patient.

Do not connect other external components to the medical system.

GB

- 12 -

5.1 Safety requirements during Installation

WARNING

The electrical installation of the premises where the device is installed and used must comply with existing legislation and the safety specifications relevant to power systems.

Install the device in a place where it will be protected from blows and from accidental splashes of water or other liquids. Do not expose the device to direct sunlight, UV light.

Do not install the device on top or close to sources of heat. When installing, make sure that there is adequate air circulating around the device. Do not install the device close to solvents or flammable liquids as they could damage the plastic casings of the device.

The plug of the Switching Power Supply (Direct Plug-In), model PSAC05R-050L6, is considered the mean of isolating of the medical equipment from supply mains. When the plug is connected, the power supply adapter must be easily accessible. Leave sufficient free space around it. The medical equipment shall not be positioned to make difficult to operate on the disconnection device (so on the plug of the Switching Power Supply).

Never short circuit the end elements of the battery with any metal object as it could cause a risk of burns, fire and explosion.

The device is transportable but it must be handled with care. Transport the device in a horizontal position.

Do not let the device be exposed to shaking or vibrations. For sizes and weight see paragraph 12 (Technical specifications).


WARNING – RISK OF EXPLOSION

Do not install the device where there is the risk of explosions. The device must not be used where an inflammable atmosphere is present (anaesthetic mixtures, oxygen, nitrous oxide, etc.). Install the device in a well ventilated area.

5.2 Connection of the charging unit to the mains power outlet

The Translux 2Wave is supplied with a separated universal-input unit (6) that accepts 100 to 240 V (), 50 / 60 Hz.

WARNING

Before connecting the charging unit (6) to the mains network, check carefully that the voltage and frequency of the mains network correspond to the values stated on the rating plate of the device. The rating plate is fitted to the bottom of the charging unit.

Use only the switching power supply adapter which is provided with the device. The use of any other power supply adapter can result in damage of the battery.

Do not connect the switching power supply adapter to the mains socket by using a multiple portable socket-outlet.

**NOTE**

The Translux 2Wave power supply unit (6) is supplied with three interchangeable adapters:

Fig. (6 a) = Europe type 1

Fig. (6 b) = UK type 2

Fig. (6 c) = US / JP type 3

- ① Select the correct adapter between the ones available: Europe – Type 1, UK – Type 2, US / JP – Type 3.
- ② Each adapter has a slot into which the retaining latch on the power supply unit fits.
- ③ With the slot oriented toward the power supply unit, place the adapter onto the power supply unit so that the slides on the adapter are aligned with their respective slots on the power supply unit.
- ④ Pressing down on the adapter and power supply unit, slide the adapter into the power supply unit until a click is heard and the adapter is locked in place.

**NOTE**

Pressing down on the adapter during installation ensures that all the slides are held by the retainers on the power supply unit. If the adapter feels loose, remove it and re-install.

- ⑤ To remove the adapter from the power supply unit, press the latch retainer button (marked with the word “PUSH”) and slide the adapter away from the power supply unit and remove the adapter.
- ⑥ Retain unused adapters for future use.
- ⑦ Connect the plug of the power cord into the socket on the underside of the charging unit.

GB

- 13 -

**WARNING**

Check the power cable, the power supply unit and the adapter regularly to make sure that they are in good condition.

If damage is apparent do not use damaged part until it is replaced.

Use original Kulzer spare parts and accessories only.

- ⑧ Plug the power supply unit with its adapter into the mains power outlet.
- ⑨ The GREEN LED on the charger unit will light up (icon).

5.3 Inserting the battery into the hand piece

The Translux 2Wave is supplied with a powerful Lithium-ion rechargeable battery.

The battery is supplied separately in the device packaging and must be inserted into the hand piece before initial operation.

**WARNING**

Carefully follow all the safety precautions described in the paragraph 3.2.3 (BATTERY – Safety precautions).

**CAUTION**

NEVER place the hand piece in the charging unit without the battery inserted in the hand piece.


Insert the battery into the hand piece in the sequence described below:

- ① Remove the battery from its packaging.

**WARNING**

Never short circuit the terminals of the battery with any metal object as it could cause a risk of burns, fire and explosion.

- ② Use a flat tool (for example a coin) to unscrew counter-clockwise the battery cap of the hand piece (Fig. ① and ②).
- ③ An arrow is marked on the battery case (Fig. ③).

**NOTE**

This arrow indicates the orientation with which the battery must be inserted into the hand piece.

- ④ Slowly slide the battery with the arrow oriented toward the housing into the hand piece until it stops (Fig. ④).

**NOTE**

If the battery is inserted in the wrong way (i.e. in the opposite direction of the arrow), the hand piece will not function. In case the battery is inserted in the wrong way, no damage is caused to the hand piece's electronics.

**WARNING**

Use original Kulzer batteries only. The use of other manufacturer batteries or non-rechargeable batteries / primary batteries is a potential hazard and can cause irreversible damage to the unit.

- 5 Use a flat tool (for example a coin) to screw clockwise the battery cap of the hand piece.

**NOTE**

The unit will not function properly if the cap of the hand piece is not screwed on completely.

5.4 Inserting the fibre optic into the hand piece

- Rotate the light guide slightly (1) and insert it into the hand piece (2) until final positioning is reached. The final position is reached when you hear a “click”.
- Mount the protective light cone (3) on the top end of the light guide.

Positioning the light guide

- Rotate the light guide into the desired position for the polymerisation.
- To make full use of the light intensity provided, place the light guide as close to the composite as possible.

Avoid direct contact with the composite material!

**CAUTION – Keep the light guide clean at all times to obtain full light intensity.**

The light intensity is considerably reduced if:

- The light guide is not in perfect condition (scratched or chipped).
- The light guide is not correctly placed in the hand piece.
- The light guide has composite stains.

**WARNING**

A damaged light guide reduces the light intensity and must be replaced with a new one immediately. Sharp edges may cause serious injury. Use original Kulzer light guides only.

**WARNING – Infection control:**

To make treatment as safe as possible for patients and operators, the light guide and the protective light cone must be cleaned, disinfected and sterilised before each treatment. Follow the instructions listed in paragraph 8 (Cleaning, disinfection and sterilisation) step by step.

5.5 Battery charging**CAUTION – NEW Battery: first charging**

The Translux 2Wave battery has been partially charged at the factory. Therefore the battery must be fully charged before using the Translux 2Wave for the first time.

**NOTE – NEW Battery: first charging**

For the first charge, this process will take about **3 hours**. New batteries or batteries stored for a long time may take longer to charge. The battery will reach its full capacity after a few full charge / discharge cycles.

**WARNING – RISK OF EXPLOSION**

Use only the charging unit, supplied with the Translux 2Wave, to charge the battery. Never attempt to charge the battery using another battery charging unit or power source, due to risk of explosion and fire.

Insert the hand piece with the battery into the Translux 2Wave charging unit.

A beep confirms that the contacts are in the correct charging position!

**IMPORTANT NOTE – Alignment key**

The alignment key (5), located at the back bottom part of the hand piece's case, allows the correct insertion of the hand piece in the charging unit. When inserting the hand piece into the charging unit, the alignment key must be aligned with the slot positioned inside the charging unit.

Battery charging process.

The battery status is indicated by the battery icon at the bottom of the hand piece display.

a) Battery is completely empty

In this case, before starting the charging process, the electronics of the hand piece carry out a pre-qualification phase during which a few functioning parameters of the battery are checked out.

During the pre-qualification phase, that can last approximately from 10 to 30 minutes, the battery icon on the display shows one single black bar blinking (Fig. a).

Ending the pre-qualification phase, without detecting any battery defect, the “regular” battery charging process starts and all the four black bars inside the battery icon blink progressively (Fig. b).

IMPORTANT NOTE – Defective battery

If during the pre-qualification phase, the hand piece electronics detects a battery malfunction, the electronics interrupts the charging process and the frame of the battery icon (without internal black bars) will blink (Fig. c).

In this mode, a **beep** signal is emitted every second for a period of 30 seconds.
Replace the defective battery with a new one!

WARNING

Use Kulzer batteries only!

b) Battery is partially charged

In case of a battery being partially charged, (one or more internal black bars are displayed inside the battery icon) the charging process starts almost immediately after having inserted the hand piece into the charging station.

All the four black bars inside the battery icon blink progressively (see Fig. b).

c) Battery is fully charged


When the charging phase is finished (battery fully charged) all the four black bars inside the battery icon are permanently displayed (Fig. d).

NOTE

When the hand piece is placed in the charging unit its functioning is inactive. The operativeness of the hand piece buttons are prevented.

The hand piece can be removed from the charging unit during the charging process; however, the battery may not be fully charged and may shorten the operational duration.

GB

- 15 -

6 Operation

WARNING – Check the condition of the device before starting the treatment.

Before each treatment, always make sure that the device is in proper working order and that the accessories (light guide and protective light cone) are efficient. If anything is noted during treatment, do not carry out the treatment and contact Kulzer or a local authorised Customer Service agent (see paragraph 14 Service). Do not use the device and its accessories if damage is apparent. Prior to each use of the Translux 2Wave ensure that the emitted light intensity is sufficient to guarantee the polymerisation. Check the light intensity by using the light meter integrated into the charging unit (see paragraph 6.3 Measurement of the Light Intensity).

WARNING – Control of infections

First use

*The accessories, light guide and protective light cone, are delivered in **NON STERILE** condition and must be sterilised before being used for the first time.*

Every use

Once used, the accessories, light guide and protective light cone, must be reprocessed prior to reuse, according to the procedures described in the paragraph 8 (Cleaning, disinfection and sterilisation).

6.1 Selecting the curing time

The Translux 2Wave can operate with four different exposure cycles:

- **Exposure with slowly increasing output: “Soft-Start”** (cycle duration 20 seconds).
The “Soft-Start” means increase of light intensity from 50% up to 100% within 2 seconds.
- **Exposure with constant output:** cycle duration 20 seconds, 10 seconds or 5 seconds.

Press the “Time / Program” button  2 seconds to activate selection and scroll through the curing modes on the display.

The selected exposure time is indicated (in the display) by a round spot next to the selection.

NOTE

When the “Soft-Start” program is selected, round spots simultaneously lighten next to each icon  and  :

Each time the “Time / Program” button is briefly pressed, the display setting advances to the next lower exposure curing mode in a loop.

NOTE



During the light emission the “Time / Program” button is inactive.

NOTE – “Freeze Mode”

10 seconds after selecting an exposure mode the “Time / Program” button is deactivated.

Press the “Time / Program” button 2 seconds to activate selection and scroll through the curing modes on the display.

GB

- 16 -

6.2 Activating and deactivating the exposure cycle

Briefly press the “STAND-BY” button . The light emission will turn on.

Once the curing cycle is complete, the unit will automatically turn the light off.

Acoustical signal:

The following acoustical signals (beep) are emitted:

- 1 beep while press the “STAND-BY” button.
- 1 beep after 10 seconds of exposure time (only 20 second Mode).
- 1 beep at the end of the exposure time.

NOTE

The exposure cycle can be interrupted at any time, before the exposure time is over, by pressing again the “STAND-BY” button.

Acoustical signal:

An acoustical signal (beep) is emitted, when the “STAND-BY” button is pressed to interrupt the cycle.

6.3 Measurement of the Light Intensity

The charging unit of the Translux 2Wave is equipped with an integrated light intensity meter (8).

CAUTION

Before each treatment always check the light intensity by using the light intensity meter integrated in the Translux 2Wave charging unit.

- Without applying pressure, place the light guide flat and directly on the testing area (8).
- ▶ Activate the unit by pressing the “STAND-BY”  button (whichever exposure time selected).

The light intensity measuring is immediately indicated by the 2-colour LED (Symbol ) placed at the bottom right of the testing area.

GREEN LED: The light intensity is sufficient for the treatment.

YELLOW LED: The light intensity is insufficient and patient/s must not be treated.


CAUTION

If the light intensity is insufficient, please carry out the following controls and measures:

- The light guide has not been inserted into the hand piece (see paragraph 5.4 Inserting the fibre optic into the hand piece).
- Check the light guide for contamination or defects.
- Clean the light guide of any contamination (see paragraph 8 Cleaning, disinfection and sterilisation) or replace the damaged light guide with a new one.

Use original Kulzer light guides only.


WARNING

If the measures described above do not result in any improvement, do not use the device for any treatment. Switch off the unit (disconnect the charging unit from the electrical power source) and ensure that unauthorised unit re-activation is not possible. Contact a Kulzer's service partner, or technically trained personnel authorised by Kulzer (see paragraph 14 Service).

6.4 Sleep Mode

The hand piece is equipped with a "Sleep Mode" to minimise the unit's energy consumption. Outside the charging unit, all internal functions of the hand piece are automatically turned off ("Sleep Mode") if it is not used for approximately 5 minutes.

In "Sleep Mode" the hand piece display is in the condition shown in Fig. e:

- To terminate the "Sleep Mode" press the "STAND-BY" button . The hand piece is ready for operation.

NOTE

Terminated the "Sleep Mode", the hand piece and the display return to the latest setting used.

6.5 Forced Sleep Mode

Outside the charging unit it is possible to switch the hand piece to the "Sleep Mode" at any time.

- To activate the "Forced Sleep Mode" press and keep pressed for 5 seconds the "Time / Program" button . The hand piece and the display switch to the same conditions described in the above paragraph 6.4 (Sleep Mode).
- To terminate the "Sleep Mode" press the "STAND-BY" button . The hand piece is ready for operation.

NOTE – Memory Function

By terminating the "Sleep Mode", the hand piece and the display return to the latest setting used.

6.6 Additional information on the unit operations

NOTE

If the hand piece is placed in the charging unit before the curing cycle is terminated, the light emission automatically switches off.

When the hand piece is placed in the charging unit its functioning is inactive. The operativeness of the hand piece buttons is prevented.

During the light emission the "Time / Program" button is inactive.

The Translux 2Wave is equipped with microprocessors which continuously check the optimum charging parameters of the battery. After every treatment, the hand piece should therefore be inserted in the charging unit for storage irrespective of the battery level.


7 Diagnostic signals

7.1 Faulty battery

In the case of a defective battery, the battery icon, without internal black bars, continuously blinks (Fig. f).

A "beep" sound is emitted every second for a period of 30 seconds.

After 30 seconds only the battery icon frame remains blinking.


IMPORTANT

This error condition is detected only when the hand piece is placed in the charging unit.

7.2 Faulty LED

When a malfunction LED condition is detected by the microprocessor, the word "LED" blinks at the top left of the display (Fig. g).

7.3 Overheat protection

If multiple curing cycles in rapid succession occur (normally after 18 consecutive cycles, each lasting 20 seconds), the built-in temperature overload sensor may activate.

Once the overheat protection activates, an acoustical signal (**3 beeps**) is emitted and the letter "T" comes on, steady, at the right top of the display (**Fig. h**).

The activation of the overheat protection temporarily stops the use of the unit for a few minutes.

Allow the hand piece to cool down for at least 4 – 5 minutes, then start the next exposure by pressing the "STAND-BY" button .

NOTE

The hand piece will automatically exit from the cooling state upon reaching operational temperature and the letter "T" on the display switches off.

7.4 BUILT-IN Temperature overload sensor defective

In case of a defective temperature overload sensor the letter "T" comes on (**blinking**) at the right top of the display (**Fig. k**).

NOTE

If the temperature overload sensor fault occurs during an exposure cycle, the software allows completing the exposure. No further exposure cycle is allowed by pressing the "STAND-BY" button  and "4 beeps" are emitted.

Any further operation with defect temperature overload sensor is prevented.

7.5 Battery-low signal

When after frequent use the battery charge drops to the minimum level, the Translux 2Wave microprocessor permits another few cycles in any curing mode, without re-charging the battery.

This condition is signalled at the end of each cycle by "2 beeps".

At the end of a few cycles the remaining battery charge does not allow any further exposure cycle and this condition is shown on the display by the battery icon without internal black bars (**Fig. m**).

Charge the battery!

8 Cleaning, disinfection and sterilisation

DISCONNECT FROM MAINS

Before carrying out the cleaning and disinfection procedures disconnect the charging base from the electrical power source!

WARNING

Control of infections

*The light guide and the protective light cone are delivered in **NON STERILE** condition and must be cleaned, disinfected and sterilised (steam) before being used for the first time as well as before every use.*

Do not sterilise the hand piece and the charging unit!

Do not sterilise the battery!

8.1 Disassembly

 Remove the protective light cone from the light guide.

 Remove the light guide: Rotate the light guide slightly and pull it out from the hand piece.

8.2 LIGHT GUIDE – Cleaning and disinfection

CAUTION

Do not use sharp or pointed instruments to clean the light guide, since they could scratch the surface of the light guide and thus decrease light transmission.

Do not use cleaning / disinfectant agents containing organic, mineral and oxidising acids, hydrogen peroxide, chloride, iodine, bromide and solvents.

- 1) Immediately after use, any possible residue, such as composite material, must be removed from the end surface of the light guide. Any polymerised composite should be removed by using ethanol and a plastic spatula or soft nylon brush may help in removing the material residue.
- 2) Disinfect the light guide with a clean, soft lint free cloth dampened with a mild disinfectant solution, featuring a mild pH (pH 7), according to the manufacturer's specifications. **Allow the disinfectant solution to dry.**

8.2.1 LIGHT GUIDE – Inspection before the sterilisation

Before proceeding with the sterilisation process, carefully check the light guide for damaged surfaces, discolouration and contamination.

WARNING

Do not use a damaged light guide. If damage is apparent, replace the light guide with a new one.
If the light guide is still contaminated, repeat the cleaning and disinfection procedure.

8.3 PROTECTIVE LIGHT CONE – Cleaning and disinfection

CAUTION

Do not use cleaning / disinfectant agents containing organic, mineral and oxidising acids, hydrogen peroxide, chloride, iodine, bromide and solvents.

- ▶ Clean and disinfect the protective light cone with a soft cloth wetted with a mild disinfectant solution, featuring a mild pH (pH 7), according to the manufacturer's specifications. **Allow the disinfectant solution to dry.**

8.3.1 PROTECTIVE LIGHT CONE – Inspection before the sterilisation

Before proceeding with the sterilisation process, carefully check the protective light cone for damaged surfaces, veins discolouration and contamination.

WARNING

Do not use a damaged protective light cone. If damage is apparent, replace the protective light cone with a new one.

A damaged protective light cone cannot provide enough protection from the unit light.

If the protective light cone is still contaminated, repeat the cleaning and disinfection.

8.4 LIGHT GUIDE and PROTECTIVE LIGHT CONE – Packaging before the sterilisation

The light guide and the protective light cone can be sterilised by using a single-use medical steam sterilisation pouch of the appropriate size.

CAUTION

The light guide and the protective light cone must be packed individually in a single pouch.

When sterilising multiple instruments in one autoclave, ensure that the autoclave maximum load is not exceeded.

WARNING

Ensure that the inner pouch is large enough to contain the single instruments, without stressing the seals or tearing the packaging.

8.4.1 LIGHT GUIDE and PROTECTIVE LIGHT CONE – Sterilisation

8.4.1.1 Sterilisation method

- Carry out the sterilisation by means of a pre-vacuum steam autoclaving.

CAUTION

Never use any other sterilisation method, because of potential incompatibility with the materials with which is made the light guide.

DO NOT USE the following sterilisation methods: Ethylene oxide sterilisation, hot air sterilisation, flash autoclaving, STERRAD sterilisation, STERIS system or comparable sterilisation systems.

DO NOT sterilise the light guide and the protective light cone by using: Hydrogen peroxide, peracetic acid system, formaldehyde and glutaraldehyde sterilisation.

8.4.1.2 Sterilisation parameters

- Minimum sterilisation parameters that provide a 10^{-6} sterility assurance level – (SAL).

Cycle Type	Minimum Temperature	Minimum Exposure Time (pouch)	Minimum Dry Time
3 times Pre-vacuum (min. pressure 60 mBar)	132°C + 3°C (270°F + 5°F)	4 minutes	20 minutes

CAUTION


Do not exceed 135°C (275°F).

8.5 HAND PIECE – Cleaning and disinfection

8.5.1 Cleaning Mode

The “Cleaning Mode” prevents the accidental activation of the light exposure during the cleaning and the disinfection of the hand piece.

To enter the “Cleaning Mode” follow the procedures described here below:

8.5.2 Cleaning and disinfection

WARNING

The hand piece is not protected against the penetration of liquids.

Do not spray liquids directly onto the surface and charging contacts of the hand piece.

Cleaning and disinfectant agents must not enter the hand piece!


NEVER immerse the hand piece in liquids.

NEVER clean the hand piece under running water.

Do not sterilise the hand piece. Risk of serious damage, electrical shock and fire!

- Clean and wipe the hand piece casing with a clean, soft lint free cloth dampened with a mild surface disinfection agent, featuring a neutral pH (pH7) according to the manufacturer's specifications.

Allow the disinfectant solution to dry or dry the residual disinfectant on the hand piece casing with a soft lint free cloth.

CAUTION

Do not use solvents such as acetone, isopropyl alcohol, and hydrogen peroxide, chloride, iodine, bromide, phenol or other aggressive solutions to clean and disinfect the plastic surface of the hand piece. These substances may be harmful, causing discolouration and / or damage to the plastic materials of the hand piece.

Do not use abrasive cleaners!

IMPORTANT: Cleaning of the charging contacts of the hand piece

After having carried out the cleaning / disinfection of the hand piece's plastic casing, always clean the charging contacts which are positioned in the bottom part of the hand piece. The charging contacts of the hand piece must always be free from possible soil, composite and cleaning / disinfectant residues.

**IMPORTANT**

The hand piece charging contacts must be cleaned regularly (at least once a week) and always after having carried out the cleaning and disinfection of the hand piece plastic casing. Dirty charging contacts or contacts exposed to liquids, like cleaner or disinfectants, can negatively affect the charging process and prevent the battery of being charged.

▶ Clean the charging contacts of the hand piece with an alcohol dampened soft lint free cloth or cotton swab.

**CAUTION**

Do not use sharp or pointed instruments to clean the hand piece contacts, since they could scratch and damage their surface and thus worsen the conductivity with the charging contacts.

Before re-using the hand piece make sure that the hand piece's surface and the charging contacts are perfectly dried. If necessary, dry the contacts by blowing air onto them with compressed air.

8.6 CHARGING UNIT – Cleaning casing**DISCONNECT FROM MAINS**

Disconnect the input power supply from the power outlet and charging unit prior to cleaning / disinfection.

**WARNING**

The charging unit casing **is not protected** against the penetration of liquids.

Do not spray liquids directly onto the surface of the charging unit casing.

**CAUTION**

DO NOT sterilise the charging unit. It will no longer function and may cause a risk of serious damage, electrical shock and fire!

Do not use solvents such as acetone, isopropyl alcohol, and hydrogen peroxide, chloride, iodine, bromide, phenol or other aggressive solutions to clean the plastic surface of the charging unit. These substances may be harmful, causing discolouration and / or damage to the plastic materials of the charging unit.

Do not use abrasive cleaners!

▶ Clean and wipe the charging unit casing with a clean, soft, lint free cloth dampened with a mild surface cleaning agent, featuring a neutral pH (pH7) according to the manufacturer's specifications.

Dry the charging unit casing with a clean, non-abrasive cloth.

**CAUTION**

Be sure to dry the charging unit after cleaning. Moisture inside the unit could cause damage.

**IMPORTANT**

Upon completion of the cleaning procedure, make sure that the charging contacts inside the charging unit are perfectly dry. If necessary, dry the contacts by blowing air onto them with compressed air.

8.6.1 CHARGING UNIT – Cleaning charging contacts**IMPORTANT**

After having carried out the cleaning of the charger's plastic casing, always clean the charging unit contacts. The charging contacts of the charging unit must always be free from possible soil, composite and cleaning residues.

The charging unit contacts must be cleaned regularly (at least once a week) and, in any case, always after having carried out the cleaning of the plastic casing. Dirty contacts or contacts exposed to liquids, like cleaner or disinfectants, can negatively affect the charging process and prevent the battery of being charged.

▶ Clean the contacts of the charging unit with an alcohol dampened soft lint free cloth or cotton swab.

**CAUTION**

Do not use pointed or sharp objects to clean the charging unit contacts!

Do not spray or flood the contacts with liquids!

Do not bend or modify the shape of the charging unit contacts during cleaning.

Make sure that the charge contacts remain dry after cleaning and are not contacted by metallic and greasy parts. If necessary, dry the contacts by blowing air onto them with compressed air.

8.6.2 CHARGING UNIT – Replaceable charging contacts

Damaged or highly soiled charging unit contacts (that cannot be cleaned by following the procedures described in the above paragraph) can prevent the contacts' conductivity and consequently the battery charging.


In this case it is possible to replace the damaged charging unit contacts with new ones.

For additional information regarding replacement of charging unit contacts and spare part supply please contact your local Kulzer subsidiary.

Working instructions:

Required tool: screwdriver for Phillips-Screws, mid-size.

1. Remove hand piece out of charger base.
2. Disconnect charger from mains.
3. Unplug connecting wire from charger.
4. Turn charger upside down.
5. Remove 2 screws from the charging contact bridge (see rear fold-out cover, **point 12**).
6. Pull out charging contact bridge.
7. Insert new contact bridge and fix / connect in reverse way.

GB

CAUTION

*Do not bend or contaminate the new contact bridge!
Do care about perfect insert and fixation! Do not over tight the screws!*

- 22 -

9 Shelf storage

Storage of the hand piece with battery during extended periods of non-use. Before or after a long inactivity period, fully charge the battery or store it in the operating charging unit. Although the built-in safety circuit of the battery prevents the full discharge, it is strongly recommended to fully charge the battery at least once every 5 months after inactivity.

10 Disposal

NOTE – IMPORTANT

The device must not be disposed off in normal domestic waste.

For the disposal of spare parts, or the unit, please contact the Kulzer subsidiary directly in your country.

Dispose unserviceable batteries and light guides according to the corresponding legal requirements in your country.

WARNING – Damaged or exhausted battery disposal

Dispose a damaged or exhausted battery after taking the necessary steps to prevent external short circuit by the following method:

After having insulated the terminals of the battery with tape dispose off in accordance with your local council.

The conditions and precautionary measures for disposing of the device are subject to the valid legal provisions similar to any other electronic device which has become unusable.

Disposal of old equipment according to WEEE (European Directive) or Electrical and Electronic Equipment Act (German Law ElektroG).

11 Troubleshooting


WARNING

If the below steps do not solve the problem, contact Kulzer or a local authorised Customer Service agent: **DO NOT USE** the Translux 2Wave and **DO NOT ATTEMPT** to service the unit as it may cause serious injury.

If the device does not seem to be working properly, read this user manual again and then check the following table:

PROBLEM	POSSIBLE CAUSE	REMEDY	
The charging unit does not power on. The GREEN LED  on the charging unit is switched off.	The input power supply is not connected to electrical supply outlet and / or to the charging unit power socket.	Make sure that the input power supply is securely connected to electrical supply outlet and to the charging unit power socket.	
	Power outlet carries no voltage.	Use a different power outlet.	
	The interchangeable adapter is not properly inserted into the input power supply casing.	Carefully read paragraph 5.2 of this user manual.	
	The power cable of the input power supply is damaged.	Replace the input power supply with a new one. (The input power supply and its cord cannot be separated from each other.)	
	The input power supply is defective.	Replace the input power supply with a new one.	
	The charging unit is defective.	Contact Kulzer or local authorised Customer Service agent.	
Hand piece NOT inserted into the charging unit. The light emission does not start when the "STAND-BY" button  is pressed and no information appears on the hand piece's display.	Hand piece without battery.	Insert the battery into the hand piece (see paragraph 5.3).	
	Hand piece in "Sleep Mode".	Press the "STAND-BY" button to bring the hand piece out of the "Sleep Mode". NOTE: Memory Function Terminate the "Sleep Mode" in the hand piece and the display will return to the latest setting used. Make reference to paragraph 6.4.	
	Empty battery. There is not enough charge in the battery to turn on the hand piece and display.	Place the hand piece in the charging unit and re-charge the battery (see paragraph 5.5).	
	Electronics malfunction.	Contact Kulzer or local authorised Customer Service agent.	
	Defective battery.	To verify if the battery is defective, place the hand piece in the charging unit. A "prequalification sequence" starts in which Translux 2Wave micro-processor checks the battery status. This prequalification sequence can last up to max. 30 minutes. If, during this check sequence, a "beep" signal is emitted every second for a period of 30 seconds and the frame of the battery icon (without internal black bars) blinks, the battery is detected to be defective. Replace the battery with a new one (see paragraph 5.3).	
	"2 beeps" are emitted at the end of the exposure cycle and the battery icon on the display appears empty.	Low battery charge.	The "2 beeps" signal that the residual charge level of the battery is adequate only for a few subsequent exposure cycles. Place the hand piece in the charging unit and re-charge the battery (see paragraph 5.5).

PROBLEM	POSSIBLE CAUSE	REMEDY
<p>The light emission does not start when the “STAND-BY” button  is pressed or during an exposure cycle the light exposure is interrupted with a beep sound and the word “LED” blinks at the left top of the display.</p>  <p>blinking</p>	<p>The LED pack could be defective.</p>	<p>Verify the LED pack functionality by means of the following procedure: Insert hand piece into the charger and press the “STAND-BY” button consecutively for 5 times. Remove the hand piece from the charger unit. The word “LED” switches-off on the display, the hand piece has been reset. If word “LED” appears again blinking in display at start of the next exposure cycle and no light emission starts, the LED pack is defective. In this case please contact Kulzer or local authorised Customer Service agent.</p>
<p>After several consecutive exposure cycles, the light emission does not turn on, when the “STAND-BY” button  is pressed, and the letter “T” comes on (steady) at the top right of the display. “3 beep” signals are emitted.</p>  <p>steady</p>	<p>Overheating protection.</p> <p>The hand piece has become overheated in the course of use and the temperature overload sensor has activated.</p> <p>The activation of the overheat protection temporarily prevents (for a few minutes) the use of the unit.</p>	<p>Allow the hand piece to cool down for at least 4 – 5 minutes, then start the next exposure by pressing the “STAND-BY” button.</p> <p>Note: The hand piece will automatically exit from the cooling state upon re-aching operational temperature and the letter “T” on the display switches off.</p>
<p>During an exposure cycle, the light emission is interrupted and “3 beeps” are emitted. The letter “T” comes on (steady) at the right top of the display.</p>  <p>steady</p>	<p>Overheating protection.</p> <p>The hand piece has become overheated in the course of use and the temperature overload sensor has activated.</p> <p>The activation of the overheat protection temporarily prevents (for a few minutes) the use of the unit.</p>	<p>Allow the hand piece to cool down for at least 4 – 5 minutes, then start the next exposure by pressing the “STAND-BY” button.</p> <p>Note: When the cooling time is over, the letter “T” on the display switches off.</p>
<p>Malfunctions during charging phase. A “beep” signal is emitted every second for a period of 30 seconds and the frame of the battery icon (without internal marks) blinks.</p>  <p>blinking</p>	<p>During the battery pre-qualification phase (battery condition: completely empty ref. paragraph 4.2.5 and 5.5).</p> <p>The electronics detects a defective or malfunction status of the battery and interrupts the charging process.</p>	<p>Replace the battery with a new one (see paragraph 5.3).</p> 
<p>Hand piece placed in the charging unit but the charging phase does not start.</p> <p>The black bars in the battery icon does not cycle progressively.</p> 	<p>Defect: interrupt of charging contacts.</p> <p>Improper contact between the charging base and the hand piece.</p>	<p>See paragraph 8.5 and 8.6.2.</p> <p>Check if the hand piece has been inserted in the correct orientation into the charging unit.</p> <p>Alignment key. Make reference to the paragraph 4.2.1, points 4 and 5.</p>
	<p>Charging and / or hand piece contacts are soiled.</p>	<p>Clean the hand piece and the charging unit contacts. Make reference to the paragraphs 8.5 and 8.6.1.</p>
<p>The light intensity measured with the meter integrated in the charging unit is insufficient.</p> <p>The LED  on the charging unit lights YELLOW.</p>	<p>Light guide not properly inserted into the hand piece.</p>	<p>Check light guide for complete insertion in metal cone of the hand piece.</p>
	<p>Light guide output contaminated by foreign material or composite residues that can interfere with light output.</p>	<p>Clean the light guide output. Make reference to paragraph 8.2.</p>
	<p>Light guide damaged or aged.</p>	<p>Replace the light guide with a new one.</p>
<p>The four round spots of the display are blinking when the hand piece is placed into the charging unit.</p>  <p>blinking</p>	<p>Charging circuit inside the hand piece defective.</p>	<p>Contact Kulzer or local authorised Customer Service agent.</p>

PROBLEM	POSSIBLE CAUSE	REMEDY
	The built-in temperature overload sensor is defective.	Contact Kulzer or local authorised Customer Service agent. Note: If the temperature overload sensor fault occurs during an exposure cycle, the software allows completing the exposure. No further exposure cycle is allowed by pressing the "STAND-BY" button and "4 beeps" are emitted.

12 Technical specifications

Classification according to Directive 93/42/EEC: Class I (one)

Safety Standards:

The product has been tested and found to comply to IEC 60601-1 (second edition) and to IEC 60601-1-1.

The product has been tested and found to comply to IEC 60601-1 (third edition) +AMD1:2012 and to IEC 60601-1-2:2014 (fourth edition – for electro-magnetic compatibility – EMC).

The hand piece is intended for use within the patient environment whereas the charger unit and the power supply adapter are not intended for use within the patient environment.

Note: The patient environment is defined as an area 1.5 m from the patient as per IEC 60601-1 (third edition) +AMD1:2012.

Classification according to IEC 60601-1:

Class II, Type B, **Applied Part:** Light guide with protective cone. IP 20 (Charging unit). IP 20 (Translux 2Wave)

Charging unit:

Translux 2Wave Charging unit

Power supply for the charging unit:

Switching Power Supply (Direct Plug-In) manufactured by: Pihong Technology CO LTD – **Model:** PSAC05R-050L6.

Input: 100 – 240 V \sim – 300 mA – 50 / 60 Hz –.

Output: 5 V \equiv – 1 A Max – \ominus \oplus .

Class II equipment

Power supply for the hand piece:

Rechargeable Lithium-ion battery.

Manufacturer: Panasonic. **Model:** NCR-18500.

Nominal Voltage: 3.6 V.

Nominal capacity (Typical): 2000 mAh

Hand piece:

Model: Translux 2Wave.

Operation: Intermittent service: 120" ON – 40" OFF – max 2 times running

Light source:

High-luminosity LED.

LED Product Risk Group 2 (IEC 62471).

Wavelength range: 385 – 510 nm.

Dominant wavelength: Min. 385 nm – Max. 470 nm

Light guide:

Lead-free multi core rod, 8 mm diameter.

Sterilisable in pre-vacuum steam autoclave:

4 minutes (minimum) at minimum temperature of 132°C + 3°C (270°F + 5°F). Max. 500 cycles

Settable exposure cycles:

– "Soft-Start" exposure

20 seconds of exposure time. (The "Soft-Start" means increase of light intensity from 50% up to 100% within 2 seconds.)

Acoustical signals:

1 beep at the start of exposure.

1 beep after 10 seconds of exposure.

1 beep at the end of exposure

– Normal exposure

20, 10 or 5 seconds of exposure time.

Acoustical signals:

1 beep at the start of exposure.

1 beep after 10 seconds of exposure (only 20 s Mode).

1 beep at the end of exposure

Charging time of empty battery:
Environmental conditions of unit:

Approx. 3 hours

	Operation	Transportation and storage
Temperature	10°C – 35°C (50°F – 95°F)	-20°C – 40°C (-4°F – 104°F)
Relative humidity	45% – 85%	45% – 85%
Atmospheric pressure	800 – 1060 hPa	500 – 1060 hPa

Weight and dimensions:

Charging unit: weight 450 g.
Dimensions: L 140 x W 58 x H 62 mm
Hand piece: weight 150 g (incl. light guide).
Dimensions: L 275 mm, max. Ø 24 mm

12.1 Electromagnetic compatibility EN 60601-1-2

WARNING

The device requires specific EMC precautions and must be installed and started up in accordance with the EMC information given in this paragraph.

Portable and mobile radio communication appliances may affect the correct function of the device.

GB

- 26 -

Guidance and manufacturer's declaration – Electromagnetic emissions		
The Translux 2Wave is intended for use in the electromagnetic environment specified below. The customer or the user of the Translux 2Wave should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – Guidance
RF emissions CISPR 11	Group 1	The Translux 2Wave uses RF energy only for its internal function. Therefore it's RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Translux 2Wave is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – Electromagnetic immunity

Enclosure port

The Translux 2Wave device is designed to work in the electromagnetic environment specified below. The client or user of the Translux 2Wave device should always make sure that it is used only in such environment.

Phenomenon	Basic EMC standard or test method	Immunity test levels	Electromagnetic environment – Guidance
Electrostatic discharge (ESD)	IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors must be wood, concrete or ceramic. If the floor is lined with synthetic materials, the relative humidity should be at least 30 %
Radiated RF EM fields ^{a)}	IEC 61000-4-3	3 V/m ^{f)} 80 MHz – 2,7 GHz ^{b)} 80 % AM at 1 kHz ^{c)}	Portable or mobile RF communication device must not be used in close proximity of the product, including its cables, except when these respect the distances of separation recommended and calculated from the equation applicable to the frequency of the transmitter.
Rated power frequency magnetic fields ^{d) e)}	IEC 61000-4-8	30 A/m ^{g)} 50 Hz or 60 Hz	Grid frequency magnetic fields should have levels characteristic of a typical commercial or hospital environment.

- a) The interface between the PATIENT physiological signal simulation, if used, and the Translux 2Wave shall be located within 0,1 m of the vertical plane of the uniform field area in one orientation of the Translux 2Wave.
- b) Translux 2Wave that intentionally receive RF electromagnetic energy for the purpose of their operation shall be tested at the frequency of reception. Testing may be performed at other modulation frequencies identified by the RISK MANAGEMENT PROCESS. This test assesses the BASIC SAFETY and ESSENTIAL PERFORMANCE of an intentional receiver when an ambient signal is in the passband. It is understood that the receiver might not achieve normal reception during the test.
- c) Testing may be performed at other modulation frequencies identified by the RISK MANAGEMENT PROCESS.
- d) Applies only to Translux 2Wave with magnetically sensitive components or circuitry.
- e) During the test, the Translux 2Wave may be powered at any NOMINAL input voltage, but with the same frequency as the test signal.
- f) Before modulation is applied.
- g) This test level assumes a minimum distance between the Translux 2Wave and sources of power frequency magnetic field of at least 15 cm. If the RISK ANALYSIS shows that the Translux 2Wave will be used close than 15 cm to sources of power frequency magnetic field, the IMMUNITY TEST LEVEL shall be adjusted as appropriate for the minimum expected distance.

Guidance and manufacturer's declaration – Electromagnetic immunity

Input a.c. power port

The Translux 2Wave device is designed to work in the electromagnetic environment specified below. The client or user of the Translux 2Wave device should always make sure that it is used only in such environment.

Phenomenon	Basic EMC standard or test method	Immunity test levels	Electromagnetic environment – Guidance
Electrical fast transients / bursts ^{a) i) o)}	IEC 61000-4-4	±2 kV contact 100 KHz repetition frequency	The quality of the power supply should be that of a typical commercial or hospital environment.
Surges Line-to-line ^{a) b) j) o)}	IEC 61000-4-5	± 0,5 kV, ± 1 kV	The quality of the power supply should be that of a typical commercial or hospital environment.
Surges Line-to-ground ^{a) b) j) k) o)}	IEC 61000-4-5	± 0,5 kV, ± 1 kV, ± 2kV	The quality of the power supply should be that of a typical commercial or hospital environment.
Conducted disturbances induce by RF fields ^{c) d) o)}	IEC 61000-4-6	3 V ^{m)} 0,15 MHz - 80 MHz 6 V ^{m)} in ISM bands between 0,15 MHz and 80 MHz ⁿ⁾ 80 % AM at 1 KHz ^{o)}	Portable or mobile RF communication device must not be used in close proximity of the product, including its cables, except when these respect the distances of separation recommended and calculated from the equation applicable to the frequency of the transmitter.
Voltage dips ^{g) p) r)}	IEC 61000-4-11	0% U _T ; 0,5 cycle ^{q)} At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % U _T ; 1 cycle and 70 % U _T ; 25/30 cycles ^{h)} Single phase: at 0°	The quality of the power supply should be that of a typical commercial or hospital environment.
Voltage interruptions ^{h) i) o) r)}	IEC 61000-4-11	0% U _T ; 250/300 cycle ^{h)}	The quality of the power supply should be that of a typical commercial or hospital environment.

a) The test may be performed at any one power input voltage within the Translux 2Wave RATED voltage range. If the Translux 2Wave is tested at one power input voltage, it is not necessary to re-test at additional voltages.

b) All Translux 2Wave cables are attached during the test.

c) Calibration for current injection clamps shall be performed in a 150 Ω system.

d) If the frequency stepping skips over an ISM or amateur band, as applicable, an additional test frequency shall be used in the ISM or amateur radio band. This applies to each ISM and amateur radio band within the specified frequency range.

e) Testing may be performed at other modulation frequencies identified by the RISK MANAGEMENT PROCESS.

f) Translux 2Wave with a d.c. power input intended for use with a.c.-to-d.c. converters shall be tested using a converter that meets the specifications of the MANUFACTURER of the Translux 2Wave. The IMMUNITY TEST LEVELS are applied to the a.c. power input of the converter .

g) Applicable only to Translux 2Wave connected to single-phase a.c. mains.

h) E.g. 10/12 means 10 periods at 50 Hz or 12 periods at 60 Hz.

i) Translux 2Wave with RATED input current greater than 16 A / phase shall be interrupted once for 250/300 cycles at any angle and at all phases at the same time (if applicable). Translux 2Wave with battery backup shall resume line power operation after the test. For Translux 2Wave with RATED input current not exceeding 16 A, all phases shall be interrupted simultaneously.

j) Translux 2Wave that do not have a surge protection device in the primary power circuit may be tested only at ± 2 kV line(s) to earth and ± 1 kV line(s) to line(s).

k) Not applicable to CLASS II Translux 2Wave.

l) Direct coupling shall be used.

m) r.m.s. , before modulation is applied.

n) The ISM (industrial scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0, 15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.

o) Applicable to Translux 2Wave with RATED input current less than or equal to 16 A / phase and Translux 2Wave with RATED input current greater than 16 A / phase.

p) Applicable to Translux 2Wave with RATED input current less than or equal to 16 A / phase.

q) At some phase angles, applying this test to Translux 2Wave with transformer mains power input might cause an overcurrent protection device to open. This can occur due to magnetic flux saturation of the transformer core after the voltage dip. If this occurs, the Translux 2Wave shall provide BASIC SAFETY during and after the test.

r) For Translux 2Wave that have multiple voltage settings or auto ranging voltage capability, the test shall be performed at the minimum and maximum RATED input voltage. Translux 2Wave with a RATED input voltage range of less than 25 % of the highest RATED input voltage shall be tested at one RATED input voltage within the range.

Guidance and manufacturer's declaration – Electromagnetic immunity			
Patient coupling port The Translux 2Wave device is designed to work in the electromagnetic environment specified below. The client or user of the Translux 2Wave device should always make sure that it is used only in such environment.			
Phenomenon	Basic EMC standard or test method	Immunity test levels	Electromagnetic environment – Guidance
Electrostatic discharge (ESD) ^{c)}	IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors must be wood, concrete or ceramic. If the floor is lined with synthetic materials, the relative humidity should be at least 30 %
Conducted disturbances induce by RF fields ^{a)}	IEC 61000-4-6	3 V ^{b)} 0,15 MHz - 80 MHz 6 V ^{b)} in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 KHz	Portable or mobile RF communication device must not be used in close proximity of the product, including its cables, except when these respect the distances of separation recommended and calculated from the equation applicable to the frequency of the transmitter.

a) The following apply:

- All PATIENT-COUPLED cables shall be tested, either individually or bundled
- PATIENT-COUPLED cables shall be tested using a current clamp unless a current clamp is not suitable. In cases where a current clamp is not suitable, an EM clamp shall be used.
- No intentional decoupling device shall be used between the injection point and the PATIENT COUPLING POINT in any case.
- Testing may be performed at other modulation frequencies identified by the RISK MANAGEMENT PROCESS.
- Tubes that are intentionally filled with conductive liquids and intended to be connected to a PATIENT shall be considered to be PATIENT-COUPLED cables.
- If the frequency stepping skips over an ISM or amateur radio band, as applicable, an additional test frequency shall be used in the ISM or amateur radio band. This applies to each ISM and amateur radio band within the specified frequency range.
- The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz; 3,5 MHz to 4,0 MHz; 5,3 MHz to 5,4 MHz; 7 MHz to 7,3 MHz; 10,1 MHz to 10,15 MHz; 14 MHz to 14,2 MHz; 18,07 MHz to 18,17 MHz; 21,0 MHz to 21,4 MHz; 24,89 MHz to 24,99 MHz; 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.

b) r.m.s., before modulation is applied

c) Discharges shall be applied with no connection to an artificial hand and no connection to PATIENT simulation. PATIENT simulation may be connected after the test as needed in order to verify BASIC SAFETY and ESSENTIAL PERFORMANCE.

Guidance and manufacturer's declaration – Electromagnetic immunity**Signal input / output parts port**

The Translux 2Wave device is designed to work in the electromagnetic environment specified below.
The client or user of the Translux 2Wave device should always make sure that it is used only in such environment.

Phenomenon	Basic EMC standard or test method	Immunity test levels	Electromagnetic environment – Guidance
Electrostatic discharge (ESD) ^{e)}	IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors must be wood, concrete or ceramic. If the floor is lined with synthetic materials, the relative humidity should be at least 30 %
Electrical fast transients / bursts ^{b) h)}	IEC 61000-4-4	±1 kV contact 100 KHz repetition frequency	The quality of the power supply should be that of a typical commercial or hospital environment.
Surges Line-to-ground ^{a)}	IEC 61000-4-5	± 2kV	The quality of the power supply should be that of a typical commercial or hospital environment.
Conducted disturbances induce by RF fields ^{b) d) g)}	IEC 61000-4-6	3 V ^{h)} 0,15 MHz - 80 MHz 6 V ^{h)} in ISM bands between 0,15 MHz and 80 MHz ⁱ⁾ 80 % AM at 1 KHz ⁱ⁾	Portable or mobile RF communication device must not be used in close proximity of the product, including its cables, except when these respect the distances of separation recommended and calculated from the equation applicable to the frequency of the transmitter.

a) This test applies only to output lines intended to connect directly to outdoor cables.

b) SIP/SOPS whose maximum cable length is less than 3 m in length are excluded.

c) Testing may be performed at other modulation frequencies identified by the RISK MANAGEMENT PROCESS.

d) Calibration for current injection clamps shall be performed in a 150 Ω system.

e) Connectors shall be tested per 8.3.2 and Table 4 of IEC 61000-4-2:2008. For insulated connector shells, perform air discharge testing to the connector shell and the pins using the rounded tip finger of the ESD generator, with the exception that the only connector pins that are tested are those that can be contacted or touched, under conditions of INTENDED USE, by the standard test finger shown in Figure 6 of the general standard, applied in a bent or straight position.

f) Capacitive coupling shall be used.

g) If the frequency stepping skips over an ISM or amateur radio band, as applicable, an additional test frequency shall be used in the ISM or amateur radio band. This applies to each ISM and amateur radio band within the specified frequency range.

h) r.m.s., before modulation is applied.

i) The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.

Test specifications for enclosure port immunity to RF wireless communications equipment

The Translux 2Wave device is designed to work in an electromagnetic environment in which RF radiated disturbances are kept under control. The purchaser or user of the Translux 2Wave device can contribute to prevent electromagnetic interference by ensuring a minimum distance between the RF mobile and cordless communication device (transmitters) and the Translux 2Wave device, as recommended below, in relation to the maximum output power of the radio-communication devices.

Test frequency (MHz)	Band ^{a)} (MHz)	Service ^{a)}	Modulation ^{b)}	Maximum power (W)	Distance (m)	Immunity test level (V/m)
385	380–390	TETRA 400	Pulse modulation ^{b)} 18 Hz	1,8	0,3	27
450	430–470	GMRS 460, FRS 460	FM ^{c)} ± 5 kHz deviation 1 kHz sine	2	0,3	28
710	704–787	LTE Band 13, 17	Pulse modulation ^{b)} 217 Hz	0,2	0,3	9
745						
780						
810	800–960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation ^{b)} 18 Hz	2	0,3	28
870						
930						
1720	1700–1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation ^{b)} 217 Hz	2	0,3	28
1845						
1970						
2450	2400–2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ^{b)} 217 Hz	2	0,3	28
5420	5100–5800	WLAN 802.11 a/n	Pulse modulation ^{b)} 217 Hz	0,2	0,3	9
5500						
5785						

NOTE If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the Translux 2Wave may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

- a) For some services, only the uplink frequencies are included.
b) The carrier shall be modulated using a 50% duty cycle square 2Wave signal.
c) As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

**WARNING:**

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Translux 2Wave device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

12.2 Charging unit identification plate



Description of the symbols on the identification plate

- 1 CE mark.
- 2 Caution, consult accompanying documents.
- 3 Consult instructions for use.
- 4 Double insulation (device complies with safety class II).
- 5 Disposal of old electrical and electronic equipment.
(Applicable in the European Community and other European countries with separate collection systems.)
- 6 Nemko mark.

12.3 Hand piece serial number

The hand piece serial number is engraved on the bottom part of its casing.

Engraved
serial number.

GB

13 Warranty

Each Kulzer unit is subject to a stringent check and final inspection prior to delivery to guarantee perfect operation of the equipment. For all new Kulzer products purchased from a Kulzer retailer or importer, Kulzer offers a warranty of 2 (two) years from the day of purchase for material and production faults. During the warranty period, Kulzer undertakes to repair (or replace, if it chooses) parts of the products deemed by Kulzer to be faulty, free of charge. No provision has been made for complete replacement of the Kulzer products.

Under no circumstances shall Kulzer be liable for any direct or indirect personal injury or damage in the following cases:

- The unit was not used for its intended purpose.
- The unit was not used according to the instructions and specifications described in this manual.
- The electrical installation in the rooms in which the unit is used does not comply with existing legislation and relevant specifications.
- The assembly, additions, setting, alteration and repair was not performed by Kulzer specialist personnel.
- The conditions for storing the unit do not comply with the specifications listed in paragraph 12 (Technical specifications).

Excluded from the warranty are transportation damage, damage caused by improper use or negligence, damage caused by connection to a power supply other than that specified for the equipment and damage to the control lamps, buttons and all accessories. The warranty shall become invalid if the unit is altered or repaired by personnel not belonging to Customer Service. When making a warranty claim, the customer forwards the damaged unit, at his own cost, to the retailer or Kulzer importer where he bought the unit. The unit, including the accessories, must be well packed before shipping (if possible, in the original packaging).

A card containing the following information must be enclosed with the unit if returning the unit:

- a) Name, address and telephone number of owner.
- b) Name and address of retailer / importer.
- c) Photocopy of the delivery note / purchase receipt for the unit by the owner which shows the date, designation and serial number of the unit.
- d) Description of malfunction.

The warranty does not cover transportation or any transportation damage.

In the event of damage caused by an accident or improper use, or which occurs after expiration of the warranty period, repairs will be invoiced according to effective material and labour costs.

14 Service

We appreciate your comments, feedback, and suggestions.

Contact in the countries and more information are available by the displayed QR code or at our website www.kulzer.com

15 Document history

- | | |
|---------|--|
| 2012-08 | Draft version for Prototypes. |
| 2013-09 | First edition. |
| 2014-03 | Editorial changes paragraph 2.1, 3.2.1, 10, 11, 12 and 12.2. |
| 2015-05 | Editorial changes paragraph 14.1, fold-out cover new fig. battery. |
| 2017-07 | Document version 11. Changed name manufacturers, new address and logo of the manufacturer.
Editorial revision paragraph 14. |
| 2017-08 | Editorial changes. |

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