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

Instructions for use:

1. Indications

The reinforcement mesh nets or grids are used to strengthen removable dentures. The reinforcement nets are placed in the thin area of the denture to prevent breakage when the denture is subjected to strong pressure. The reinforcement nets are available as prefabricated elements which can be customized to the required size by the dental technician when processing the denture. The reinforcement nets are products which are only used in conjunction with a dental prosthesis. The reinforcement nets and the prosthesis together compensate for the lost teeth in the patient's mouth.

When used correctly, the reinforcement nets are completely covered by denture acrylic (with the exception of item no. 232 1100), so that there is no direct contact with the soft tissue, unless the product has been used incorrectly or in the event of damage.

2. Contraindications

The items are to be used exclusively for the purpose of their application. Contraindications are not known. Material tolerances see point 5.

3. Directions for use

Item nos. 222 / 223 / 232 can be easily cut to shape using suitable scissors. Please do not cut these fine-meshed nets with a separating disc, in order to avoid possible oxidation. It is necessary to cut item no. 235 with a separating disc due to this product's higher material thickness. We recommend using the Renfert separating disc, item no. 570322. Take care not to apply too much pressure and only use a cutting speed of 15,000 1/min. Possible discoloration (oxidation) on the processed parts must be removed with a suitable rubber polisher. In some unfavorable situations (exposed metal, cracks in the acrylic) it is possible that chloride induced corrosion may occur. It is essential that these areas are opened, cleaned and polished, then resealed again. Before the finished restoration is inserted into the patient's mouth, it should be visually checked to ensure it is in good condition. Processing residues can be removed by cleaning.

The nets and mesh grids are not sterile. When using disinfection liquid it is important to test the grids for compatibility prior to use.

4. Hazard Warning

Sharp edges may occur when custom fitting the nets/grids. This could cause problems in the finished denture (e.g. cracks/flaws), but may also be an injury hazard for the patient, dentist or dental technician. Sharp edges must always be rounded off with a suitable rubber polisher. It is advisable to test the compatibility of the mesh net and denture acrylic or opaque prior to actual use.

It is essential to avoid a gap between the gold edge and the denture acrylic, as residual food could collect in the crevice and lead to infection within the mouth. Areas of net or mesh which has become exposed due to surface damage during processing, could cause ions from the base material to be set free into the body. This can lead to an allergic reaction (see 5) or headaches. For this reason, this information must be passed on to the dentist.

5. Allergic reaction

Constituents from the net or mesh could cause an allergic reaction.

For this reason it is essential to check the alloy constituents before use:

Item nos. 222/223 contain X5CrNiMo 17-12-2 (1.4401), with AuCo coating (max. 0,3% Co).

Item no. 2350004 contains X5CrNi 18-10 (1.4301).

Item no. 2350104 contains X5CrNi 18-10 (1.4301), with AuCo coating (max. 0,3% Co).

Item no. 232 contains X5CrNiMo 17-12-2 (1.4401) and X5CrNi 18-10 (1.4301), with AuCo coating (max. 0,3% Co).

All non-precious steel alloys contain nickel!

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