WAX Disc

Directions for Use

Mode d'emploi
Istruzioni per l'uso
Gebrauchsanweisung
Instrucciones de uso

 Responsible Manufacturer:
DENTSPLY Prosthetics
570 West College Avenue
York, PA 17401 U.S.A.
1-800-243-1942
prosthetics.dentsply.com

DP-0000177 Rev.1 (09/16)

THE DENTAL
SOLUTIONS
COMPANY™

GB

Technical specifications:

Description	Requirement ISO 15854:2005	Dentsply Sirona Wax Blank
Burning residual (1h at 700°C (1292°F))	< 0.1%	<0.05%

Indications:

Dentsply Sirona Wax discs are indicated for milling patterns from which custom made dental prostheses, such as crowns and bridges are cast.

Contraindications:

- > Waxing-up technique

Warnings:

Possible cross-reactions or interactions of this product with other products or material already present in the oral environment must be taken into consideration by the professional when selecting this product.

Precautions:

- > Keep product dust away from eyes
- > Do not inhale dust particles during grinding
- > Use local vacuum suction and suitable mouth/face protection during manual machining at the workplace
- > Keep away from heat and light
- > Maximum storage temperature: 35°C (95°F)
- > The product colors may vary due to the use of natural wax

The safety and warning notes listed here describe how to use our product in a safe and risk-free manner. Notify the dentist in charge of all factors described above if you use this product for a custom design and make sure to comply with the pertinent Safety Data Sheets (SDS)

Symbols on product labels:

	Product code or reorder number
	Batch or Lot number
	Expiration date
	Date of manufacture
	Consult instructions for use

Use Instructions:

Dentsply Sirona Wax should be fixed in the milling machine according to the instructions from the device manufacturers.

Important Note: Avoid overheating during machining and finishing.

Recommended operating temperature: Between 20°C (68°F) to 30°C (86°F)

Please use the suggested furnace cycle/program shown below:

Step heating	Heat rate	X1	X3	X6	X9
Room temperature to 260°C (500°F)	3°C (37.4°F)/min				
Holding time at 260°C (500°F)*		40 min.	60 min.	70 min.	90 min.
260°C to 580°C (500°F to 1076°F)	6°C (42.8°F)/min				
Holding time at 580°C (1076°F)*		20 min.	30 min.	40 min.	50 min.
580°C (1076°F) to end temperature**	9°C (48.2°F)/min				
Holding time at end temperature*		30 min.	40 min.	50 min.	60 min.

Important Notes:

* If more rings are inserted together in a furnace, each holding time should be extended by 10 minutes.

** Use the end temperature indicated by the investment and restoration material being used, as metal (casting), ceramic pressing, etc.

Limits:

- Posteriors: Do consider indications of material to be done, as metal (casting), ceramics pressing, etc.
- Minimum thickness: 1mm occlusal and cervical 0.6mm
- Minimum diameter connection: Anteriors: 3.6mm and Posterior: 4.4mm

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